LACTOSCAN COMBO Operation Manual

Accurate Somatic cells counting & Ultrasonic milk analyzer in one device Simple. Fast. Accurate. Reliable. Portable.



Available online at http://lactoscan.com/editor/ufo/manuals/COMBO/IM_COMBO_EN.pdf

Operation Manual of LACTOSCAN Combo is made to help better understanding operation of the device.

The Manual consists of step-by-step instructions for setting, use and maintenance of the device, preparation of samples, collection and processing of the received data.

The information in this manual is described as precise as possible and is applicable for the latest hardware and software versions of the LACTOSCAN COMBO. However, it can be changed without any preliminary consultation or notification.

Last edited: 20.06.2019

Safety measures

- 1. Always make sure that the power supply of the incoming current matches the supply in your region.
- 2. For "Working environment conditions" see page 11.
- 3. Always make sure the main switch is in position "OFF" before connecting it to the electrical network.
- 4. Never touch the power cable with wet hands.
- 5. Always make sure that the grounding clamp on both the device and the wall outlet are well connected. The power cable must be connected to the grounding.
- 6. In order to prevent a possible electrical shock, make sure, that the supply is well grounded.
- 7. In case of damaged device, turn it off and contact the authorized service. Do not disassemble the device. If you do so, the warranty will not be considered valid.
- 8. Use only the authorized accessories (USB flash drive, LACTOCHIP, wireless keyboard and mouse).
- 9. Use this device only the way that is indicated in this Operation Manual and in the documents related to its components. Each inappropriate use of the device is greatly concerning and may cause damages as the mentioned in the warning notes.
- 10. In case of smoke coming out of the device, turn in off and contact the authorized service.

2

Safety symbols

The following symbols can be found on the device and in this document. Please, give them special attention and always use the device in the safest possible way.

Symbol	Meaning
	Attention and Warning
	Device is covering all the provisions of the EU
	Protection from change of the pole of the supply
⊕ <u></u>	Plus and minus of the power receptacle

Warnings

Object	Warning
Embedded tablet battery	There is a risk of explosion if the battery is being changed improperly. This battery not to be changed by the customer. Please contact the authorized service.
Box	Do not disassemble the box. There are not adjustable components inside the box. If there is an issue, please contact the service.
Guidance	Do not try to service the device, unless this instruction is perfectly understood and discussed.
	If you do not pay enough attention to this instruction, this may lead to harm to the service supplier, electric shock during operation, mechanical or other risks.

Waste	After using the LACTOCHIP it is necessary to throw it in the appropriate bin.
Operator	Must have the basic knowledge on the procedure on cell counting.

Content:

Unpacking of LACTOSCAN COMBO	
Content of the kit	
LACTOSCAN COMBO	
I. Somatic cell counter (SCC) of LACTOSCAN COMBO	
Introduction of the Somatic cell counter	
Working principle of SCC	
Description of LACTOSCAN SCC KIT	
Description of LACTOCHIP	
Description of the software displays	
Preparation for work with LACTOSCAN SCC	Error! Bookmark not defined.
Preparation of sample for analysis	Error! Bookmark not defined.
Sample analysis	Error! Bookmark not defined.
Data samples entry during the analysis	Error! Bookmark not defined.
Reference sample testing	Error! Bookmark not defined.
Saving image from the analysis as .BMP file	Error! Bookmark not defined.
Setting the analysis' parameters	Error! Bookmark not defined.
Deleting results from the data base	Error! Bookmark not defined.
Cleaning of the Database	Error! Bookmark not defined.
Printing the results from the data base	Error! Bookmark not defined.
Printing the results on an external printer:	Error! Bookmark not defined.
Saving the results in PDF file:	Error! Bookmark not defined.
Saving the results in Excel file:	Error! Bookmark not defined.
Focusing of LACTOSCAN SCC	Error! Bookmark not defined.
Work in Advanced mode:	Error! Bookmark not defined.
Exit and switch off the LACTOSCAN SCC	Error! Bookmark not defined.
Maintenance and cleaning	Error! Bookmark not defined.
LACTOSCAN SCC - determination the number of the some	atic cells in raw milk - method
	EITUI: DUUKIIIAIK IIUT UETIIIEA.

Detection of Yeast mastitis in raw milk samples using LACTOSCAN SCC / LACTOSCAN COMBO
II. ULTRASONIC MILK ANALYZER135
Specification of LACTOSCAN COMBO's Ultrasonic milk analyzer
1. FUNCTION
2. TECHNICAL PARAMETERS
2.1. Working modes characteristics:
2.1.1. Measurement mode milk / dairy product – first type 139
2.1.2. Measurement mode milk / dairy product – second type 139
2.1.3. Measurement mode milk / dairy product – third type 139
2.1.4. Cleaning
2.2. Measuring range:140
2.3. Accuracy:
2.4 Correct ambient conditions:
2.5. Continuous working time:
2.7 Milk sample volume per one measurement:
3. Qualification of raw milk, thermally treated milk, other dairy products and derivatives
3.1. Taking samples and preparation for analysis143
3.2. Making the measurement143
3.2.1. Preparing the analyzer for working mode143
4. CLEANING THE ANALYZER147
4.1. Periodically cleaning (rinsing) the analyzer147
4.1.1. Periodical cleaning frequency
4.1.2. Making the rinsing
4.2. Complete cleaning
4.2.1. Complete cleaning frequency148
4.2.2. Cleaning
5. Menus' description153
5.1 Menu File 153
System Setup

Reset System	153
Exit	153
ShutDown	153
5.2 Menu Supplier	153
Add Supplier	153
Edit Supplier	153
5.3 Menu Tools	154
LSCal	154
Calibration	154
Program control:	154
5.4 Menu LS-MENU	156
Special Modes	156
5.5 Menu Correction	156
Temperature Correction	156
All Corrections	157
Corrections	157
5.6 Menu Settings	158
Set Calibration Names	158
Save Calibration	159
Restore Calibration	159
Autoprint Enable/Disable	159
Format Results:	160
High Fat Speed For Calibration	160
5.7 Menu Tests	161
Test Pump	161
Ultrasound	161
Set Amplitude	162
5.8 Menu Cleaning	163
5.9 Menu pH & Co Meter	163
pH Meter En/Dis	163
pH Test	163
pH Measure	163
pH Meter Calibration	163
Co Meter En/Dis	164
7 Operational manual © 2017 "Milk All rig	otronic" Ltd hts reserved

Co Meter Test
Co Meter Calibration 16
Sample calibration
5.9 Menu Pump Control 16
Pump In
Pump Out
Stop 16
5.11 Menu Measurement
Discard Last measurement 16
5.12 Menu Help 16
LS Identity 16
About16
Admin Mode 16
5.13 Menu Admin Mode 16
Restore Factory Settings 16
6. Peristaltic pump service
Peristaltic pump
7. Possible malfunctions and error messages, troubleshooting
8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
8. Making corrections and recallibration of the device 17 8.1. Taking samples and preparation of samples for checking the accurracy of the milk analyser, making corrections and recalibration 17 8.2. Determination the type of the discrepancy: 17 8.2.1. Making measurements 17 8.2.2. Analysing the measurement results 17 8.3. Making corrections, limits and changing steps 17
8. Making corrections and recallibration of the device 17 8.1. Taking samples and preparation of samples for checking the accurracy of the milk analyser, making corrections and recalibration 17 8.2. Determination the type of the discrepancy: 17 8.2.1. Making measurements 17 8.2.2. Analysing the measurement results 17 8.3.1. Possible corrections, limits and changing steps 17 9.Basic printer characteristics Error! Bookmark not defined
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device
 8. Making corrections and recallibration of the device 17 8.1. Taking samples and preparation of samples for checking the accurracy of the milk analyser, making corrections and recalibration 17 8.2. Determination the type of the discrepancy: 17 8.2.1. Making measurements 17 8.2.2. Analysing the measurement results 17 8.3.1. Possible corrections, limits and changing steps 17 9.Basic printer characteristics 9.1 Printer control panel 9.3 Consequence of actions when printing 18 10.1. Connecting to 12 V DC power supply.

Appendices	186					
Appendix 1: Preparation of samples for milkanalysers' calibration						
Appendix 2. Freezing point determination191						
1. Methods for determination	191					
2. The basic freezing point	191					
Appendix 3. pH measuring	194					
1. General information	194					
2. pH Electrode	194					
2.1. Electrode part	194					
2.2. Electrode care & Electrode maintenance	194					
2.3. Storage	194					
2.4. After Use	194					
2.5. Electrolyte Replacement (for refillable electrode only)	195					
2.6. New electrolyte preparation:	195					
2.7. Re-use the electrode	195					
2.8. Electrode cleaning	195					
2.9. Electrode activation	196					
2.10. Rejuvenation Procedure	196					
2.11. Electrode Lifespan	196					
3. Buffer Solutions	197					
4. pH Electrode Calibration	197					
5. PH helpful hints	201					
6. Preparation for pH measurement	202					
7. PH measuring	207					
Appendix 4. Conductivity measuring	207					
1. Method of determination	207					
2. Conductivity measurement	208					
3. Corrections in conductivity measurement.	208					
4. Conductivity calibration buffer preparation	209					
Warranty card	212					

Service card	
Contacts	

Unpacking of LACTOSCAN COMBO

Unpack carefully the LACTOSCAN Combo following the sequence illustrated below:



5. Automatic pippets

Content of the kit

The Somatic Cell Counter LACTOSCAN Combo is sent with the following items:

Object	Amount
LACTOSCAN COMBO	1
LACTOSCAN COMBO STARTER KIT	1
LACTOSCAN COMBO 14V power cable	1
LACTOSCAN COMBO Operation Manual	1
Automatic pipette from 10 to 100 μ L – set at 8 μ L	1
Automatic pipette from 20 to 200 μL –set at 100 μL	1
USB wireless keyboard with batteries	1
USB wireless mouse with batteries	1
Mini Vortex	1
Plastic sample holder	2
Spare Pipes	1
Alkaline cleaning solution Lactodaily	1
Acidic cleaning solution Lactoweekly	1

When receiving the device, please make sure it is not damaged during transportation. Also make sure that all the parts, including the above mentioned, are in the box. In case of damage/missing parts, contact the carrier company. The warranty does not cover damages caused during transportation.



Disregarding and not taking off all the transport brackets, Styrofoam may lead to a damage of the device.

The transport brackets, Styrofoam, must be placed before transporting, in order to prevent a damage.

Device registration:

Each customer can register the LACTOSCAN COMBO by visiting <u>www.lactoscan.com/device registration</u>. The serial number is needed to be filled in, name of the owner and contacts data. After registration, notifications will be received about software updates and information about new analysis, which may be done with the LACTOSCAN COMBO.

Purpose of the device:

LACTOSCAN COMBO is designed to be used only for analysis of milk in milk collecting centers, laboratories, etc. It is not manufactured for testing human milk or animal but with therapeutic or diagnostic usage.

Environment conditions

In order to work properly and stable for a long period, LACTOSCAN COMBO needs to be installed in a room meeting the following environmental conditions:

- Indoor temperature 10-35 °C. It is not recommended to place it in a room where the temperature is below 4 °C.
- Prevent exposing the device to direct sun light.
- Prevent exposing to a continuous vibration.
- Relative humidity 0-95%.
- In a place free of corrosive gazes or other corrosive substances.
- In a place free of dust.
- 10 cm (4 inches) is the minimum required distance for normal air flow.
- Do not place any heavy object on the device.

LACTOSCAN COMBO

The front and upper panel, as well as the external accessories, are shown below.

Front panel:



Entrance of the carriage. It opens automatically and ejects a stand for the LACTOCHIP with the analyzed sample.

Touch-screen display. It operates the device. The necessary functions and results from the analysis are displayed there.

Embedded printer. It allows printing the results after analysis or the results of previous tests.

pH probe place. To place the pH probe inside.

ON/OFF Button. It is used to turn on and off the display.

Wireless keyboard and mouse. Make entering data easier.

Back panel:



Switch ON/OFF for the SCC. It is used to power on / off the SCC part of the device.

12 ÷**14** V **input, min 3** A. Power supply. Always connect only the powering adapter, supplied with the device.



Attention!

Do not connect the power supply of the device if the toggle switch ON/OFF is not on "Off" position.

USB ports. For connection of USB flash-drive and wireless keypad, mouse and external printer.



Attention!

Do not disconnect USB flash drive while the device is working.

Serial interface RS232. Used only in service center. Please do not connect something to it.

Switch ON/OFF for the milk analyzer. It is used to power on / off the SCC part of the device.

pH probe connector. It is used to connect the pH probe cable towards the device.

HDMI. It is used to connect an external monitor towards the device.

Specifications of LACTOSCAN COMBO:

CHARACTERISTICS	Туре:	Somatic Cells Counter & Ultrasonic milk analyzer					
	Dimensions (HxWxL):	39 cm x 30 cm x 26 cm					
	Weight:	6.000 kg					
	Working power:	100-230 VAC, 2.5A, 120V					
	Frequency:	50/60 Hz					
	Power supply:	12 VDC, 5A, 60W					
	Working place:	Inside only!					
	Working temperature:	15-35 ℃					
	Working humidity:	0-95%					

LACTOSCAN COMBO Somatic cell counter and ultrasonic milk analyzer in one device.

The Operation manual consists of 2 separate descriptions for the both devices. Here's how to switch between the Somatic cell counter and the Ultrasonic milkanalyser:

If you are in the SCC, by pressing the icon in the red rectangle with picture of the device will switch to the milk analyser:

	Main	🔟 Graph	Options	Advanced				S × 18-	10 17:04
			No Image				Start		
						ID			
						Milk type			
						Average cells per			
						me			
	ID	Chip	Date	Client	Milk	Dilution	Cells per mL	Avg. cell size [µm]	
	Empty								Details Summary
wv	vw.Lactoscan.com								2016-08-10 17:04

If you are in milk analyzer's mode:

File	Suppliers	Tools	LS-Menu	Measurement	SCC	Help	53 🗸	
	Supplier Number:				^			
	FAT 00.	00	SN	NF 00.00		DEN 00.00	LAC 00.00	

I. Somatic cell counter (SCC) of LACTOSCAN COMBO

Introduction of the Somatic cell counter

The LACTOSCAN COMBO's SCC is based on fluorescent microscope technique of counting cells. Thanks to the fluorescent dye, LED optics and CCD filming Technologies, the analysis of the milk is precise, reliable and fast.

In order the somatic cell to be counted by with SCC, the sample is mixed with SOFIA GREEN dye. It is only 8 μ L that is needed to be pipette onto the single LACTOCHIP. After that, the CHIP is loaded in the device. The analysis is being conducted during a period between 10 seconds and 2 minutes and the duration depends on the number of filmed fields. The system of SCC is focusing automatically on the CHIP and the dyed cells are filmed by the sensitive CCD camera. This algorithm of analysis of digital images determines the number of the fluorescent cells and counts their concentration and size. The result is automatically displayed.



Working principle of SCC

Unique, 3D, multi-image, patent application protected, sequential scanning process, based on a precise fluorescent optics and low magnification, images analysis software, LACTOSCAN SCC is fast, precise and reliable counter of somatic cells. Via automatic displacement of the mechanism on axles X – Y and liquid lens Z, the device is capturing maximum 60 images. After capturing, the images are being processed by the embedded software and the average result, calculated by using the formula from IDF/ISO 13366, of all the filmed images is displayed. The whole process, after placing the LACTOCHIP in the cartridge, is automatic.



SCC compared with a standard methods for Somatic Cells Counting

The below described data are based on validation and check of SCC done by the "Biotechnologies" body of University "Prof. AsenZlatarov", Bulgaria

1. Accuracy of analysis – comparison between SCC and Direct microscopic counting of somatic cells (DMCSC)

Images of somatic cells (SCC) Images of somatic cells (DMCSC)



19 Operational manual

© 2017 "Milkotronic" Ltd All rights reserved



SCC

Standard method

NSC (cells/ml)	*CV %	NSC (cells/ml)	*CV %
100 000	5%	100 000	7%
500 000	3%	400 000	5%
1 000 000	2%	600 000	4%

*Coefficient of variation

2. Compatibility – ratio between data from SCC and DMCSC



Description of LACTOSCAN SCC KIT

Types of LACTOSCAN SCC KIT:

- LACTOSCAN SCC STARTER KIT 100 tests consists of 25 disposable LACTOCHIP x4, 100 micro test tubes with SOFIA GREEN lyophilized dye and 200 tips for automatic pipette. The expiration date of the set is one year after its manufacturing date.
- LACTOSCAN SCC KIT x4 400tests-consists of 100 disposable LACTOCHIP x4, 400 pcs micro test tubes with SOFIA GREEN lyophilized dye and 800 tips of automatic pipette. The expiration date of the set is one year after its manufacturing date.

Description of LACTOCHIP

LACTOCHIP is a disposable chip, specially designed to be used with LACTOSCAN SCC. It is made of ABS and PMMA.



LACTOCHIP x4 has four separate closed chambers (A, B, C μ D), that allow analysis of four different samples. Each chamber's capacity is 8 μ L. LACTOSCAN SCC makes 16 pictures of the sample in each chamber and then it performs the analysis via a specific algorithm.



LACTOCHIP 4L



Attention!

Before starting work with LACTOSCAN SCC KIT always check LACTOCHIP letter. It can be L or R. Letter L indicated LACTOCHIP 4L and letter R indicated LACTOCHIP 4R. The LACTOCHIP letters match with the already put in the program LACTOCHIPs data. One LACTOSCAN SCC KIT included only LACTOCHIPs with letter L or only with letter R.

Specifications:

LACTOSCAN COMBO'S SCC			
	Туре:	Somatic Cells Counter	
	Dimensions (HxWxL):	25.5 cm x 38 cm x 30 cm	
CHARACTERISTICS	Weight:	7.300 kg	
	Working power:	100-230 VAC, 2.5A, 120V	
	Frequency:	50/60 Hz	
	Power supply:	14 V, 5A, 60W	
TECHNICAL SPECIFICATION	Working place:	Inside only!	
	Working temperature:	15-35 ℃	
	Working humidity:	0-95%	
	Time of analysis:	From 10 sec to 2 min in dependence of the number of pictures taken	
LACTOSCAN SCC STARTER KIT 100 tests			
PHYSICAL CHARACTERISTICS	Dimensions (HxWxL):	355 mm x 255 mm x 135 mm	
	Weight:	1.500 kg	
	Number of LACTOCHIP x4:	25 pcs	
	Number of SOFIA GREEN lyophilized dye:	100 pcs	

LACTOSCAN SCC KIT x4 400 tests

	Dimensions (HxWxL):	355 mm x 255 mm x 135 mm
	Weight:	1.500 kg
Physical characteristics	Number of LACTOCHIP x4:	100 pcs
	Number of SOFIA GREEN lyophilized dye:	400 pcs
	Number of tips for the automatic pipette:	800 pcs
LACTOCHIP x4		
	Material:	ABS and PMMA
	Dimensions (HxWxL):	2 mm x 25 mm x 75 mm
Physical characteristics	Weight:	0.005 kg
	Volume of the camera:	8µL

Description of the software displays

LACTOSCAN SCC software consists of the following displays:

• Display "Initial"



Device's ID	Information about the software versions of the device
Exit button	By pressing this button the operator switches to menu Switch off/Operational system, from where whether to turn off the device or the program can be chosen.
	By pressing this button the User manual for LACTOSCAN SCC is opened in PDF file.
Summary	Information about the condition of the device after self- diagnosing test.
Continue	By pressing this button the operator continues to the main display of the software.
🛫 сомз	COM 3 must be selected, in order the device to work

properly.

PCB Device PCB Device must be selected, in order the device to work properly. For demonstrations choose «Simulator», this way only the software will be working, without the device itself.



Microscope Device must be selected, in order the device to work properly. For demonstrations choose «Simulator», this way only the software will be working, without the device itself.



Main display

Main buttons	Switch the main displays of the software
Functional buttons	Allows the user to mark and unmark cells, to select cells, to zoom the image by pressing on the buttons.

Shows software edited image on which mastitis cells are colored in green, yeast cells in blue and the background is presented completely smooth.







Software marks all the counted cells. Mastitis cells are marked with red square, yeast cells in blue square.



Shows the real image





Allows user to select whether to see only mastitis cells, yeast cells or all cells and to zoom the image:

Marked only mastitis cells

Software edited image



Marked only yeast cells

Software edited image



Real Image



Real Image



Marked all cells

Software edited image



<u>Zoom</u>

Real Image





100% Original size 200% Enlarged 2 times 300% Enlarged 3 times 400% Enlarged 4 times 500% Enlarged 5 times



Opens the image in full screen and returns the image to Main mode.

	Transfers the image in Advanced mode where the user can examine it in details.
Start	Starts the sample analysis
٩	Turns off the device
	By pressing this button the User manual for LACTOSCAN SCC is opened in PDF file.
	Switches to print page
	Shows a report from the made analysis in PDF format.

Calla Tast Ltcc - Lateau forms Cer County 1.2.05	
Reference	
5000011822-4451 PM Stri 10 System Cree mit Gladen 1.0	
Result	
16 snapshots, covering 1.87 ut, Average cet size 7.6 um Cells per mt, 294e+3	
Example picture	
000000001-0-05.11#	
Cell size histogram	
Number of sales	

Image check	Shows the density of the somatic cells in the selected sample
Results field	Summary of the selected analysis
	The flags show the following animal health status:
	1. 🖳 - Sub-clinical mastitis
	2. 🖳 - Clinical mastitis
	3. 😲 - Presence of yeast cells in milk sample. Possible yeast mastitis. A test for yeast mastitis must be performed.
	Opens the search panel
Summary	Switches to a table with summarized information about the performed analysis
Details	Switches to a table with detailed information about the already performed analysis
30 Operational ma	nual © 2017 "Milkotronic

• Display "Options"

Allows the operator to set options of analysis according to the local or regional requirements for the different types of milk.



Options display	
Main buttons	Switch main displays of the software
Parameters of measurement	Sets the min. and max. size of the counted cells for the different milk samples
Precautious levels	Sets a precautious level of sub-clinical and clinical mastitis according to the animal species
Laboratory offset	Coefficient for corrections
Apply and store	Applies the set options and their storage
Set default options	Sets default options of the measurement parameters and precautious levels
Advanced option	Opens menu with options needed for service purposes only
Focusing	Pressing the button switches to Setting focus menu

• Display "Graphics"



Main buttons	Switch main displays of the software
Summary	Summary of the data used to create the graphics
	Switches the source of information used to create the graphics
Min and Max size of the cell fields	Sets the border size of the cells included in the final analysis result by pressing the arrow buttons
Graphics	Graphical representation of the analysis data. "Size of the cell graphic" presents the data for the number of the cells according to their size. "Sample graphic" shows the size and number of cells that are taken into consideration during counting the analysis results.
Average size	Average size of the cells
Green line	Shows the average size of the cells
	© 2017 Milliotropie" Ite

Shows an image which is not included when calculating the final result. In the table of details for the corresponding result, the image is marked with NO in column VALID, because of unacceptable deviation from the average of results derived from other images.



Display "Advanced" • Main buttons Full screen button Advanced 📕 Graph 22.07.2017 13:31 ч. Options Functional buttons Test sectors 3D mode sD 🔿 Image check Open_button Save button 🗖 Open. Sav 🕒 Stream 3 11 2 10 Stream button Probe type Preserved cow milk Park button ime 0.234 seconds Eject button or 1.000 ÷ LACTOSCAN 4R 50µm 4x16 LACTOCHIP type and 2nd lights power 50.0% parked 454 Microfluid camera position Absolute cell count Select microfluid camera sectors and capture image Brightness settings Milk sample data

Red column

Advanced display	
Main buttons	Switch the main displays of the software
	<text></text>
Functional buttons	Allows the user to mark and unmark cells, to select cells, to zoom the image by pressing on the buttons.
	Shows software edited image on which mastitis cells are colored in green, yeast cells in blue and the background is presented completely smooth.



Software marks all the counted cells. Mastitis cells are marked with red square, yeast cells in blue square.



Shows the real image





Allows user to select whether to see only mastitis cells, yeast cells or all cells and to zoom the image:

Marked only mastitis cells



Mark only mastitis cells Mark only yeast cells Mark all cells Fit to window 100% Original size 200% Enlarged 2 times 300% Enlarged 3 times 400% Enlarged 4 times 500% Enlarged 5 times

Software edited image



Marked only yeast cells Software edited image



Real Image



Real Image


Marked all cells

Software edited image



Real Image



Zoom



Test sectors	Allow user camera whi mode.	to selec ich he v	t the ex wants to	kact sector test and	from exami	the mi ne in a	crofluid Idvance
3D	Examine	the	test	image	in	3D	mode.

	Main Graph Options Advanced 23.7.2017 r. 08:21 4. Image: Construction of the second sec
Image check	Shows the density of the somatic cells in the selected sector.
P Open	Opens an image from previous testing for further examination in Advance mode.
Save	Saves the shown image in the Image check field.
	Parks the cartridge with the loaded LACTOCHIP.
$\overline{\bigcirc}$	Ejects the cartridge with LACTOCHIP.
LACTOSCAN 4R 50µm 4x16	Selects the LACTOCHIP type.
Microfluid camera parked	Shows the position of the cartridge with the microfluid camera / LACTOCHIP.
Absolute cell 454 count	Shows the counted cell on the image from the selected test sector.
	By using the arrows, the user can select the test sector of the microfluid camera and then by pressing capture button to make a photo of the sector.
Exposure time 0.234 seconds	Allows making manual or automatic changes of brightness, light power and focus.
Dilution Probe type 1 Preserved cow milk	Entering information about the dilution and probe type of the examined sample.

Preparation for work with LACTOSCAN SCC

Working sequence when preparing the LACTOSCAN SCC for work is described below:

1. Place the device on a flat, dry surface with nothing around it (4 inches radius)

2. Make sure the power supply is equal to the necessary for the equipment and connect the cable to the inlet on the back panel of the analyzer. Using improper power supply may lead to serious damage of the device.

3. Place the plug into the electricity network.

4. Switched the ON/OFF key, placed at the back panel of the analyzer and wait a minute until it charges.

5. Turn on the display by using the ON/OFF key placed at its upper side.

6. LACTOSCAN SCC launches the software and carries out an independent diagnostic test, including all the optic components. If an issue is found during the check with any of the device's components, the program would not allow continuing to the Main screen menus and indicates which component is problematic.



The issues may be caused by disconnection of cables, lack of power supply or software failure during the booting process. To eliminate them it is necessary to update the program by pressing the Retry button (try again).

7. Press button Help to open an instruction describing how to check the connections.



If the issue is not resolved, please contact your local distributor or <u>www.lactoscan.com/service</u>.

8. Switching to Main display menus, by pressing the button "Continue".



Preparation of sample for analysis

To prepare a sample for analysis are needed:

- Raw milk;
- LACTOSCAN SCC KIT;
- 2 pcs. automatic pipettes.



Attention!

To obtain results with reference accuracy, the requirements for proper sampling specified in ISO 707 | IDF 50 must be met.



Attention!

Measurement accuracy depends on the correct and consistent implementation of all stages of sample preparation and good mixing of the sample. To minimize the difference in reporting the results of several tests on the same sample of milk, always stir well before taking the sample. The sample is representative only when stirred well.

Stages:

1. Preparation of the raw milk:

It is mandatory a raw milk, just milked or preserved with room temperature 15-25°C. The necessary min. volume is 30 mL.



Attention!

According to **INTERNATIONAL STANDARD ISO 13366-1 IDF 148-1:2008**, to obtain the best results, you should observe the following principles:

If the samples are without preservative, they should be measured within 6 hours after milking.

If samples cannot be measured in the course of these 6 hours, they must be preserved with Bronopol ($C_3H_6BrNO_4$), Potassium Dichromate ($K_2Cr_2O_7$) or Formalin (CH_2O) in amounts specified in the standards for sampling for analysis. The final concentration of Bronopol shall not exceed 0,05 g per 100 ml of test sample. The final concentration of potassium dichromate shall not exceed 0,1 g per 100 ml of test sample. They may be stored in a refrigerator at 4 °C ± 2 °C for no longer than 6 days.

We recommend using Bronopol as a milk preservative!

Before measuring samples, they must be heated up to 40 °C and cooled down to 20 °C and then stirred thoroughly with Vortex mixer. When samples are kept in a refrigerator, the fat globules float to the top and majority of the leukocytes adhere to them, and therefore the somatic cells go up together with fat globules. Often if the sample is not heated up to 40 °C and cooled down to 20 degrees, it is not possible to be mixed thoroughly, leading to uneven distribution of somatic cells in the sample volume. Then the measurements will vary.

Preserved and stored in a fridge sample is suitable for measuring no longer than 5-6 days.

The sample must be no more than 50 ml, and must not fill the bottle with the sample to the cap in order to allow easier mixing with Vortex mixer or by hand.

Attention!

/!\ If the analysis is not conducted within 3-4 hours after milking, it is necessary to preserve the milk. When the raw milk is preserved is recommended to use preservatives formalin, bronopol or potassium dichromate.

Attention!

If the preserved milk is chilled below 10 °C, it must be tempered naturally to room temperature 15-25 °C. Freshly milked milk is not necessary to be chilled or heated.



Attention!

Do not use for analysis raw or preserved milk with acidity above:

- 18 °T (Therner) for cow milk
- 17 °T for buffalo milk
- 16 °T for goat milk
- 22 °T for sheep milk



Attention!

In case of measurement of milk with Fat more than 5 %,

for

35 -

25

15 10

example buffalo milk, it is necessary the milk sample to be diluted with water in ratio 1:1. Then 100 µL of it is taken and added to the lyophilised dye. Adding water prevents difficulties in milk samples entry into the microfluidic chamber.

Using the Mini Vortex mixer stir the raw milk sample. For stirring, place the tip of the container in the stirrer, press and keep it pressed for 1-2 seconds, then remove it (see 1, 2, 3.). Repeat it 3-4 times paying attention during the stirring process the sample not to reach the cap of the container.



2. Pipetting 100 μL raw milk in micro-tube with SOFIA GREEN lyophilized dye:

Take one micro tube containing SOFIA GREEN liquid dye, open it, and place it on the rack.

Attention!

Before start working with the automatic pipettes, read carefully their Instruction Manual. Make several trials sucking and pipetting water in order to understand when exactly the first and second stop of the working button is reached.

Take preliminary set to 100μ L automatic pipette. Make sure that the front cone of the pipette is clean. Place it vertically over one of the tips on the working rack and place the cone of the pipette in the opening of the tip by slightly pressing it.

From initial position (see 2.1), press the working button of the pipette till the first stop is reached (see 2.2), keep it pressed and dip 2-3 mm of the tip in the milk (see 2.3). Smoothly release the working button and take out the tip from the liquid. Touch the walls of the bottle to remove the excessive milk (see 2.4, 2.5).



Pipette milk in the opened micro-tube on the rack by smoothly pressing the working button of the pipette from the initial position to the first stop (see 2.6, 2.7). After a short period press the button to the second stop (see 2.8). In this way you'll empty the tip and will guarantee precise pipetting. Always pipette the milk without dipping the tip in the lyophilized dye.

Release the button to its initial position (see 2.9).



3. Stirring the sample:

Close the micro-tube containing SOFIA GREEN dye and milk sample. Take it from the rack and place the tip of the micro-tube in the opening of the stirrer Mini Vortex. Press and hold it pressed for 1-2 seconds and remove. Repeat 8-9 times being careful while stirring the solution not to reach the cap of the micro-tube. (see 3.1, 3.2 and 3.3)



4. Interaction of milk with dye:

1 minute is needed for this interaction. If it is less 1 minute or more than 20, the analysis result may be with deviation 2-3%.

5. Repeated stirring the sample:

Take the micro-tube containing the sample from the rack and place its tip in the opening of the stirrer Mini Vortex. Press and keep it pressed for 1-2 seconds, remove. Repeat 2-3 times, paying attention place the tip of the container in the stirrer, press and keep it pressed for 1-2 seconds, then remove (see 3.1, 3.2, 3.3.). Repeat it 3-4 times paying attention during the stirring process the sample not to reach the cap of the container.

Attention!

If more than 5 minutes elapse after the milk was placed inside in micro-tube with SOFIA GREEN dye, stir with Mini Vortex mixer the sample once again, before filling in the LACTOCHIP.

6. Pipetting 8 µL sample in the micro-fluidic camera of the LACTOCHIP x4:

Open one LACTOCHIP x4.

Attention!



Do not touch the upper surface of the LACTOCHIP x4. Always hold its side edges.

To eject the solution into the micro-fluidic camera of LACTOCHIP x4, take preliminary set to 8 μ L automatic pipette.

Make sure that the front cone of the pipette is clean. Place it over one of the tips on the rack and place the cone of the pipette in the opening of the tip by slightly pressing it.

Open the micro-tube containing the solution.

From initial position (see 6.1), press the working button of the pipette till the first stop is reached (see 6.2), keep it pressed and dip 2-3 mm of the tip in the solution (see 6.3). Smoothly release the working button to the initial position. Take out the tip from the liquid by touching the walls of the bottle to remove the excessive solution (see 6.4, 6.5).



Now there's 8 μL solution in the tip.

Take the LACTOCHIP x4 by holding its side edges.

Pipette the solution at an angle of approximately 80[°] to the filling opening in semicircular shape. Pipetting is done by smooth pressing the working button of the pipette from the initial position to the first stop (see 6.6, 6.7). Hold the button at the first stop, remove the pipette from the LACTOCHIP and smoothly release the button to the initial position (see 6.8).

In this way you'll empty the tip and precise pipetting will be assured. Release the button to its initial position.

Attention!

Do not use the second stop in order to avoid air entering into the micro-fluidic camera.

Attention!

Avoid forming bubbles in the micro-fluidic camera and splashes during pipetting the sample.



In order to load the rest of the micro-fluidic cameras of the LACTOCHIP x4 / LACTOCHIP x2, repeat the procedure described in points 1 to 6 by consecutively filling the micro-fluidic cameras from A-D.



Attention!

It is recommended to use the all microfluidic cameras at once. If you use only 1 or 2, store the LACTOCHIP x4, paying attention not to contaminate it with dust or other pollution as it will lead to false results of the analyses.



Attention!

It is recommended to place the loaded LACTOCHIP x4 in the device and to start analyses within 1 minute. Delay may lead to inaccurate results due to evaporation of the sample and air entering it.

7. Starting analysis:

Place the loaded with sample LACTOCHIP x4 in the cartridge of the LACTOSCAN SCC. Using the software, start the analysis.

8. Disposal:

Using the button for removing the tip, leave the tip inside the micro-tube with the sample. Dispose the micro-tube with the sample residue, the tip and used for analysis LACTOCHIP x4 in suitable container.

Sample analysis

1. Press on the Main screen in order to proceed with entering data of the sample.



2. Choose the type of chip from the drop down menu

LACTOSCAN 4R 50µm 4x16

$\mathbf{\Lambda}$

Attention!

Before starting work with LACTOSCAN SCC KIT always check LACTOCHIP letter. It can be L or R. Letter L indicated LACTOCHIP 4L and letter R indicated LACTOCHIP 4R. The LACTOCHIP letters match with the already put in the program LACTOCHIPs data. One LACTOSCAN SCC KIT included only LACTOCHIPs with letter L or only with letter R.



In our example we choose $\boxed{\text{LACTOSCAN 4R 50 \mu m 4x16}}$ for LACTOCHIP 4X16 with thickness of the four micro fluidic chambers from 50 µm.



3. All the micro fluidic chambers of the LACTOCHIP are selected. If you want

to use only 1 or 2, using button not be using in this test. deselect the chambers which you will

In our example we deselect Sample C and D.



4. For a quick analysis, press button and select a test from the dropdown menu. You can select from test "Scanning 25% of positions", test Scanning 50% of positions", test "Scanning 100% of positions".



5. Complete the value of the dilution's multiplier in field 1 under Dilution.



Attention!

The value of the multiplier for dilution depends on the fat content of the milk.

For milk with fat:

- Up to 5 %, the value of the multiplier for dilution is 1, as the milk is added directly to the lyophilized dye in the micro tube.

- Above 5 %, the value of the multiplier for dilution is 2, because, for example, in case of measurement of milk with Fat more than 5 %, for example buffalo milk, it is necessary the milk sample to be diluted with water in ratio 1:1. Then 100 μ L of it is taken and added to the lyophilised dye. Adding water prevents difficulties in milk samples entry into the microfluidic chamber.



Attention!

Don't forget to change the multiplier for dilution in the field "Dilution" when working with milk, containing fat over **5%** and add distilled water to the lyophilized dye.



6. Choose the type of milk from the drop-down menu.



7. Enter the name of the supplier or choose from the drop-down menu.

CORE)A
Dilution	Milk	
1	Cow milk	÷
Client		
0000-000	1014	-
0000-0001	1014	^
Iva Penev	0131	
Iva Penev	01310	E
lva Penev	01311	
lva Penev	01312	
lva Penev	01313	
Iva Penev	01314	
Iva Penev	01315	
Iva Penev	01316	
Iva Penev	01317	
lva Penev	01318	
Iva Penev	01319	
Iva Penev	0132	
Iva Penev	01320	
lva Penev	01321	
lva Penev	01322	
lva Penev	01323	
Iva Penev	01324	-

8. After the information for all the samples is entered, press in order to proceed with the analysis. On the screen a message will be displayed, showing that the cartridge for placing the LACTOCHIP comes out of the door.

Pressing will revert to the Main screen.







Image



10. The device is making the analysis of the samples.



11. If you want to see test image in real time, press



12. In order to print the results, press

13. SAVE der to sav Cancel sults and continue the analysis of new samples,

press or press without saving the results. , to continue with the analysis of new samples



14. Take the LACTOCHIP out.



15. In order to go back to the Manual screen, press



Data samples entry during the analysis

The software of LACTOSCAN SCC allows the entry of data for new samples during the process of analyzing the samples.

In order to do this:



The data for the samples which are being analyzed at the moment are saved in the form for facilitating the user's work, in case analysis of identical samples for different clients or identical species of animals have to be made.

2. Enter the data for the new samples.



If

is grey, this means that the previous analysis has not been finished yet. In this

case wait for the analysis to be finished and the *screen* to become active, or go back to the screen with the current analysis.





4. Press , i

, in order to proceed directly to the new analysis.

By pressing , the results are automatically saved in the data base of the device. And the software prompts to pull out the LACTOCHIP with the samples from the previous analysis.



5. Pull out the LACTOCHIP with the samples from the previous analysis.



6. Press . The software prompts for placing the LACTOCHIP with the samples for which data are entered.





7. Place the LACTOCHIP with the samples with entered values.



8. Press in order to start the analysis.



Reference sample testing

To receive results with reference accuracy from somatic cell counting using LACTOSCAN SCC the following steps must be done:

1. Samples taking:

Samples must be taken with accordance to the recommended sampling method given in ISO 707|IDF 50.

2. Samples storage:

Samples must be stored with accordance to the recommendations in ISO 13366-1 | IDF 148-1:2008, page 4:

"7.1 Storage

Prior to testing or preservation, store the test samples at a temperature of 4 °C \pm 2 °C.

Analyse the test samples within 6 h after sampling. In the case of longer storage, add chemical preservatives such as boric acid, bronopol or potassium dichromate. The final concentration of boric acid shall not exceed 0,6 g per 100 ml of test sample. The final concentration of bronopol shall not exceed 0,05 g per 100 ml of test sample. The final concentration of potassium dichromate shall not exceed 0,1 g per 100 ml of test sample. Store the thus preserved test samples at a temperature of 4 °C \pm 2 °C for no longer than 6 days.

For environmental reasons, it is recommended to restrict the use of potassium dichromate to samples that require a long shelf life only."

Note! Bronopol is the recommended preservative for preparing samples for analysis with LACTOSCAN SCC.

3. Sample Preparation procedure recommended by IDF Reference laboratory:

Heat the sample in a water bath with the temperature between 38 and $40^{\circ}C$ ($100 - 104^{\circ}F$) for a period of 10 minutes. Shake vigorously for 2 minutes and then emerge the vial in the water bath for an additional 10 minutes. Cool the sample to the temperature of 20 °C.

4. Testing:

To do this the following consumables will be required: $1 \times X4$ cartridge, $4 \times eppendorf$ tubes with Sofia Green dye and $8 \times pipette$ tips.



• After sample preparation take 100μ L of sample and put it into an eppendorf tube. Do this a total of 4 times for the 4 cameras of the chip.

• Incubate for 5 minutes while mixing frequently.

• After incubation do a final mixing and take 8 μ l of the sample and put into the X4 cartridge – do this a total of 4 times, each time with sample from a new eppendorf tube.

The procedure for preparing single sample for testing is described in section **<u>Preparation</u>** of sample for analysis.

Note! A new tip must be used for preparing each of the 4 samples and a new tip must be used for putting the sample in each chamber A, B, C, D.

• Let the cartridge sit for 30 seconds before testing in order the cells to stop moving inside the chambers.

<complex-block>

 Mark
 Mark

In testing mode, press button

to select "Reference testing mode".

• Fill in the data for the tested milk sample.







• The device is making the analysis of the sample.



LACTOSCAN SCC performs the test and gives final result. The final result is an automatically calculated average result from the results of the tests of chambers A, B, C, D.

5. Coefficient of variation (CV) %:

Coefficient of variation (CV) % is measured in Reference mode following the instruction for Reference testing. It is determined by analysis of standard reference milk samples from **certified dairy laboratories**. To determine the reference coefficient of variation (CV) % only such standard reference milk sample must be used. If nonstandard milk samples are used, a deviation of the CV% may be received. The standard reference milk samples do not include a large amount of cell fragments and do not include yeast cells. The presence of yeast cells in milk can easily be determined by examining the graph of the analyzed milk sample. Below you will find examples for graph of normal milk without yeast cells and graph of milk with yeast cells:



Normal milk

Milk with yeast cells

6. Calculating Coefficient of variation (CV) %:

To calculate the Coefficient of variation (CV) %, following the above procedure, you have to test 10 times milk sample with somatic cell concentration 50000-200000 cell/ mL, to test 10 times milk sample with somatic cell concentration 400000-600000 cell/mL and to test 10 times milk sample with somatic cell concentration 800000-1400000 cell/mL. Then to fill in the received results in the table below.

Sample tables:

Empty	ed in						
Raw data	Sample 50000- 200000 cell/mL 1	Sample 400000- 600000 cell/mL 2	Sample 800000- 1400000 cell/mL 3	Raw data	Sample 50000- 200000 cell/mL 1	Sample 400000- 600000 cell/mL 2	Sample 800000- 1400000 cell/mL 3
Measurement 1				Measurement 1	142 000	571 000	1 153 000
Measurement 2				Measurement 2	161 000	563 000	1 164 000
Measurement 3				Measurement 3	147 000	548 000	1 146 000
Measurement 4				Measurement 4	157 000	530 000	1 142 000
Measurement 5				Measurement 5	142 000	543 000	1 146 000
Measurement 6				Measurement 6	145 000	584 000	1 122 000
Measurement 7				Measurement 7	151 000	558 000	1 100 000
Measurement 8				Measurement 8	150 000	584 000	1 154 000
Measurement 9				Measurement 9	140 000	539 000	1 120 000
Measurement 10				Measurement 10	145 000	570 000	1 155 000
Mean value				Mean value	148000	559000	1140200
Standard deviation				Standard deviation	6815	18708	19904
Coefficient of				Coefficient of			
variation %				variation %	4,60	3,35	1,75

Formulas for calculating:

Mean value – The all 10 measurements were summed and then divided into 10 (number of samples).

Standard deviation in statistics, typically denoted by σ , is a measure of variation or dispersion (refers to a distribution's extent of stretching or squeezing) between values in a set of data. The lower the standard deviation, the closer the data points tend to be to the mean (or expected value), μ . Conversely, a higher standard deviation indicates a wider range of values.

$$s = \sqrt{\frac{1}{N-1} \sum_{i=1}^{N} (x_i - \overline{x})^2},$$

Where:

xi is one sample value

 $\bar{\boldsymbol{x}}$ is the sample mean

N is the sample size

Coefficient of variation: CV, % = Standard deviation x 100/mean value.

The above table can be downloaded from: <u>https://autocellcount.com/page/coefficient-of-variation</u>

7. LACTOSCAN SCC setting:

For measuring standard reference milk samples from certified dairy laboratory and calculating the Coefficient of variation (CV) %, LACTOSCAN SCC must be set as follow:

- Minimum cell size in micrometer : 6
- Maximum cell size in micrometer : 20

The Menu Options allows the user to set the analysis parameters for different types of milk.

Saving image from the analysis as .BMP file

LACTOSCAN SCC allows Users to save an image from already done analysis as .BMP file. The .BMP file is easy for further process and can be opened by or imported in any Windows program.

To do this:

1. Select the analysis from the data base:

	Graph Options	Advanced			17.1.2017 r. 16	T Main	🕼 Graph 👘 O	ptions Advanced			17.1.2
			© SI © SI	ev2					in 282C Client Var	SHA213 rved cow milk	
									· · · · · · · · · · · · · · · · · · ·		
			Average cells per mL 3443	3·10°	Ģ				Average cells per 952-	10 ³	
10	Chip Date	Client	Average cells per ml. 3443 Milk Dilution	Cells per mL	Avg. cell size (pm)		Chip D	ate Client	Average cells per mL 952- Milk Dilution	10 ³ Cells per ml.	Avg. cell
10 2728	Chip Dee LACTOSCAN 48 Styr 163.2016 r. 11:	Client 141 - WAROSHA21	Average cells per ml. 3443 Milk Olivrion Preserved cow milk 1	3-10 ³ Cells per ml. 554-10 ⁴	Avg. cell size (pm) 12.12		Chip Da A LACTOSCAN 4R Silpr 14	ate Client 5.3.2116 r. 11:12:41 • VARD SHA21	Average cells per ml. 952- Milk Dilation Preserved cow milk 1	Collis per ml. 554-10 ⁹	Avg. cell 12.12
10 2728 2748	Chip Date LACTOSCAN 48 Styr 16.3,2016 r. 11: LACTOSCAN 48 Styr 16.3,2016 r. 11:	Client 141 - VAROSHA21 26 - VAROSHA217	Average cells per ml. 3443 Milk Dilution Preserved cow milk 1	Cells per mL 554-10° 271-10°	Avg. cell size (µn) 12.12 12.46		Chip D 3 LACTOSCAN 4R Shir 10 5 LACTOSCAN 4R Shir 10	ate Client 5.2016 r. 11:12:41 • VURO SHA21 5.32016 r. 11:28:26 • VARO SHA217	Average cells per mL 952- Milk Dilation Preserved cow milk 1 Preserved cow milk 1	10 ³ Cells per ml. 554-10 ⁴ 271-10 ⁴	Avg. cell 12.12 12.45
10 2728 2748 2748 276A	CNP Dee LACTOSCAN 48 Stay 16.3.398 s. 11 LACTOSCAN 48 Stay 16.3.398 s. 11 LACTOSCAN 48 Stay 16.3.398 s. 112	Client 511 - VAROSHA21 26 - VAROSHA217 511 - VAROSHA214	Average cells per ml. 3443 Milk Dilution Preserved cow milk 1 Preserved cow milk 1	Cells per mL 554-10' 271-10' 1141-10'	Avg. cell size (µn) 12.12 12.46 12.36		Chip D B LACTOSCAN 4R Skyr M B LACTOSCAN 4R Skyr M A LACTOSCAN 4R Skyr M	ate Client 5.1.2916 (r. 11:12:41 - VARO SHA21 5.1.2916 (r. 11:21:25 - VARO SHA217 5.1.2916 (r. 11:21:11 - VARO SHA214	Average cells per mL 952* Milk Dilution Preserved cow milk 1 Preserved cow milk 1	Calis per ml. 554-10° 271-10° 1742-59°	Avg. cell 12.12 12.46 12.36
10 2728 2748 2754 2754 2754	Chip Oute LACTOSCAN # Styr #13,3046 ± 111	Client 141 - VAROSHA21 28 - VAROSHA217 11 - VAROSHA217 16 - VAROSHA213	Average cells per mt. 3443 Milk Ditution Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1	Cells per ml. 554-10° 271-10° 1142-10° 918-10°	Avg. cell nize (µm) ^ 12.12 12.46 12.36 12.25 12.25		Chip D B LACTOSCAN 4R Sign 11 S LACTOSCAN 4R Sign 11 A LACTOSCAN 4R Sign 11 A LACTOSCAN 4R Sign 11 A LACTOSCAN 4R Sign 11	ate Client 3.2916 r. 11:12:41 VARO SHA21 5.2916 r. 11:2026 v VARO SHA21 5.2916 r. 12:026 v VARO SHA213 5.32916 r. 11:37:06 v VARO SHA213	Average cells per mt. 952- Milk Dilation Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1	Cells per ml. 554-10° 271-10° 1747-59° 918-10°	Avg. cel. 12.12 12.46 12.36 12.25
10 2728 2748 2748 276A 278A 2784	Chp Dee LACTOSCAN 48 Styr 16.3,095 E 111 LACTOSCAN 48 Styr 16.3,095 E 111	Client 241 - VAROSHA21 28 - VAROSHA217 11 - VAROSHA214 6 - VAROSHA214 22 - VAROSHA214	Average cells per mt 3443 Milk Ditution Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1	3-103 Cells per mL 554-10 271-10 1147-10 918-10 787-10	Avg. cell tize (µm) 12.12 12.46 12.26 12.26 11.98	0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	Chip D 3 LACTOSCAN 4R Skyr M 4 LACTOSCAN 4R Skyr M 4 LACTOSCAN 4R Skyr M 5 LACTOSCAN 4R Skyr M 8 LACTOSCAN 4R Skyr M	ate Client 5.2916 r. 11:12:41 · WARO SHA21 5.2916 r. 11:20:42 · WARO SHA27 5.2916 r. 11:20:42 · WARO SHA27 5.2916 r. 11:27:32 · WARO SHA274 5.2916 r. 11:27:32 · WARO SHA274	Average cells prr ml. 952 - Mitk Dilation Preserved cow mitk 1 Preserved cow mitk 1 Preserved cow mitk 1 Preserved cow mitk 1	Cells per ml. 55410° 27510° 1747-59° 918-10° 787-10°	Arg. cell 12.12 12.46 12.36 12.25 11.90
10 2728 2748 2748 276A 276A 2778 2778 2778 2778 2778 2778 2778 277	Chip Device LACTOSCAN 48 Spy H3.2846 - 111 LACTOSCAN 48 Spy H3.2846 - 111	Client 241 - VAROSHA21 25 - VAROSHA214 26 - VAROSHA214 26 - VAROSHA214 22 - VAROSHA214 22 - VAROSHA218	Average cells per mt 3444: Mik Preserved cow mik 1 Preserved cow mik 1	Cells per mL 554-10° 221-10° 1147-10° 918-10° 787-10° 588-10° 588-10°	Avg. cell size (µm) 12.12 12.46 12.35 12.25 11.38 10.50 11.30		Chip D 8 LACTODICAN 48 SHIP 19 9 LACTODICAN 48 SHIP 19 0 LACTODICAN 48 SHIP 19	ate Client 5.2364 ; 11:241 · WADOHA21 5.2364 ; 11:2265 · WADOHA21 5.2016 ; 11:2265 · WADOHA21 5.2016 ; 11:2167 · WADOHA21 5.2016 ; 11:216 · WADOHA21 5.2016 ; 11:216 · WADOHA21 5.2016 ; 11:216 · WADOHA21 5.2016 ; 11:4627 · WADOHA21	Average cells pr. 952 -	Colls per ml. 55510° 27110° 1747-59° 918-10° 787-10° 5282-59° 508-100	Avg. cell 12.12 12.46 12.36 12.25 11.90 16.50
10 2728 2748 2748 276A 276A 276A 276B 276B 276B 276B 276B 276B 276B 276B	CMp Down LACTOSCAN 49 Styr H3.298 F. H1 LACTOSCAN 49 Styr H3.298 F. H1	Client 41 - VAROSHA21 28 - VAROSHA217 11 - VAROSHA214 46 - VAROSHA214 22 - VAROSHA218 22 - VAROSHA218 45 - VAROSHA218	Average cells per ml. 3445 Witk Dilution Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1 Preserved cow milk 1	Cells per ml. 554-10° 271-10° 1747-10° 174	Avg. cell tize (µm) 12.12 12.46 12.25 11.38 11.38 11.70 11.70		Chip D 3 LACTOSCAN 4R Spir 10 4 LACTOSCAN 4R Spir 10 4 LACTOSCAN 4R Spir 10 5 LACTOSCAN 4R Spir 10 6 LACTOSCAN 4R Spir 10	ate Client 5.32916 / 1112/41 · WARDSHA21 5.32916 / 1126/2 · WARDSHA22 5.32916 / 1126/2 · WARDSHA214 5.32916 / 1127/2 · WARDSHA214 5.32916 / 1127/2 · WARDSHA214 5.32916 / 1152/5 · WARDSHA216	Average cells pg Mith 9522- Preserved cow mitk 1 Preserved cow mitk 1	Calls per ml. 554-10° 271-10° 1441-19° 918-10° 707-10° 5082-19° 309-10°	Arg. cel 12.12 12.46 12.35 12.25 11.97 18.59 11.70
10 2728 2748 2748 276A 276A 2778 2778 2778 2778 2778 2778 2778 277	Chip Deve (A-CTOSCAN # 59gr H3.238* 1.11 (A-CTOSCAN # 59gr H3.238* 1.12)	Client 141 - VAROSHA21 28 - VAROSHA217 21 - VAROSHA214 22 - VAROSHA214 22 - VAROSHA214 22 - VAROSHA214 23 - VAROSHA216 24 - VAROSHA216 16 - VAROSHA216	Average cells per mit. 3444: Preservet cow mill: 1 Preservet cow mill: 1	3-103 Cells per ml. S54 107 1147-107 918-107 787-107 S382-107 309-107 289-107	Avg. cell tize (an) 12.12 12.46 12.25 11.98 10.59 11.76 11.76 12.93		Chip D S LACTOSCAM 4E Sign D LACTOSCAM 4E Sign L LACTOSCAM 4E Sign	nto Client 5.2846 r. 11-5241 - WADDIAL21 5.2846 r. 112242 - WADDIAL21 5.2846 r. 112241 - WADDIAL21 5.2846 r. 112241 - WADDIAL21 5.2846 r. 112712 - WADDIAL21 5.2846 r. 112712 - WADDIAL21 5.2846 r. 112712 - WADDIAL21 5.2846 r. 11261 - WADDIAL21 5.2846 r. 11261 - WADDIAL21	Average cells per 952 - Mit 952 - Preserved cov milk 1 Preserved cov milk 1	103 Calls per ml. 554-109 1747-509 1747-509 1747-509 1747-509 1747-509 1747-509 1747-509 1747-509 1747-509 1747-509 1749-509 1709-509	Arg. cel 12.12 12.46 12.25 11.98 16.50 11.70 12.93
10 2778 2778 2778 2778 2778 2778 2778 277	CLU LACTOSCAN REMOVE HL. JUNC HL LACTOSCAN REMOVE HL LACTOSCAN REMOVE HL. JUNC HL LACTOSCAN REMOVE HL LACTOSCAN	Client 41 - VAROSNA21 28 - VAROSNA21 11 - VAROSNA214 36 - VAROSNA214 32 - VAROSNA214 32 - VAROSNA218 45 - VAROSNA218	Average cells per discussion of the second core mills 1 Preserved core mills 1	3-10 ³ Cells per ml. 554-10° 271-10° 1147-10° 918-10° 787-10° 5820-10° 5820-10° 5820-10° 5820-10°	Avg. cell size (pre) 12.12 12.26 12.25 11.28 11.28 11.28 11.29 1		Chip Di 8 LACTOSCAN 4R Sbp. 19 9 LACTOSCAN 4R Sbp. 10 9 LACTOSCAN 4R Sbp. 10 9 LACTOSCAN 4R Sbp. 10	ate Client 5.3286 (+ 1124/2 + WADDSHA27 5.3286 (+ 1124/2 + WADDSHA27 5.3286 (+ 1126/2 + WADDSHA27 5.3286 (+ 1137/6 + WADDSHA27 5.3286 (+ 1137/6 + WADSHA27 5.3286 (+ 1136/2 + WADSHA27 5.3286 (+ 1146/2 + WADSHA27 5.3286 (+ 1146/2 + WADSHA27) 5.3286 (+ 1146/2 + WADSHA27)	Average cetts per 952 - bits Deserved over milk 1 Preserved over milk 1	Colls per ml. 554-10° 215-10° 1747-59° 918-10° 717-50° 918-10° 717-50° 309-10° 6475-50°	Ang. cel 12.12 12.46 12.25 11.25 11.50 11.70 12.93 11.99
10 2728 2748 2744 275A 2796 2796 2796 2796 2796 2796 2796 2796	Chip Date LCCIDICAL # Sty 11.2374 TH LCTOCAL # Sty 11.2374 TH	Client H1 - VAROSHA21 28 - VAROSHA217 11 - VAROSHA213 36 - VAROSHA213 22 - VAROSHA214 22 - VAROSHA214 23 - VAROSHA216 24 - VAROSHA215 36 - VAROSHA213 36 - VAROSHA213 36 - VAROSHA213	Average cells per 11 3444 Billit Deserved corr mills 1 Preserved corr mills 1	3-103 Cells per ml, 554-107 1747-107 1747-107 538-107 339-107 289-107 289-107 289-107 249-107 542-107 542-107	Aug. cell size (sm) 12.12 12.26 12.26 12.26 12.25 12.25 10.59 10.59 10.79 12.20 11.74 12.20 11.74 12.75 12		CMp D 3 LACTOSCAM 48 Pbp 41 4 LACTOSCAM 48 Pbp 41 4 LACTOSCAM 48 Pbp 41 4 LACTOSCAM 48 Pbp 41 5 LACTOSCAM 48 Pbp 41 6 LACTOSCAM 48 Pbp 41	ata Client 2014 - Cl	Average calls are by 52.	Colls per ml. 554-10° 271-10° 1747-10° 1747-10° 309-10° 209-10° 209-10° 209-10°	Avg. cel 12.12 12.46 12.36 12.25 11.98 10.59 11.70 12.93 10.99
10 9 2728 9 2948 9 278A 9 278A 9 278A 9 278A 9 278A 9 282A 9 284A 9 284A	Chip Date LaClostical et 80, 11,3141, 11,11 11,1141,114,114 LaClostical et 80, 11,2341,117 11,114,114,114,114 LaClostical et 80, 11,2341,117 11,114,114,114,114,114,114,114,114,114,	Client 41 - WAROSHA21 22 - WAROSHA21 23 - WAROSHA214 45 - WAROSHA214 23 - WAROSHA214 23 - WAROSHA214 24 - WAROSHA214 25 - WAROSHA214 26 - WAROSHA213 28 - WHROSHA213 28 - WAROSHA213 28 - WAROSHA213 28 - WAROSHA213 28 - WAROSHA213 28 - WAROSHA213 28 - WAROSHA213 29 - WAROSHA213 29 - WAROSHA213 20 - WAROSHA213	Average cells per difference of the second core mills 1 Preserved core mills 1	3-103 Cells per ml. 554-107 211-107 1147-107 918-107 787-107 5582-107 368-107 5582-107 5582-107 552-107 247-107 522-107 247-107	Asg. cell tize (pri) 12.12 12.26 12.25 12.25 12.25 11.28 10.59 11.76 12.93 10.99 11.76 1.55 11.76		CL0 CL0 B LACTOSCAM 48 Stop 11 S LACTOSCAM 48 Stop 11 A LACTOSCAM 48 Stop 11 A LACTOSCAM 48 Stop 11 A LACTOSCAM 48 Stop 11 S LACTOSCAM 48 Stop 11 B LACTOSCAM 48 Stop 11 CACTOSCAM 48 Stop 11 LACTOSCAM 48 Stop 11 LACTOSCAM 48 Stop 11 LACTOSCAM 48 Stop 11 LACTOSCAM 48 Stop 11 LACTOSCAM 48 Stop 11	ata Client 3.394 - 11224 - 140205427 3.394 - 11224 - 14005427 3.394 - 11224 - 14005427 3.394 - 11224 - 14005427 3.394 - 11226 - 14005427 3.394 - 11266 - 14005427 3.394 - 11266 - 14005427 3.394 - 12666 - 14005427 3.394 - 12667 - 14005427 3.394 - 1400547 - 14005427 3.394 - 14005427 3.394 - 1400547 - 14005427 3.394 - 1400547 - 14005427 3.394 - 14005427 3.395 - 14005427	Average cells are 952-	Colls per ml. 554-10° 215-10° 1142-10° 918-10° 717-10° 215-10° 2142-10° 209-10° 6475-10° 6475-10° 6475-10° 5247-10°	Arg. cell 12.12 12.46 12.36 12.25 11.96 11.70 12.99 11.76 13.99 11.76 9.35
10 2728 2748 2758 2768 2768 2708 2708 2708 2708 2708 2708 2708 270	Cop Desc LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111 LaCTORCAL #109 bits 123.595 + 111	Client 41 - WAROSHA21 21 - VAROSHA21 21 - VAROSHA214 46 - VAROSHA214 46 - VAROSHA214 22 - VAROSHA214 22 - VAROSHA216 23 - VAROSHA212 24 - VAROSHA212 26 - VAROSHA212 26 - VAROSHA212 27 - VTellmpeM022 22 - VTellmpeM022 22 - VTellmpeM022 23 - VTellmpeM022 24 - VTellmpeM022 24 - VTellmpeM022 25 - VTellmpeM022 25 - VTellmpeM022 26 - VTellmpeM022 27 - VTellmpeM022 27 - VTellmpeM022 28 - VTellmpeM022 29 - VTellmpeM022 20 - VTellmpeM02 20 - VT	Assesse only 3444 Mile Obtion Preserved corr mill 1 Preserved corr mill 1	3-103 Cells per ml. 554-107 271-107 144-107 5582-107 5582-107 306-107 306-107 306-107 522-107 222-107 222-107 222-107 222-107 222-107 222-107 221-107	Aug. cell tite (and) 12.12 12.26 12.25 12.25 10.59 10.59 10.59 11.76 8.75 8.71 1.71 1.75 8.71 1.71 1.75 1.71 1.75 1		Chip D B LACTOICAM #R Step 19 LACTOICAM #R Step 19 LACTOICAM #R Step 19 LACTOICAM #R Step 10 LACTOICAM #R Step 10	ate Client 1.30% (11524) WAROSHA21 3.0% (11524) WAROSHA21 4.0% (11524) WARO	Average calls are going of the second	103 Colls per mt. S54-10° 221-10° 1447-19° 548-10° 548-10° 209-10° 6475-10° 6475-10° 6475-10° 6475-10° 6475-10°	Arg. cell 12.12 12.46 12.36 12.25 11.94 15.50 11.70 12.93 10.99 11.76 9.35 9.31

2. Go to the image and press the right button of the mouse on it.

1.00	under Cont	one Advanced			47 4 2047	- 46.27		Main	Ch Graph Ontions Ad	vanced			47 4 2017
1 101	ipn opr	Auvanceu			17.1.2017	1. 10.27 4.		15 main		Farren			17.1.2017
						(O)					<u> </u>		
				Start				1				Start	
				Start				. <u>18</u> .08				otant /	
								1.4	Save As				
								- 14 - ¹⁷		1 de			
			D										
				10 282C								0 282C	
			Clie	err VAROSHA213						190 - P	Clier	VAROSHA213	
			Milk typ	pe Preserved cow mill			//				Milk typ	 Preserved cow milk 	
			Average cells p	52·10 ³				1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 - 1999 -			Average cells pr	952·10 ³	
	0	Clive		01.0 C.B.					China Data China			Walan Calle second	Anno collision for
CTOSCAN 4R 50ur 1	6.3.	2015 r. 11:12:41 • VAROSHA21	Preserved cow milk	1 554-10 ⁴	12.12			2728	LACTOSCAN 4R 50ur 16.3.2016 r. 11:12:41 • VARO	ISHA21 P	Preserved cow milk 1	554-10 ^o	12.12
ICTOSCAN 4R 50pr 16	3.	2016 r. 11:20:26 • VAROSHA217	Preserved cow milk	1 271-104	12.46			2748	LACTO SCAN 4R 50pr 16.3.2016 r. 11:20:26 • VARO	ISHA217 P	Preserved cow milk 1	271-10*	12.46
CTOSCAN 4R 50pr	16.3.	2016 r. 11:29:11 VAROSHA214	Preserved cow milk	1 1747-10	12.36	La		276A	LACTOSCAN 48 50gr 16.3.2016 r. 11-20-11 • VARO	SHA214 P	Preserved cow milk 1	1747-10*	12.36
ICTOSCAN 4R SOJA	16.3	2016 r. 11:37:06 • VARUSHALT3	Preserved cow milk	1 910-10	12.25	Ĕ		2738	LACTOSCAN 48 50pr 15.3.2016 / 11.37.06 / VARO	SHA213 P	Preserved cow milk 1	787-10*	12.25
CTOSCAN 4R 50pr 16	Ĵ.	2016 r. 11:46:27 - VAROSHA218	Preserved cow milk	1 5282-10*	10.50	15		2798	LACTOSCAN 4R 50pr 16.3.2016 r. 11.46-27 · VARO	SHA218 P	Preserved cow milk 1	5282-10*	10.50
ICTOSCAN 4R 50pr 16	3.	2016 r. 11:58:45 - VAROSHA210	Preserved cow milk	1 309-10*	11.70	(O)		281D	LACTOSCAN 4R 50pr 16.3.2016 r. 11:58:45 • VARO	ISHA210 P	Preserved cow milk 1	309-10 ⁴	11.70
ICTOSCAN 4R 50pr 16.3.		2016 r. 12:00:16 · VAROSHA211	Preserved cow milk	1 289-10*	12.93			282A	LACTO SCAN 4R 50pr 16.3.2016 r. 12:00:16 · VARO	SHA211 P	Preserved cow milk 1	289-10*	12.93
ICTOSCAN 4R 50µr (16.3.20	30	16 r. 12:00:43 • VAROSHA212	Preserved cow milk	1 6475-10	10.59			282B	LACTOSCAN 4K SQF 15.3.2016 r. 12:00:43 · VARO	PSHAZIZ P	Preserved cow milk 1	6475-10*	10.39
CTOSCAN 48 5hr 17.3.2016 r. 11	2016 c. 11	E1/28 · Vladimosk052	Preserved now milk	1 247-107	9.95			5 DEL 2010	LACTOSCAN 48 Shar 17.3.2015 r. 11.03.28 v Viada	mosk057 P	Preserved cow milk 1	207-10	9.95
CTOSCAN 4R 50pr 17.3.2	2	016 r. 11:07:22 • Vladimpek055	Preserved cow milk	1 62-10	9.21	i i i i i i i i i i i i i i i i i i i		5 304C	LACTOSCAN 4R 59µr 17.3.2016 r. 11.07.22 · Viada	mpek055 P	Preserved cow milk 1	62-10*	9.21
CTOSCAN 4R 50pr 17.3.20	20	16 r. 11:23:50 · Vladimpek055	Preserved cow milk	1 84-10*	11.49			\$309C	LACTOSCAN 4R 50pr 17.3.2016 r. 11.23:50 • Vladi	mpek055 P	Preserved cow milk 1	84-10*	11.49
ICTOSCAN 4L 50µr 29.2		2015 r. 09:40:09 + K.Rusev2	Preserved cow milk	1 3443-10*	11.68			a 322D	LACTOSCAN 4L 50pt 29.3.2016 r. 09:40:09 · K.Rus	IEV2 P	Preserved cow milk 1	3443-19*	11.68
			<u></u>			17 N 2817 E NB 27 N	an	e e Ladracan con					
		b		Save As		17 1207 : W274		ww.Ladoargan.com					
5	s the	butto	n 🗌	Save As		17 1 2017 x 1627 u		wwiaddacae.com					
59	s the	butto	n 🗌	Save As		12307/4224		en Listensen ese					
1	s the	butto	n	Save As	17.1.2017	et 1300 / 4620 v		e Sare inoge file	Vicus Colors and				17.1.201
S	the	ons	n	Save As	17.1.2017	11:2072-10214		Save image fite	tertanti e KCleopticeje e kopi	• Hellower	na knojež – P		17.1.201
S	the	butto	n	Save As	17.1.2017	асанската (г. 16:28 ч.	R	Sove image file Copressionage in EN + C	transmits a KC bageloops a long 2	a a	al logit2 P Si + 0		17.1.201
S	the	ons (Advanced	n	Save As	17.1.2017	r. 16:28 ч.	F	See image file Cyrrenessan - Die Baatssimon - D	type marks a % (* Hogeliforgis a Hoge) Taxa mark	• •	2010 21 - 0		17.1.201
5 1	the	butto	n	Save As	17.1.2017	r. 16:28 ч.		See inage file	ternole + XC tegeldoge + keg2	n= 42 Tuotera	्र स्वार्थ स्वार स्व	Start	17.1.201
st	he	butto	n	Save As	17.1.2017	7 : 15:28 ч.		Sole image file Sole image file Cyransogaet * N Stappenson X Styress N Styress N Styress	tter mette a ScChapeldinge a Image Nati ante	in 17 Marca	e sta	Start	17.1.201
s t	he	butto	n	Save As	17.1.2017	r. 16:28 ч.		Silve image file Cyranogenet Bagesparsen Kopynen Kopynen Kopynen Kopynen Miganog	Yana manta a SC Yangadi Singka a Hanga 2 Inte a mata 1	io 17 Manua	na langa 2 – p sa langa 2 – p sa langa 2 – p	Start	17.1.201
s tł	o opt	butto	n	Save As	17.1.2017	7 : 16:28 4.		Size image for	Ter same a kChopelogie a kopst na mor 1	in 2 Datase	21 • Ø	Start	17.1.201
s th	e opt	butto	n	Save As	17.1.2017	Г. 16:28 ч. (7. 16:28 ч.) (2)		Since image file Parameters	type names a SC Charged Single & Honge 2 to at same 1	namen al anti-	er hergel – P Si + O	Start	17.1.2012
s the	Ipt	butto	n	Save As	17.1.2017	7 : 16:28 ч.		Storinge for Participation of the second Participation of the second Pa	Yzeronata w KCongologija w Ingez Dan sava 1		n top? Ri * Ø	Start	17.1.2017
s tł	ne • opt	butto	n	Save As	17.1.2017	стано, цато (г. 16:28 ч. (б) (Д)		Since Hange for U Dial Consensation Conse	tter mette + SC Stepphöret + Imp/J htt anna 1	10 22 Matter	at map) Ci + Q	Start	17.1.201
s th	Opt	butto	n	Save As	17.1.2017	7 - 15:28 ч.		Sterings for Participation Partici	trat montes & K2 (transpillinguis & longe) in a survi I		a top2 P	Start	17.1.2012
s ti	ne	butto	n	Save As	17.1.2017	r. 16:28 ч.		See Hage for U Dia Constraint Constraints	THE South & KCCoopeRings's & Soup J and The South Sout	10 27 Macana	alwy) R + 0	Start Paccasta	17.1.2012
s the		butto	n	Save As	17.1.2017	7.15:28 4. (2)		Size stage for e to a stage for e to a stage for e to a stage for definition of the stage for defi	taranata a (Chaptiling) a Inagri Ina Inat I	namenal na 22 Danner	a topz - P k) + Q	Start 282C Wreasyvel cow milk	17.1.2017
s the	ic	butto	n Cicilia type	Save As Start	17.1.2017	сант, цтэ г. 16:28 ч. (2) (2) (2) (2) (2) (2) (2) (2) (2) (2)		See ange for Decemped of the second	Ter same a kCoopelogie a loop J ne neu		21+ 0	Start 24000HA213 Wreserved cow milk	17.1.2017
s th	e opt	butto		Save As	17.1.2017	 7: 16:20 ч. (b) (c) (Store image for Store	trans marks + 1 (Changel Bings + Honge 1 Inst marks 1	in 2 Manual	20 mg 27 p 21 + 0	Start 7450 GMA213 Preserved cow milk 952-10°	17.1.2017
s th	ne	ons Advanced	n	Start	47.1.2017	7. 1623 4. 8.		Stars range for = 80 + 80 + 80 + 80 + 100	tree mante a XCCompetitionale e lange 2 en anna 1 dess Billing III. con on Reconst P. Jangi	• 22 Dames		Start Preserved cow milk 952-10 ³	17.1.2017
s th	Deter	Advanced		Save As Start	• 17.1 2017	(1. 16/28 vs. (1. 16/28 vs.	⇒	Steerings for Background and Steerings for Background and Steerings for Steerings for Steering for Steerings for Steerings for Steering for Steering for Steering for Steering for Steerings	trae name + 20 (Hangeldingels + Hangel Han Hangel 1 dam Blang II, om or Ma Kana ("Bagi	• 2 Tomes		Start Z82C VAROBHA213 Preserved cow milk 952-103 zmr	17.1.2017
s th				Save As Start	• 17.1.2017	r, 19:20 ч. (6) (7) (7) (7) (7) (7) (7) (7) (7) (7) (7	⇒	Schemagner S	Anna Sheng to gene of the conte (" Seg)	20 Danner 20 Danner 20 Danner		Start VAROSHA213 Preserved com milk 952-103	47.1.2017
s th	Deter Styr 16.1. Styr 16.1.		n R R R R R R R R R R R R R R R R R R R	Save As Start	• 17.1.2017	(f. 1823 4. (f) (f) (f) (f) (f) (f) (f) (f) (f) (f)		State stage for United States and States States and States States and States State	The create a SC Supplicing a Page 7 Total create 1 Medical Biology II. Jacobson (* Marg) UCCTOSCAN III Sign (* 1.3.291 + 111-22 * 1.4000) UCCTOSCAN III Sign (* 1.3.291 + 111-22 * 1.4000) UCCTOSCAN III Sign (* 1.3.291 + 111-22 * 1.4000)	• 27 Transv 2000		Start 282C VAROSHA213 Preserved cow milk 952-10 ² Immediate State 254-10 ²	17.1.2017
en e	001	butto	N N N N N N N N N N N N N N N N N N N	Save As	17.1.2017 4 4 4 4 4 4 4 4 4 4 4 4	r. 1623 4.		Concerning for	Transmiss & SC Standardings & Jungs 2 Tax mass 1 1 International Standard Standa	• 22 17200000 	restanting Presented care milli	Start	17.1.2017
s the b ph Copion Copie de Copion Copie de Copie de Copie de Copie Copie de Copie de Copie de Copie de Copie Copie de Copie de Co	0015 - 111 2016 - 111 2016 - 111 2016 - 111 2016 - 111		n	Save As Start ** 2220 *** 952-10 ² *** 952-10 ² *** 952-10 ² *** 952-10 ²	* 17.1.2017	(7. 16:23 4. (7. 16:23 4. (8) (9) (9) (9) (10		the mage for	The second as a SCC impediately a stage 1 and an analysis of the second stage 1 and	• 22 Talana - 22 Talana - 23 - 24 - 24		Start Preserved cow milk 952-10 ³	47.4.2017
	pt			Save As Start 282C VAROBHAZ13 Preserved contained VAROBHAZ13 Preserved contained Start Chart Start Chart Start	* 17.1 2017	Image: Control of the second	⇒	The same type of the sa	The encoder as a Competitive as a surged a	• 22 Promotion SNU27 P SNU27 P SNU27 P		ZB2C VAROBHA213 Mreserved cov milk S52-003 Mre Citie gent U299 U299 With Citie gent U299 U299 With Citie gent	17.1.2017
s t	Control 1000 (1999) (19			Save As Start Start Start Preserved cow mile 952-10 ² There 952-10 ² Start 952-10 ²	* 17.1.2017 17.1.2017	r. 19:20 ч.		Care image for	And Sheng is one of the control (Seg)	• 12 Talana SSU277 P SSU278 P SSU278 P SSU278 P SSU278 P SSU278 P		Start VAROBHA213 Preserved com mik 952-10 ²	17.1.2017
S ph			n Caracteristication Car	Save As Start	• 17.1.2017 17.1.201	7. 19:28 4.		The angle is a set of the second set of the seco	Not create > 5/C Supplicing is + Image? Not create 1 Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image: Image is create of the create? Image is create of the create? Image is create of the create? Image is create of the create? Image is create of the create? Image is create of the create? Image is create? Image is create? <td></td> <td></td> <td>Start 282C VARDSHA213 Preserved cow milk 952-10° Inter Start</td> <td>17.1.20</td>			Start 282C VARDSHA213 Preserved cow milk 952-10° Inter Start	17.1.20

LACTOSCAN 48 59µr 17.32915 r. 11.03.28 v Vladimpek052 Preserved cow milk 1 LACTOSCAN 48 59µr 17.32915 r. 11.97.22 v Vladimpek055 Preserved cow milk 1 LACTOSCAN 48 59µr 17.32915 r. 11.97.25 v Vladimpek055 Preserved cow milk 1 LACTOSCAN 48 59µr 17.32915 r. 11.92.59 v Vladimpek055 Preserved cow milk 1 LACTOSCAN 48 59µr 17.32915 r. 11.92.59 v Vladimpek055 Preserved cow milk 1

LACTOSCAN 4R 59µ 17.3.2015 r. 11:02.28 · Vladimpek052 Preserved cow milk 1 LACTOSCAN 4R 59µ 17.3.2016 r. 11:07.22 · Vladimpek055 Preserved cow milk 1 LACTOSCAN 4R 59µ 17.3.2016 r. 11:07.219 · Vladimpek055 Preserved cow milk 1 LACTOSCAN 4R 59µ 17.3.2016 r. 11:07.219 · Vladimpek055 Preserved cow milk 1 247-10⁴ 62-10⁴ 84-10⁴ 9.95 9.21 11.49 11.68

5 300 304C 309C 5 322D

247-10⁴ 62-10⁴ 84-10⁴

9.95 9.21 11.49 3. Select the folder where you want to save the image, write the name of the image, select the .BMP file format from the dropdown menu and press button SAVE to finally save the image as .BMP file.



NOTE! The image can be saved in .TIFF file format as well. The .TIFF file is a RAW image format and the saved images can be processed in special software without a losing the quality of the initial taken image.

Before saving the image, the User can select a visualization mode such us Software processed image, Marked cells image, Raw image and to save in .BMP file the image as it is seen in the software.

Examples:



Setting the analysis' parameters

The Menu Options allows the user to set the analysis parameters for different types of milk.



1. Setting the minimum and maximum size of the cells:

1.1. Choosing the sample of milk, for which the set parameters have to be applied.

With the assistance of the drop-down menu, choose the sample of milk for which change of the pre-set parameters is needed.



1.2. Filling in the required parameters.

When the milk sample is chosen, pre-set values of the minimum and maximum size of the cells, which will be counted appear below it. With up and down arrows or by manual entry of value, the range of the counted cells may be changed according to the certain type of milk.

1.3. Setting the limits of the clinical and subclinical mastitis.

According to the breed of the animal, the local and state requirements in the different countries, the limits of the clinical and subclinical mastitis may differ.

For example in Germany it is accepted that the healthy animals are these, which milk contains somatic cells up to 100 000 pcs/ml. Over 100 000 pcs/ml up to 400 000 pcs/ml the animal has a subclinical mastitis. And over 400 000 pcs/ml the animals has a clinical mastitis.

By coloring the analysis result in green, orange and red, the device signalizes for the health condition of the animals.



Setting the limits of the subclinical mastitis is made by changing the value in the field "Warning level of cells count". Changing of the value is made with up and down arrows or by setting the value manually.

Setting the limits of the clinical mastitis is made by changing the value in the field "Critical level of cells count". Changing of the value is made with up and down arrows or by setting the value manually.



1.4. Saving the entered settings

By press the button APPLY Save changes, all entered settings are saved in the program. After being saved, the button changes its color from green to grey.



2. Setting the Laboratory offset in percentages:

It is accepted that there is a difference of the results from analysis of one and same milk sample between two laboratories. According to ISO 13366-1:2008 "Milk – Enumeration of somatic cells – Part 1: Microscopic method (Reference method) », Annex A Collaborative trail, Table A. 1 – Results from interlaboratory test, the deviation between two laboratories can be up to 20%. If LACTOSCAN SCC is checked by using reference samples from local laboratory, a difference in the results between the reference samples value and the results from LACTOSCAN SCC can be received. In this case, a correction using the option Laboratory offset can be done by filling in a number indicating the percentages of correction.

The Laboratory offset number can be:

0 (zero) – meaning that no correction is done.

+N (positive number) – meaning that positive correction in percentages % is done. The final result of the analysis is increased with the set percentages of correction.

N (negative number) – meaning that negative correction in percentages % is done.
 The final result of the analysis is decreased with the set percentages of correction.



Attention!

To determine the need of correction and to do a correction, the instructions below must be strictly followed:

- For determining the deviation, to be used **only reference milk samples from IDF accredited laboratory**.

- To be used minimum 3 different reference samples with the following ranges of cells concentration:

 1^{st} Reference sample cell concentration range: **50 000 – 200 000 cell/mL** 2^{nd} Reference sample cell concentration range: **400 000 – 600 000 cell/mL** 3^{rd} Reference sample cell concentration range: **800 000 – 1 200 000 cell/mL**

- Reference samples must be measured in **Reference mode** of LACTOSCAN SCC following the procedure for Reference sample testing.

- The received results from the analyses to be filled in the Excel table for calculation of Laboratory offset available at:

http://lactoscan.com/editor/ufo/files/Laboratory_offset_table.xlsx

Or to be written in a table shown below:

Raw data	Sample 50000- 200000 cell/mL 1	Sample 400000- 600000 cell/mL 2	Sample 800000- 1400000 cell/mL 3
Measurement 1			
Measurement 2			
Measurement 3			
Measurement 4			
Measurement 5			
Measurement 6			
Measurement 7			
Measurement 8			
Measurement 9			
Measurement 10			
Mean value			
Reference sample			
value			
Difference in			
percentages			
Average deviation			

- The Mean values (Average results) for each sample to be calculated using the formula below:

$MV = \frac{M1 + M2 + M3 + M4 + M5 + M6 + M7 + M8 + M9 + M10}{10}$

Where:

 $\mathsf{MV}-\mathsf{Mean}$ value (Average result) from the results of the 10 measurements of one and same reference sample

M1,M2,M3,M4,M5,M6,M7,M8,M9,M10 – Results from the ten measurements of the one and same reference sample

- The received Mean values (Average results) from the calculation of all three samples to be compared with the reference sample values from the labels of the analyzed reference samples. The difference in percentages between each the Mean value and the Reference sample value to be calculated using the formula below:

Where:

D1% - the difference in percentages for 1^{st} reference sample RSV1 – 1^{st} Reference sample value from the label MV1 – Mean value for 1^{st} sample

- The difference in the received results for the three samples must be only positive (+) value or only negative (-) value for all the three of them. If you receive positive (+) and negative (-) values, this means that there is an inaccuracy in the measurement, according to the value of the reference sample, written on the label of the sample. In such cases <u>no corrections</u> should be done, until the correctness of the laboratory samples is proven.

- If the results from the three samples are all positive (+) or all negative (-), then the deviation in percentages has to be calculated. This is done using the formula below:

$$D\% = \frac{D1\% + D2\% + D3\%}{3}$$

Where:

D% - the average deviation of the three samples

D1% - the difference in the received results for the first sample in percentages

 $\mathsf{D2\%}$ - the difference in the received results for the second sample in percentages

D3% - the difference in the received results for the third sample in percentages

After calculating the average deviation for the three samples, the received number can be filled in the field for Laboratory offset.

Sample table with results from calculation of Laboratory offset:

Raw data	Sample 50000- 200000 cell/mL	Sample 400000- 600000 cell/mL 2	Sample 800000- 1400000 cell/mL 3
Measurement 1	142 000	571.000	1 153 000
Measurement 2	161 000	563,000	1 164 000
Measurement 3	147 000	548 000	1 146 000
Measurement 4	157 000	530,000	1 142 000
Measurement 5	142 000	543 000	1 146 000
Measurement 6	145 000	584 000	1 122 000
Measurement 7	151 000	558 000	1 100 000
Measurement 8	150 000	584 000	1 154 000
Measurement 9	140 000	539 000	1 120 000
Measurement 10	145 000	570 000	1 155 000
Mean value	148000	559000	1140200
Reference sample			
value	156000	565000	1230000
Difference in			
percentages	5,13	1,06	7,30
Average deviation		4,50	

Excel formulas for calculating:

Mean Value: =AVERAGE(Measurement 1: Measurement 10) Reference sample value: fill in the value from the label of the reference sample Difference in percentage: =100*(Reference sample value-Mean value)/Reference sample value

Average deviation: =(Difference in percentages for sample 50000-200000 cell/mL + Difference in percentages for sample 400000-600000 cell/mL + Difference in percentages for sample 800000 -1400000 cell/mL)/3

To do this:

2.1. From Main menu go to Options menu:



2.2. Enter the calculated deviation number.

Note! If you have received 5% deviation, you should enter number 5. If you have received -5% deviation, you should enter number -5.



To restore the factory settings, number 0 must be filled in the Laboratory offset field.

Deleting results from the data base

The User can delete results from the data base table. To do this:

1. Select the result that you want to delete and press the right button of the mouse to appear button Delete.

				Start		Mairi	Graph Options	Advanced		Sta	rt) (23.02.2018
• 5			Average cells	per 72-10 ³					Average cell	mL 7377.1	10 ³	
10	Chip Date 1	Client	Average cells	pir 72.10 ³ Dilution Cells per mL	Avg. cell size [sm]	ID OI	Chip Date	Client	Average cell	mL 7377-1	10 ³ Cells per ml.	Arg. cell size [arr]
ID 280	CMp Date 1 LACTOSCAN 48 Spyr 19.02.2016 11-48-45 v	Client Varosha101	Average cells	per mL 72.10 ³ Dilution Cells per mL 1 1 433410°	Avg. cell size [pm] 10.22	10 29C	Chip Date LACTOSCAN 4R 6kp 19.02.2016 11:53-43	Client Varosha103	Average cell Milk Preserved cow mill	Dilution	Cells per ml. 184-10 ⁴	Avg. cell size [are] 10.47
10 10 280 290	CNp Date LACTOSCAN 4R 50p 15/82.2016 11.4046 v LACTOSCAN 4R 50p 15/82.2016 11.53.04 o	Client Varosha101 Varosha103	Average cells y Milk Preserved cow milk Preserved cow milk	Per 72-10 ³ Dilution Cells per mL 1 4334.10' 1 164.10'	Avg, cell size [pm] 10.22 10.47	10 29G 31C	Chip Date LACTOSCAN 4R Rige 19.82.2016 11:53-43 LACTOSCAN 4R Rige 19.82.2016 12:16:24	Client Varosha193	Average cell Milk Preserved cow milk Preserved cow milk	Dilution	Calls per ml. 184-10 ² 801-10 ²	Avg. cell size [arr] 10.47 11.12
10 2 20C 2 25C 3 35C	Chip Date LACTOSCAN 48 Styr 15/82.2016 11.48-64 v LACTOSCAN 48 Styr 15/82.2016 11.53-03 v LACTOSCAN 48 Styr 15/82.2016 12:03-20	Client Varosha101 Varosha103 Varosha103	Average cells y Milk Preserved cow milk Preserved cow milk	Per 72-103 Difution Cells per mL 1 4334-10° 1 184-10° 1 801-10°	Avg, cell size (pm) 10.22 10.47 11.12	ID 29C 31C 328	Chip Date (ACTOSCAN 48 Skyr 19.02.2016 116.53.43) (ACTOSCAN 48 Skyr 19.02.2016 126.526.52) (ACTOSCAN 48 Skyr 19.02.2016 52.105.2	Client Varosha103 Varosha101 Varosha104	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilution	Calls per ml. 184-10 ⁹ 001-10 ⁹ 355-10 ⁹	Arg. cell size [µrs] 10.47 11.12 12.04
10 9 280 290 310 5 328	CHp Dete LACTOSCAR 48 by 19822016 11.48/4 v LACTOSCAR 48 by 19822016 11.53/0 v LACTOSCAR 48 by 19822016 11.63/0 v LACTOSCAR 48 by 19822016 11.61/0 v	Client Vanosha101 Vanosha103 Vanosha101 Vanosha101	Average cells Milk Preserved cow milk Preserved cow milk Preserved cow milk	Per 72-103 Ditution Cells per mt. 1 133-10° 1 184-10° 1 881-10° 1 335-10°	Avg. cell size (pre) 10.22 10.47 11.12 12.04	ID 29C 31C 328 34A	Chip Date LACTOSCAN 48 Sky 19.2.2016 11.53.43 LACTOSCAN 48 Sky 19.2.2016 12.95.20 LACTOSCAN 48 Sky 19.2.2016 12.210.52 LACTOSCAN 48 Sky 19.21.2015 12.10.52	Client Varosha193 Varosha191 v Varosha194 v JMI OOD19	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilution	Calls per ml. 184-10° 001-10° 355-10° 390-10°	Arg. cell size [pre] 10.47 11.12 12.04 9.15
10 29C 29C 31C 328 31A	CMp LACTOSCAR 48 Skp 16 02 2166 11 AU-04 v LACTOSCAR 48 Skp 16 02 2166 11 AU-04 v LACTOSCAR 48 Skp 10 22 2166 12 AU-04 LACTOSCAR 48 Skp 10 22 2016 12 AU-04 LACTOSCAR 48 Skp 10 22 2016 12 AU-04	Client Vansha101 Vansha103 Vansha104 JAN OOD9	Average cells Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Per 72-103 Dilution Cells per mL 1 433410° 1 88410° 1 88410° 1 335410° 1 335410°	Avg. cell size (µn) 10.22 10.47 11.12 12.04 9.15	10 29C 31C 328 344 210	Chip LACTOSCAN 48 Styr: 19.02.015 11:53:43 LACTOSCAN 48 Styr: 19.02.2016 12:552 LACTOSCAN 48 Styr: 19.02.2016 12:552 LACTOSCAN 48 Styr: 10.23.2016 12:5529 LACTOSCAN 48 Styr: 10.23.3016 10:5529	Client Varosha183 Varosha181 Varosha181 Varosha184 Varosha184 VarOsh9184 JMI OC0292	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilution 1 1 1 1 1 1 1 1 1	Calls per ml. 184-10° 001-10° 355-10° 290-10° 2023-10°	Arg. cell size [rm] 10.47 11.12 12.04 9.15 12.43
10 24C 25C 31C 328 34A 338	Chip Date LACTOSCAN el Sigu 19/02/2016 11:34:04 v V LACTOSCAN el Sigu 19/02/2016 11:35:04 v V LACTOSCAN el Sigu 12/02/2016 11:35:04 v V	Client Varosha101 Varosha101 Varosha101 Varosha101 JIN OOD93 JIN OOD92	Average cells Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Per 72-103 Ditution Cells per mt. 1 4354 10° 1 184-10° 1 881-10° 1 336-10° 1 336-10° 1 3562 10°	Avg. cell size [jm] 10.22 10.47 11.12 12.04 0.15 12.43	ID 29C 31C 328 34A 328 428	Chip Date LACTOICAN 46 Styr 1102.2015 11:53.47 LACTOICAN 46 Styr 1102.2015 11:53.47 LACTOICAN 46 Styr (2013)1111253 LACTOICAN 46 Styr (2013)1111253 LACTOICAN 46 Styr (2013)1111253	Client Varosha193 Varosha194 Varosha194 JM OCD99 JM OCD92 JM OCD924	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilution 1 1 1 1 1	Calls par rd. 18410° 00110° 35510° 39010° 5925-19° 12-10°	Arg. cell size [are] 10.47 11.12 12.04 9.15 12.42 11.35
10 24C 29C 37C 328 336 338 233	Chip LaCTOCCAN et Sign 1942 2016 1146-014 1 LACTOCCAN et Sign 1942 2016 1120-01 LACTOCCAN et Sign 1942 2016 1120-01 LACTOCCAN et Sign 1942 2016 1116-010 LACTOCCAN et Sign 1243 2016 1116-01 LACTOCCAN et Sign 1243 2016 1116-01 LACTOCCAN et Sign 1243 2016 1116-01	Client Varosha101 Varosha103 Varosha104 JIM COD09 JIM COD092 JIM COD092	Average cells p Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Pert 72:103 Ditution Cells per m. 1 433-19 1 18-107 1 35-107 1 350-107 1 550-19 1 550-19	Avg. cell size [sm] 10.22 10.47 11.12 12.04 9.15 12.43 11.3	10 29C 31C 328 34A 380 428 43	Che Date LACTOSCAN 48 Styr: 19.23/05 1953/40 LACTOSCAN 48 Styr: 19.23/01 1955/40 LACTOSCAN 48 Styr: 19.23/01 1915/2 LACTOSCAN 48 Styr: (20.33/04 1915/27 LACTOSCAN 48 Styr: (20.33/04 1915/27) LACTOSCAN 48 Styr: (20.33/04 1915/27) LACTOSCAN 48 Styr: (20.33/04 1915/27)	Client • Varcoha113 • Varcoha114 • Varcoha114 • JM OC029 • JM OC0294 • JM OC0294 • JM OC0294	Average cell Mile Preserved cow mills Preserved cow mills Preserved cow mills Preserved cow mills Preserved cow mills Preserved cow mills	Dilution 1 1 1 1 1	Calls per ml. 184-10° 801-10° 335-10° 300-10° 5923-10° 72-10° 727-10°	Avg. cell size [im] 10.47 11.12 12.04 9.15 12.43 11.33 11.5
10 28C 29C 31C 328 31C 328 310 310 310 310 314 316 314 316 316 316 316 316 316 316 316 316 316	Chip Date LACTOSCAN 44 Skyr (1982-2016 11-48-44 v) LACTOSCAN 44 Skyr (1982-2016 11-38-44 v) LACTOSCAN 45 Skyr (1982-2016 11-38-34 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 46 Skyr (1982-314 v) LACTOSCAN 47 Skyr (1982-11 v) LACTOSCAN 46 Skyr (1992-314 v) LACTOSCAN 47 Skyr (1992-11 v) LACTOSCAN 47 Skyr (1992-11 v) LACTOSCAN 47 Skyr (1992-11 v) LACTOSCAN 47 Skyr (1992-11 v) LACTOSCAN 47 Skyr (1992-11 v)	Client Vanosha101 Vanosha103 Vanosha104 JIM OOD92 JIM OOD92 JIM OOD92 JIM OOD92	Average cells Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Per 72-103 Ditution Cells per mL 1 4334 10 1 164-10 1 555 10 1 355 10 1 355 10 1 552 10 1 727 10 1 727 10	Aug. cell size (jm) 10.22 10.47 11.52 12.04 9.55 12.43 11.57	ID 29C 31C 328 34A 236 428 430 44A	Chip Date LACTOSCAN 45 Sep 19.82,2915 11.53.47 LACTOSCAN 45 Sep 19.82,2915 12.512 LACTOSCAN 45 Sep 19.82,2915 13.512	Client Varosha183 Varosha181 Varosha181 JM OC099 JM OC0992 JM OC0992 JM OC0992 JM OC0992 JM OC0991	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilation 1 5 5 1 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	Calls per ml. 104-10° 001-10° 335-10° 335-10° 5023-10° 72-10° 72-10° 72-10° 565-10°	Ang. cell size [pn] 10.47 11.12 12.24 9.15 12.24 9.15 12.43 11.39 10.51
10 28C 28C 32C 32B 33B 33B 43C 43C 44A	Chip Date LAC 10542AB 4E Natr 1982.2016 11.488-44 LAC 10542AB 4E Natr 1982.2016 11.388-44 LAC 10542AB 4E Natr 1982.2016 11.338-44 LAC 10542AB Natr 1982 11.382.2016 11.358-94 LAC 10542AB Natr 1982 11.382.2016 11.358-94 LAC 10542AB 4E Natr 1982.2015 11.3126-94 LAC 10542AB 4E Natr 1982.2015 11.3126-94 LAC 10542AB 4E Natr 1982.2015 11.3126-94 LAC 10542AB 4E Natr 1982.2015 11.3126-94	Client Vansha101 Vansha103 Vansha103 Vansha104 Vansha104 JIM 00009 JIM 000093 JIM 000093	Average cells Milk Preserved cow milk Preserved cow milk	Print, 72-103 Difution Calls per int, 1 Calls for 1 Store of the second	Avg. cell size (pre) 10.22 10.47 11.42 11.42 11.51 12.04 9.55 12.43 11.55 10.57 10.57	ID 29C 31C 328 34A 298 429 429 428 45 44A 48A	Chip Chip Cartoscan es Ser Lactoscan es	Client Varosha193 Varosha193 JM OCD99 JM OCD99 JM OCD99 JM OCD99 JM OCD99 JM OCD99 JM OCD99	Average cell Milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk Preserved cow milk	Dilation 1 1 1 1 1 1 1 1 1 1 1 1 1	Cells per ml. 184-10° 001-10° 335-10° 5023-19° 72-10° 72-70° 565-10° 565-10° 565-10°	Arg. cell size [art] 16.47 17.24 17.24 17.24 17.24 17.25 17.43 17.55 17.43 17.55 17.55 17.55 17.55 17.55
10 23C 25C 37C 37C 37C 37C 37C 37C 37C 37C 37C 37	CMI Chain LACTOCAL # Styp: TRAZINE THE #### LACTOCAL # Styp: TRAZINE THE #### LACTOCAL # Styp: TRAZINE THE ##### LACTOCAL # Styp: TRAZINE THE ##################################	Client Varusha101 Varusha101 Varusha101 JiM codbi JiM co	Average cells Milk Preserved cow milk Preserved cow milk	Print. 72-103 Ditution Cells per mt. 1 143-19 1 161-19 1 355-19 1 356-19 1 356-19 1 356-19 1 356-19 1 356-19 1 352-19 1 352-19 1 352-19	Aug. cell size (ani) 10.22 10.47 11.52 12.04 9.55 12.43 11.55 10.57 10.57 10.57 10.57 10.55	ID 29C 31C 34A 298 428 435 44A 44A 48A 451 51C	Clap Date LACTOSCAM 65 (w) 19.02,2945 (19.52,294 LACTOSCAM 65 (w) 19.02,2945 (19.55,294 LACTOSCAM 65 (w) 19.02,2945 (19.55,194 LACTOSCAM 65 (w) 19.02,2945 (19.52,194 LACTOSCAM 65 (w) 19.02,2945 (19.24,194	Client • Varosha193 • Varosha191 • Varosha191 • JBI OCD92 • JBI OCD92 • JBI OCD92 • JBI OCD91 • JBI OCD91 • JBI OCD91 • JBI OCD91 • JBI OCD91	Average cell Milk Preserved cow milk Preserved cow milk	per 7377-1 Dilution 1 5 5 5 5 0 0 0 0 0 0 1 1	Calls par ml. 18410° 80110° 33510° 30910° 5023-10° 5023-10° 5023-10° 5023-10° 56510° 85510° 85510° 85510°	Avg. cell size [ars] 16.47 11.12 12.24 9.15 12.24 9.15 12.24 11.25 11.25 11.25 11.55
0 26 26 26 36 36 36 36 40 40 40 40 40 40 40 50	Chip LLC TOTAL ARE Sign 1982 2981 14 68-14 LLC TOTAL ARE Sign 1982 2981 14 68-14 LLC TOTAL ARE Sign 1982 2981 14 58-34 LLC TOTAL ARE Sign 1982 2981 14 58-39 LLC TOTAL ARE Sign 128.23 2981 14 58-39 LLC TOTAL ARE Sign 128.23 2981 14 58-39 LLC TOTAL ARE Sign 128.23 2981 14 58-34 LLC TOTAL ARE Sign 128.24 2981 14 58-34 LLC TOTAL ARE Sign 128.24 2981 14 58-34 LLC TOTAL	Client Varisha101 Varisha101 Varisha101 JM od0092 JM od0092 JM od0092 JM od0092 JM od0092 JM od0092	Average cells Milk Preserved cow milk Preserved cow milk	Performance of the second seco	Avg. cell size (pre) 10.22 10.47 11.12 12.04 9.15 12.43 11.57 10.57 10.57 10.51 10.59 9.42	10 29C 31C 328 34A 398 428 428 428 428 428 428 428 428 548	Chep Date LACTOSCAM 65 (apr. 19.2.344) 19.2.314) LACTOSCAM 65 (apr. 19.2.344) 19.2.3341	Client Varosha193 Varosha191 Varosha191 Varosha191 Varosha194 Varo	Average cell Mile Preserved cow mile Preserved cow mile	per 7377-1 Dilation 1 5 5 0 Delete 1 5	Calls per ml. 164-10° 001-10° 335-10° 335-10° 5023-53° 72-10° 72-70° 565-10° 555-10° 555-10° 12-10° 73-10° 555-10° 12-10	Arg. cell size [are] 10.47 11.12 12.04 9.15 12.43 11.55 11.55 10.51 10.51 10.51 9.42 9.42 9.42
10 226 226 226 326 344 435 444 455 455 444	Chip Den LCPCIDCAL 41 (Reg. 10.2224 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 10.2224 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 10.224 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 10.224 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 10.224 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 10.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1 LCPCIDCAL 41 (Reg. 12.24 HT 11.14 at 1 1	Client Varisha101 Varisha101 Varisha101 Jilli OOD03 Jilli OOD032 Jilli OOD032 Jilli OOD032 Jilli OOD031 Jilli OOD031 Jilli OOD031 Jilli OOD031 Jilli OOD031	Average cells Mile Preserved cow milk Preserved cow milk	Pric. 72-103 Ditution Cells per ret. 1 U8-10 1 U8-10 1 350-10 1 350-10 1 350-10 1 550-10 1 7277-10 1 650-10 1 650-10 1 7277-10 1 850-10 1 850-10 1 7277-10 1 850-10 1 850-	Asg. cell size (and) 10.22 10.47 11.52 12.04 9.55 12.43 11.57 10.57 10.57 10.57 10.55 9.42 9.60	ID 29C 37C 328 34A 278 428 428 428 428 457 44A 48A 57C 548 550	Choir Date LACTORIAN & Style RE32(39) RE32(39) LACTORIAN & Style RE32(39) RE32(39) LACTORIAN & Style RE33(39) RE33(39) LACTORIAN & Style RE33(39)	Client Varosha131 Varosha131 Varosha134 JBI OCOD92 JBI OCOD92 JBI OCOD92 JBI OCOD92 JBI OCOD92 JBI OCOD92 JBI OCO122 JBI OCO122	Average cell Mai Preserved cow mili Preserved cow m	per 7377-1 Dilution 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	Calls per ml. 164 10° 001 10° 355 10° 355 10° 355 10° 355 10° 355 10° 355 10° 172 10° 7377 30° 565 10° 622 10° 73 10° 120°	Arg. cell size (an) 10.47 11.12 12.24 3.15 12.43 11.55 12.43 11.55 14.55 14.55 14.55 14.55 14.99 9.42 9.49 11.56
0 286 286 376 377 378 378 444 450 450 450 450 450 450 450 450 450	Cally Date LCTOCICAL #1 499, 195.2014 150.2014 LCTOCICAL #1 99, 105.2014 150.2014 LCTOCICAL #1 99, 105.2	Client Varisha101 Varisha101 Varisha101 Varisha101 JIM 000092 JIM 000092 JIM 000092 JIM 000092 JIM 000092 JIM 000092 JIM 000122 JIM 000123	Average cells Mile Preserved cour mile Preserved cour mile	Print, 72-103 Ditation Cells per ret. 1 Cells per ret. 2 Cel	Aug. cell size janij 10.22 10.47 11.52 12.44 12.55 12.43 10.57 10.57 10.57 10.57 10.57 10.57 10.57 10.57 10.57 10.57 10.57 11.56	10 29C 31C 328 34A 238 428 428 428 428 428 428 548 548 548 558 550 786	Chip Date ULXTOTICAM 69 (m) R.2.291 (m) 53.4 ULXTOTICAM 69 (m) R.2.291 (m) 53.4 LXTOTICAM 69 (m) R.2.991 (m) 53.4	Client Varolha193 Varolha191 Varolha191 Varolha191 Varolha191 Vall OCD92 Vall OCD92 Vall OCD92 Vall OCD91 Vall	Average cell Mile Preserved cow mill Preserved cow mill	Per 7377-1 Dilation 1 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	Colls per rel. 154-10° 355-10° 306-10° 306-10° 3028-90° 7276-10° 555-10° 555-10° 555-10° 555-10° 522-10° 7376-10° 555-10° 523-10° 555-10° 523-10° 555-10° 523-10° 555-10° 55	Avg. cell size jare] 10.47 11.12 12.24 1.15 17.5 17.5 17.5 17.5 17.5 17.5 17.5

2. Press the Delete button and confirm that you want to delete the selected result

	Graph Options Advan		23.02.2018 16:53
		Start	
		ID 43C Client JIM OOD099	
		Average cells per 7377.10 ³	
	Chip Date Client	Milk Dilution C	alls per mi. Avg. cell size (jun) *
	CMp LACTOSCAN 4R Styp 19.82.2016 11:33-83 v Varosha10	Milk Dilution C Preserved cow milk 1 1	ells per ml. Avg. cell size (prn) *
0.2	CNp Date Client LACTOSCAN 48 Styr 19.92.2016 11:33:33 v Vanosha10 LACTOSCAN 48 Styr 19.92.2016 12:05:24 v Vanosha10	Milk Dilution C Preserved cow milk 1 1 1 Preserved cow milk 1 8	Cells per mil. Avg. cell size [am] Aug. 84-10* 10.47 101-10* 11.12 62-10* 10.10
	Chip Date Chient LACTOSCAN 45 Styr 19.02.2016 11:53:43 v Varesha16 LACTOSCAN 48 Styr 19.02.2016 12:156:28 v Varesha16 LACTOSCAN 48 Styr 19.02.2016 12:156:28 v Varesha16	Milk Diffusion C Preserved cow milk 1 1 Preserved cow milk 1 3	ells per ml. Avg. cell size (jml) 410° 10.47 91-0° 11.12 5510° 12.04
	CMp Date Client LACTOSCAN #8 Step: 19.02.2016 153.43 Vacenball/L LACTOSCAN #8 Step: 19.02.2016 152.65.20 Vacenball/L LACTOSCAN #8 Step: 19.02.2016 152.65.20 Vacenball/L LACTOSCAN #8 Step: 19.02.2016 152.65.20 Vacenball/L LACTOSCAN #8 Step: 12.02.2016 159.25.90 180.0000 LACTOSCAN #8 Step: 12.02.3016 159.25.90 180.0000	Mile Citution Preserved cow milk 1	alfa per mt. Avg. cell size juni 64-10° 10.47 81-10° 11.12 55-10° 12.04 59-10° 3.15 50-47 12.4
	Chip LAC (195CAM 46 Spir (19.2) 246 (19.5) 4: Versenheit LAC (195CAM 46 Spir (19.2) 246 (19.6) 4: Versenheit LAC (195CAM 46 Spir (19.2) 246 (19.6) 5: Versenheit LAC (195CAM 46 Spir (19.6) 5: Versenheit LAC (195CAM 46 Spir (19.	Mili Difution C Preserved cow milk 1 1 8 Preserved cow milk 1 8 Preserved cow milk 1 3 Preserved cow milk 1 3 Preserved cow milk 1 3	Alla per mt. Avg. cell size (and 84-10° 10.47 01-10° 11.12 55-10° 12.04 80-10° 25.5 829-10° 12.43 3247 11.28 11.29 11.20
	CMp Daw Clinet LACTOSCAR 459 19.22616 13.32.4 Neunshoft LACTOSCAR 459 19.22616 13.32.4 Neunshoft LACTOSCAR 459 19.22616 13.52.4 Neunshoft LACTOSCAR 459 19.22616 13.53.4 Neunshoft LACTOSCAR 459 12.22616 13.53.4 NM CODD LACTOSCAR 459 12.22616 14.64.9 NM CODD LACTOSCAR 459 12.22616 14.64.9 NM CODD LACTOSCAR 459 12.22616 14.69.9 NM CODD LACTOSCAR 459 12.22616 14.69.9 NM CODD LACTOSCAR 459 12.22616 14.69.9 NM CODD	Dilution Citation C Preserved cow milk 1 1 1 Preserved cow milk 1 3 3 Preserved cow milk 1 5 3 Preserved cow milk 1 5 3	Alls per ml. Aug. cell size (an) 64-10° 10.47 61-10° 11.12 55-10° 12.04 59-10° 3.15 592-10° 12.43 24.19° 12.43 24.19° 12.43
	Chie Date Chief LACTOCICA 48 59 (19.2):2115 12.01-21 LACTOCICA 48 59 (19.2):2115 12.01-21 LACTOCICA 48 59 (19.2):22115 12.01-21 LACTOCICA 48 59 (19.2):2215 12.01-21 LACT	Hit Cludion C Personnel Coar milli 1 1 1 Personnel coar milli 1 1 Personnel coar milli 1 1 Personnel coar milli 1 1 Personnel coar milli 1 7 Personnel coar milli 1 7 Personnel coar milli 1 7	Alls par rd. Ang. cell size (pr.) 44-10° 10-21 47-10° 11-12 47-10° 12-24 49-10° 12-24 49-10° 12-24 49-10° 12-24 49-10° 12-24 49-10° 12-24 49-10° 12-24 49-10° 10-51 40
	Chip Daw Cliner LACTOSCAR # 59 (m.22316 133.0 c) // worshold LACTOSCAR # 59 (m.22316 133.0 c) // worshold LACTOSCAR # 59 (m.22316 13163 c) // worshold LACTOSCAR # 59 (m.22316 13163 c) // worshold LACTOSCAR # 59 (m.23316 1556 c) // with OODS LACTOSCAR # 59 (m.23316 1556 c) // with OODS	Mik Dilation C Preserved core mik 1 1 1 Preserved core mik 1 8 Preserved core mik 1 3 Preserved core mik 1 3 Preserved core mik 1 5 Preserved core mik 1 5 Prese	Alth per md. Any. cult size (mr) 44 10° 10.47 55 10° 12.04 5910° 12
	Chip Date Claim LAC1052AR 48 Sty 19.22 214 33.34 */ Venoshaft LAC1052AR 48 Sty 19.22 214 33.44 */ Venoshaft LAC1052AR 48 Sty 19.22 214 51.14 */ Venoshaft LAC1052AR 48 Sty 19.22 214 51.14 */ Venoshaft LAC1052AR 48 Sty 19.22 214 51.14 */ Venoshaft LAC1052AR 48 Sty 19.23 214 51.14 */ Venoshaft LAC1052AR 48 Sty 19.23 214 51.15 */ Venoshaft LAC1052AR 48 Sty 19.23 214 51.15 */ Venoshaft LAC1052AR 48 Sty 12.33 Venoshaft Venoshaft	Mh Diblos Preserved core min Preserved core min Pre	alli per nt., Ang. cult size (pr.) 144 - 97 - 16.47 147 - 97 - 16.27 159 - 97 - 16.28 159 - 97 - 16.29 159 - 97 - 16.29 159 - 179 - 16.51 152 - 179 - 1
	Chip Date Chart LLC ¹ (S)CLM Step 101.23 Step 101.23 Step 101.23 Step 101.23 LLC ¹ (S)CLM Step 102.23 St	Mih Dhulon C Present oru mih 1 1 8 Present oru mih 1 8 Present oru mih 1 8 Present oru mih 1 3 Present oru mih 1 3 Present oru mih 1 7 Present oru mih 1 7 Present oru mih 1 8 Present oru mih 1 8	Alli par nt. Any, cell size [and 44-10" 10.47 61-90" 11.52 55-97" 12.04 55-97" 12.04 55-97" 12.03 525-11" 12.04 525-11" 12.04 52
	Date Date Chart LC COSCAR 10 (%) State 100 (110 (147) (147	His Diudios C Preserved core mills	Alls par ref. Any, cell size (pref. 54 47 16.47
C C 8 8 8 8 8 8 8 8 8 8 0	Cop Data Data Latic DOLLAR ID 50 (2014) 50 (20	Kin Charlon Charlon Preserved coru mil 1 Preserved coru mil Pr	Alls per nt. Asy. cell size greg 44 197 to 4.2 110 07 to 12.2 300 197 to 12.2 300 197 to 12.2 300 197 to 12.2 210 10 10 10 10 10 301 10 10 10 301 10 10 10 301 10 10 301 10 10 301 10 10 301 10 10 301 10 10 301 10 30

3. The selected result is deleted.



Cleaning of the Database

The software allows in accumulation of a large number of results of analyzes, the entire database to be deleted to make space for storing the results of the upcoming analysis.

To clean the Database:

1. Go to menu Options



3. From the drop-down menu select Clean Database

LSCC - Lactoscan Somatic Cell Counting, v1.2.0.4	2	LSCC - Lactoscan Somatic Cell Counting, v1	20.42
User options User	Processing Optics Brightness Focusing	User options	Processing Optics Brightness Pocusing
Probe type for which should apply parameters over nik, most with lots firest	All values on this page are associated to satisfied proba type. The other proba types may have their own different values.	Probe type for which should apply paramete Criv rill, read with lots times Cow milk	All values on this page are associated to selected probe type. The other prote types may have their own different values.
Minimum cell size in micrometers The task better this size without the countert	Warning level of cells count (-10%mL) Pates with non-cells are colored in page E00 E00 E00 E00 E00 E00 E00 E00 E00 E0	Minimum cell size in micrometers The one tenso the size without the counted	Warning level of cells count (-101mL) Pates with work tabled to campe ECO
Maximum cell size in micrometers The only ensuring size without the country	Critical level of cells count [101/mL] Pistes with new cells are celled in red	Maximum cell size in micrometers The cells size with ut be counted	Critical level of cells count (101mL) Prides with more calls are called in ref
20 .	500	20 -	500
The upper limit of cell size in micrometers The instance cell size (and surpass the user 20	Laboratory offset in percentage the seaso-se positive function, or zero for normal instit	The upper limit of cell size in micrometers The noninent set sort server burgest this reter 20	Laboratory offset in percentage One register stronger, at zon for normal insut 0
	About Activate Window Caption bar Maintina		Alanat Animata Ulindran Caption har Matinina
	feat Care Database		Reat
	Opinist Database results Yead Calls Options		Update Extension enable Year California

4. Press Clean Database





Attention!

Once deleted, the results cannot be restored!


Printing the results from the data base

LACTOSCAN SCC gives the opportunity to print every single result in the data base of the device.

Printing the results from an integrated printer: 1. Choose the result which needs printing.

It is marked by pushing onto the information for the relevant analysis in the field "Summary".

			art (23.02.2018 16	re: Main	Graph) 🕒 Options			Sta	rt	23.02.2018 1
		ID 42B Client JIM OOF Milk type Preserv Average cells per 72-10	2094 ed cow milk						CI Mik t Average colls	235B T.UZUNO Preserve	v146 ed cow milk	
	Data Client		Collegend	And and also first						mil to be		
I LACTOS	5C 3N 42 50 m 10 02 2016 11 40 46 m Maxweb 101	Desenand con milk 1	ATTA-10	10.22	10	LACTORCAN	U809	Chenz	Descent of a second ball	Difusion	Cells per mil	Avg. cen size juni
C LACTOS	SCAN 48 50w 19.02 2016 11:53:43 v Varosha103	Preserved cow milk 1	154-10*	10.47	1210	LACTOSCAN	48 50pr 10.03.2016 09.41.53	* Agropanto	Preserved con mik		85.10	12.05
C LACTOS	CAN 4R 50ar 19.02.2016 12:05:20 v Varosha101	Preserved cow milk 1	801-10*	11.12	1250	LACTOSCAN	4R 50ar 10.03.2016 09:42:16	* Agropartn	Preserved coar milk	1	3393-104	12.13
8 LACTOS	SCAN 4R 50ur 19.02.2016 12:10:52 v Varosha104	Preserved cow milk 1	355-10"	12.04	1448	LACTOSCAN	4R 50pr 10.03.2016 11:19:40	* Agropartré	Preserved coar milk	1	1631-104	11.89
A LACTOS	JCAN 4R 50pr 02.03.2016 10:35:59 v JIM OOD09	Preserved cow milk 1	390-10*	9.15	2166	LACTOSCAN	4R 50ar 15.03.2016 11:00:13	* T.UZUNOV	Preserved coar milk	1	24-10*	7.98
B LACTOS	JCAN 4R 50pr 02.03.2016 10:56:09 ~ JIM OOD092	Preserved cow milk 1	5929-10*	12.43	216C	LACTOSCAN	4R 50µr 15.03.2016 11:01:01	N T.UZUNOV1	Preserved cow mllk	1	129-10"	9.90
B LACTOS	CAN 4R 50pr 02.03.2016 11-15:21 ~ JIM OOD094	Preserved cow milk 1	72.58	11.35	2160	LACTOSCAN	4R 50pr 15.03.2016 11:01:23	 T.UZUNOV2 	Preserved cow milk	1	239-10"	11.05
IC LACTOS	ICAN 4R 50pr 02.03.2016 11:22:14 + JIM OOD099	Preserved cow milk 1	7377-104	10.57	219C	LACTOSCAN	4R 50pr 15.03.2016 11:14:54	 T.UZUNOV3 	Preserved coar milk	1	1379-104	11.95
IA LACTOS	ICAN 4R 50pr 02.03.2016 11:28:34 + JIM OOD091	Preserved cow milk 1	565-10*	10.51	234C	LACTOSCAN	4R 50pr 15.03.2016 12:58:00	* T.UZUNOV143	Preserved cow milk	1	3719-10*	10.90
LACTOS	sCAN 4R 50pr 02.03.2016 11:48:18 v JIM OOD097	Preserved cow milk 1	832-10*	14.99	2358	LACTOSCAN	4R 50pr 15.03.2016 13:00:43	< TUZUNOV146	Preserved cow milk		1055-10*	10.08
LACTOS	SCAN 4R 50µr 02.03.2016 12:17:28 v JIM 000121	Preserved cow milk 1	79-107	9.42	₽ 11 2410	LACTOSCAN	4R 50pr 15.03.2016 13:28:04	v VAROSHA11	Preserved cow milk	1	3716-10*	11.69
8 LACTOS	SCAN 4R 50pr 02.03.2016 12:32:23 + JIM 000122	Preserved cow milk 1	326-10*	9.00	t: 1 241D	LACTOSCAN	4R 50pr 15.03.2016 13.31.18	v VAROSHA11	Preserved cow milk	1	631-10*	11.37
d0 LACTOS	ICAN 4R 50pr 02.03.2016 13:33:27 - JIM 000164	Preserved cow milk 1	6580-10*	11.56	2448	LACTOSCAN	4R 50pr 15.03.2016 13:45:40	+ VAROSHA119	Preserved cow milk	1	485-10*	9.35
IA LACTOS	sCAN 4R 50pr 02.03.2016 13:35:16 + JIM 000165	Preserved cow milk 1	6633-10*	11.97	2 2440	LACTOSCAN	4R 50pr 15.03.2016 13:46:26	v VAROSHA111	Preserved cow milk	1	190-10*	12.02
90 LACTOS NA LACTOS	CAN 48 56µ (22.03.2016 13.33.27 × J.M. 00.0144 CAN 48 56µ (22.03.2016 13.35.16 × J.M. 00.0145	Preserved cow milk 1 Preserved cow milk 1	6500-10 ⁹ 6633-10 ⁴	11.56 11.97 22.07	 244D 244D	LACTOSCAN	i 48 50µr 15.03.2016 13.43.40 i 48 50µr 15.03.2016 13.46:26	 VAROSHA119 VAROSHA111 	Preserved cow milk Preserved cow milk	1	485 10° 190 10°	9.35 12.02



3. It is not necessary to choose a printer or form of printing. By default it is accepted that LACTOSCAN SCC will print the results on the integrated printer in a short form.

4. Push the arrows from the right side of the results ID if changing or adding new results for printing is needed.



5. If printing the results from the analysis, made for a certain period of time is needed, then, with the assistance of a calendar the initial and the final date of the period can be set.





Printing the results on an external printer:

1. Before start printing, make sure that the external printer is connected and installed.

2. From the drop-down menu, choose the name of the external printer.



3. Choose Wide form for printing the results in a table.



4. Press the arrows from the right side of the results ID, if changing or adding new results for printing is needed.



Saving the results in PDF file:

1. From the drop down menu, choose MICROSOFT PRINT TO PDF.



2. Before saving, check the orientation of the page.

If you want to save the results in WIDE form, the orientation must be landscape and the paper size A4. To check the orientation, press Document properties.



3. Choose Wide form for printing the results in a table.



4. Add the results that you want to include in the PDF file.

Utered Plate In PC Decement Properties Utered Plate In PC Decement Properties Utered Plate In PC Decement Properties Utered Plate In PC Decement Plate In PC	tton	<u>ра</u> 7272331534 1944	Parts	eve the fil	Image: constraint of the second of the se	0 0 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3 3	Dah 2005/141 12:03.04 2005/041 12:03.04 2005/041 12:03.04 2005/041 12:03.04 2005/041 12:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:03.04 2005/041 10:04.04 2005/041 10:03.04 2005/041 10:04.04 2005/041 10:03.04 2005/041 10:04.04 2005/041 10:03.04 2005/041 10:04.04 2005/041 10:04.04 2005/041 10:04.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 2005/041.04 2005/041 10:04.04 <	Cited Spann Spann Spann Spann Spann Spann Spann Spann Spann	Page 11
Picture Contract to PCF 2 Contract Projections Contract Projections Output Source Contract Projections Contract Projections Contract Projections Projections Contract Projectio	0 10 14 16 16 16 16 16 16 16 16 16 16 16 16 16	Dat Open 7.06218 11.20.4 Faces 7.07218 11.40.4 Saces 7.07218 11.40.4 Saces	Page 15	⇒	Sechard Dolgshi → →	× P E + 0 V V V	Extel 7000000 7000000000000000000000000000000000000	Claint System System System System System System System System System	Page 11

6. Select the folder where you want to save the file, write the name of the file and press button SAVE to finally save the results in PDF file.

The Save Print Durgut Ac X				Printer					
+ +				CutePDF Writer -					
Organia • Navehidar III • 0		Page 1	14	Document Properties			Page 1/1		
This PC * No forms match your angests	finite.	Class	24		it:	Date	Chert	Probe	De
E Desites	201/2316 13 25 44	funtum.		International Control of Control	3250	22.02.2018 15	4 0003-0001014	Cow mile	281
S Decuments		address of the second sec		Output format					
- Downlands	7/26/2018 10:35 AM	System		() Short form					
J Mac	7/26/2016 10:51 AM	System							
Witzen	7/26/2015 10:52 AM	System		Vince form					
Windows (C)	7/25/2015 10:57 AM	System							
(P)	7/26/2016 10:53 AM	System		Selected probe(s)					
_ (P)	7/27/2016 9-13 AM	Sutem		Probe ID					
the same	7/27/2018 0 14 AM	System		From: 325D					
Serve as type: (PDF Document (* adt)	7/27/2016 9 14 AM	Suttern							
	7.97(2314 B 15 AM	System		To: 325D					
n HitsTukles Censil				Durba data					
Probe date				Probe date					
TOP 201 1 10 20 20 20 20 20 20				From: 22.02.2018 • 15:19:23 ч. ÷					
PTON: 7/20/2010 10.30.20 PM									
To: 7/26/2016 + 10:30:20 AM				10: 22.02.2018 • 16:19:23 v. 🗧					
				1 proba(s) selected, 9 page(s)					
in protect and the property				Comments and a second s					
Cancel Exit				Cancel Cancel					
Calicer			-		4 million				
			Prov.	and the second se					

Exit

Saving the results in Excel file:

1. Go to printing mode by pressing the button





2. Add the results you want to be included in the Excel file.





Printer	Thi Sove Print Diagonal Sc → → → → → → → → → → → → → → → → → → →	×
Microsoft Print to PDF	The Part of the Second Se	E - O Papita
Output former One One One Output former Image: Status Sta	Remarking devices and a second device and a se	Data Operation 20202111 02384 Apples 20202111 02384 Apples 20202111 02384 Apples 20202111 02284 Apples 20202111 02284 Apples 20202111 02284 Apples 20202111 02284 Apples 20202111 02384 Apples 2020211 02384 Apples 2020211 02384 Apples 2020211 02384 Apples

4. Select the folder where you want to save the file, write the name of the file and press button SAVE to finally save the results in Excel file.

Save Print Dutyot & X		The Save Print Durget As	
	-		
legenia = New folder 📰 = 😖	Date 17	Dyprine - Nav folder III - O	Ence 17
This PC. No thema match your search.	Date Client	The PC No Item and A year parts.	Date Client
Desktop	7/26/2316 10:35 AM System	Custop .	7/26/23/16 10:35 AM System
Development	7/96/2016 18/38 AM Sustem	L Constants	7/96/2016 10:36 444 Sustem
J. Main	Topological states and the states	J. Marine	TOLOGIA INC. INC.
Pictures	Transaction to at the second second	all Pictures	All and a second s
Witem	7/26/2016 10:52 AM System	Vitem	7/26/2015 10:52 AM System
L Windows (C)	7/25/2015 10:53 AM System	🛓 Windows (C)	7/25/2016 10:53 AM System
= P)	7/26/2016 10:53 AM System	= P)	7/26/2016 10:53 AM System
(P)	7/27/2016 0 13 AM System	(D) v	7/27/2016 E 13 AM System
Henane	7/27/2016 0 14 AM System	Titegane Reula	7/27/2016 0 14 AMI Eystem
Seve as type: PDF Decument (*.pdf)	7/27/2016 0 14 A88 System	Serve as type: IPDF Decument (*gelt)	7/27/29/16 0 14 A88 System
Hide Folders Cancel	7/07/02198-0-15-AAR Bystem	n Hele Fulders Servi Censel	7/07/0316 0.15 AM System
Probe date		Probe date	
From: 7/26/2016 10:30:20 AM		From: 7/26/2016 • 10:30:20 AM ÷	
To: 7/26/2016 • 10:30:20 AM ÷		To: 7/26/2018 • 10:30:20 AM 🚊	
10 probecký selectad. O pagedij		10 probacity solucional. D pageopi	
Cancel Exit		Cancel Exit	
	- F		

Focusing of LACTOSCAN SCC

Focusing of the device must be done every 3 months.

In order to focus the device:

1. Preparing cow milk sample with 300-500x10³ somatic cells per milliliter.

It may also be used a cow milk sample from a previous analysis with a result between 300×10^3 and 500×10^3 somatic cells per milliliter.

If such a cow milk sample is not available, it must be used such cow milk sample for which it is supposed that the presence of somatic cells is between 300×10^3 and 500×10^3 somatic cells per milliliter. The cow milk sample must be prepared for analysis as per the instructions from section "Preparing of the sample for analysis".

2. Focusing of LACTOSCAN SCC



2.1. From the screen "Options" press the button **Sector Methods**, in order to proceed to the screen "Focusing". In the four squares the last found focuses are seen.



2.2. From the drop-down menu "Type of sample", choose cow milk.



2.3. Choose the type of the microfluidic camera which will be loaded with the sample.



2.4. Press the button 💛, so that the cartridge can be shown.



2.5. Place the loaded with the milk sample LACTOCHIP in the cartridge.







The device takes the cartridge in and focusing can be started.





The device starts searching the focus of the four main points.





ACCEPT DISCARD Reject this result

ACCEPT Save this result

DISCARD Reject this result

Work in Advanced mode:





2. Place the LACTOCHIP with the sample in the cartridge and choose its type.



3. Press, to close the cartridge with the LACTOCHIP in it



4. Choose a sector in the microfluidic camera, which is needed to be displayed.

In the sample it is microfluidic camera A, section 4



5. To make a picture and see the image of the chosen section, press



6. To see 3D graphics of the image, press . The 3D graphics shows the cells in form of columns, as the height of the column depends on its brightness.



7. In order to see a 3D graphics of group of cells from the image, move the grey square from the upper left corner on the needed group of cells. The mouse could also be used, position it on the grey square, press and keep pressed the right button and move the mouse over the image.



8. All the sides of the 3D graphics could be seen by pressing on it and moving up/down/left/right until the required field is seen. The mouse can be used by positioning on the 3D graphics; pressing and keeping the right button pressed and moving the mouse.



9. To see all cells of the image marked and to see the threshold level for mastitis cells, yeast cells or all cells, press button .

Mastitis cells





Yeast cells





All cells





10. For magnification, press 📙

and choose % of magnification.









The User can select the format of the saved image. Possible formats are .TIFF and .BMP. The .BMP format is easy for further process and can be opened by or imported in any Windows program.



12. To open an image from previous analysis, press **Copen...**. Open the directory where the image was saved, click on it and press <u>Open</u>.



13. Button Stream – after pressing the button Stream, mode Movie is activated and LACTOSCAN SCC starts making consecutive pictures as the time for changing the image depends on the set exposition time. For example, if it is set to 1.945 sec. exposition time, images on the display will be changed in approximately 1.945sec. Mode Movie is used when a number of sections of

the microfluidic camera need to be seen. Instead of pressing the button after each new section chosen, LACTOSCAN SCC is automatically filming the corresponding section. While the mode Movie is active, the user can not only choose a section of the microfluidic camera, but also change the intensity of the LED, exposition time, change the focus and gain sensitivity of the optical sensor (CCD GAIN).



14. By pressing , different sections of the microfluidic camera can be chosen. Another section can be chosen by pressing its corresponding number. To make a picture and see the image of the

newly chosen section, press



In the example: section 10, microfluidic camera A.



15. For setting the intensity of the LED 1 and 2 (if they are available) and exposition time, press



16. To switch off/on one or both LEDs or to locate best brightness, press



17. To see the buttons used for setting the focus of the microscope,



Main	🖆 Graph	Options Advanced	23.02.201	8 10:10 ч.	Main Graph Options	Advanced	23.02.2018 10:36 ч.
Caracterization Concernation Dialon Probe 1 Cow Concernation Concer	10. 2 A COLO 20 AND A COLORADO	Contraction of the second			Protect Line And And Andrews And Andrews And Andrews A	A B B C C C C C C C C C C C C C C C C C	Parkana 3.05 Ju







manually set the autofocus



Attention!

In the tab "Step" is shown the step with which the position of the focus will be changed. The value of the step can be changed by writing the corresponding number in the tab.



20. In order to see the graphics of the results of the current image, before passing to menu Graph, the dilution and type of the sample must be entered.



Check of LACTOSCAN SCC magnification using a graduated scale

The check of LACTOSCAN SCC magnification is done with the use of a special LACTOCHIP with a graduated scale.

1. In order to perform the check, in menu"Advanced" press the button





2. Press the button which to open the menu for positioning of LACTOSCAN SCC microscope on the graduated scale.





3. Press

to position the microscope on the graduated scale.





3.1. Use Motor Y and Мотор X for a more precise positioning of the microscope on the graduated scale.

🖙 Flash

3.2. Button Flash opens a menu, which is used only by a trained technician. All changes, made in this menu, may lead to malfunctions of the device.



4. Press the buttor

to open the menu for setting of the intensity of LED.



For an automatic finding of the best intensity for the LED, press the button W. If the image doesn't seem light enough or the graduated scale seems too dark, use the sliding button to set the exposition.



in order to proceed to the menu for check.

5. Press the button



The values for Focus rate are different at a different intensity of the LED.

Position the mouse on the left horizontal scale, on position 0. Press the right button of the mouse and drag it to position 6 of the right horizontal scale. This way you will draw a measuring line.

For example: The Line size must be 1.400 mm $\pm 0,01 \div 0,03$ mm. When the Line size is different from 1.400 mm $\pm 0,01 \div 0,03$ mm, then the LACTOSCAN SCC magnification is not set. It is necessary to contact LACTOSCAN service center.



5.2. In order to be sure that the line is correctly positioned on position 0 and position 6 from the scale, press the button to zoom.



5.3. The button ^{Tune magnification factor} opens a menu, which is used only by a trained technician. All changes, made in this menu, may lead to malfunctions of the device.

Updating the information for the types of LACTOCHIP available in the software

Updating the types of LACTOCHIP available in the software is made with the update file. If updating the information is needed, please contact our service.

Exit and switch off the LACTOSCAN SCC

1. In order to proceed to the screen Switch off, press the button



2. In order only to exit the device's program, proceed to the operational

system Windows, and press

3. In order to turn off the program and the operational system of the device,







Attention!

Always turn off the device at the end of the working day or after finishing work with the SCC.

4. To go back to the program, press



Replacing LACTOSCAN SCC filters

The filters, situated at the rear panel and at the bottom of the LACTOSCAN SCC need replacement each month on regular basis or in case it is noticed that they are obviously dirty.

Replacing the filters of the rear panel:

1. Unscrew the screws of the filter's grid.



2. Remove it.



3. Remove the filters.



4. Place the new filters.



5. Place the grid back.



6. Screw the screws back.



In order to replace the filter at the bottom of the device, follow the instructions about replacement of filters at the rear panel, described in part "Replacement of LACTOSCAN SCC filters.

Maintenance and cleaning

LACTOSCAN SCC needs a periodical auto focus.

Change the filters of LACTOSCAN SCC once a month.

Clean the box of LACTOSCAN SCC with a soft cloth and alcohol or distilled water paying attention the cleaning liquid not to get into the power port, USB ports or through the door of the slide. Do not spill over or splash any kind of liquid directly on the device in order to avoid a power failure, when the device is on. Wipe the screen with a dry cloth right after cleaning.

Clean the sensor display of LACTOSCAN SCC with a soft lint-free cloth, moistened with a cleaning liquid for LCD displays without applying any force and being careful the liquid not to get into the buttons placed on the sides of the screen. Do not spill over a cleaning detergent directly on the screen and do not use any abrasive products, which may scratch the screen. Wipe the screen with a dry cloth right after cleaning.

Disinfect LACTOSCAN SCC by using a soft cloth and 70% alcohol paying attention cleaning liquid not to get into the power port, USB ports or through the door of the slider. Avoid using bleach because it may leave white leftovers on the device. Do not pour or splash any kind of liquid directly on the device in order to avoid a power failure when the device is turned on. Wipe the screen with a dry cloth right after the cleaning.



Attention!

Throw out the cloths in an appropriate waste container.

LACTOSCAN SCC - determination the number of the somatic cells in raw milk - method



The number of somatic cells (SCC) is one of the internationally recognized standards for milk quality control and is also a useful indicator for mastitis presence.

There are several different methods for determination the total number of somatic cells in milk (SCC). Each method is based on different feature and exploitation characteristics, but none of them is capable of determining the total number of somatic cells. In case of perfect conditions it is possible to determine the SCC of a milk sample, but there are several theoretical and practical problems, that make it impossible, especially when automation and speed of analysis are required. Because of the random positions of the cells in the sample, each result, unless total sample is tested, may only indicate a part of the SCC.

When determining the SCC it is necessary to maintain the requirement of ISO 13366-1:2008 "Milk – Enumeration of somatic cells – Part 1: Microscopic method (Reference method) » standard. Since the cells in milk are located based on the Poason law, there is a determined minimum quantity of cells, which needs to be counted in order to achieve the needed level of exactness (page 6 of the standard). For authentic determination of SCC it is necessary to be counted not less quantity of somatic cells than the mentioned in the following table:

Concentration of somatic cells thousands/ml	Coefficient of variation, CV%	Min. quantity of counted cells
< 150	10	100
150-250	7	200
250-400	6	300
≥ 400	5	400

The Poisson distribution is based on the formula: M = V = s2,

Where:

M – average value – (number of counted cells)

V – Dispersion

s – Standard deviation

The coefficient of variation (CV) would be equal to $CV = \frac{s}{M} * 100\%$

In order to achieve the needed precision, with the different methods of counting, some limitations appear, linked with the volume of the effectively measured milk. This way, if the measured sample contains 100 000 somatic cells in 1 ml, in order to achieve the needed precision (CV<10) there must be counted not less than 100 cells, in no less than 1 μ L undiluted milk. In order to achieve higher precision of measurement – for example CV<5%, it is necessary the amount of measured milk to be substantially greater – to count at least 400 cells that may be contained in 4 μ L milk.

An important moment is the correct identification of the cell, if it occurs, and rejection of each object or particle in the milk, which look like but is not a somatic cell. This is an important factor for determination the precision of each method for counting SCC. A key moment is the possibility of measuring the size of the mentioned cells. The question of selectiveness is further complicated by the fact that most milk samples are analyzed in a place, away from the place of milking, after a period of time. Since the milk is a biological medium with active enzymes and microorganisms, this may lead to a change in the number of cells and their morphology, even if they're preserved.

As a solution to the above, there is general agreement as the reference method for determining the "true" number of somatic cells in milk to accept direct microscopic analysis (DMSCC) using methylene-blue or fluorescent dyes. The method is protected in IDF / ISO standard 148. There are some aspects of this method, though, that limit its feasibility from practical point of view, like long procedure and intensive training of the operator, in order to provide objective selectiveness or precision.

Significant disadvantage of direct microscopic analysis (DMSCC) is the tendency to color artifacts and with the potential problem of cell aggregation and limited sample volume gives rise to uncertainty in the number of cells. Especially clearly this is expressed in the milk with a low content of the somatic cells (up to 300 000 / ml). The accuracy, precision and repeatability of the results of the DMSCC method depend a lot on the operator's training and abilities, notwithstanding the used equipment or protocol.

There are other methods for determining the SCC. Typical SCC methods are counting the SCC by

Similar Standard Method, based on Direct Fluorescent image low magnification microscopic recognition (Lactoscan SCC by Milkotronic Ltd., C-Reader "ADAM" by Digital Bio Technology; NucleoCounter SCC 100 by Chemometec; DCC by DeLaval);

Fluorescent flow cytometry ISO 13366-2:2006 (IDF 148-2: 2006) Milk --Enumeration of somatic cells -- Part 2: Guidance on the operation of fluoro-opto-electronic counters (Somacount 150 by Bentley Instruments; SomascopeTM by Delta Instruments). Based on disc cytometry – fluorine-optical counter with spinning disc (FossoMatic4000TM by Foss Electric). In milk testing laboratories the number of somatic cells is usually determined by automatic electronic machines, which may be exact and reliable instruments.

With them the SCC is quickly and cheaply determined.

There are also indirect methods, where the SCC is determined after the interaction of the milk with certain reagents, for example the use of California mastitis tests. They're with lower accuracy and repeatability.

Milkotronic Ltd. has developed a portable Lactoscan SCC based on direct fluorescent, low magnification microscopic somatic cell counting. Lactoscan SCC uses a very sensitive fluorescent dye Sofia Green, LED optics μ CMOS technologies for capturing in order to make the cells analysis more accurate, reliable and fast. Mostly, the objective selectiveness of Lactoscan SCC is combined with the permanent high stability of the different mechanical, electronic, optical and chemical components of the system, which provides almost identical results during the whole duration of machine's use. On top of that, the production of the machines provides high level of uniformity between the devices. This is a unique feature of Lactoscan SCC which actually offers identical results when we consider the measurement of one and the same sample, notwithstanding which machine is used, what is the location, who is using the machine, and in what moment the analysis is being made, in case that the characteristics of the samples have not changed.

In order to count the somatic cell with Lactoscan SCC, the milk sample is mixed with the dying reagent, containing fluorescent dye Sofia Green. Only 8 μ L from the dyed sample is pipetted on the measuring chamber of disposable LACTOCHIP. The chip is loaded into the device and for a period between few seconds and 2 minutes, depending on the measuring mode, the analysis is done. Lactoscan SCC system focuses automatically on the chip and the dyed cells are captured by the sensitive CMOS camera. The analysis algorithm of digital images determines the number and dimension of the fluorescent cells and counts their concentration. The results are automatically shown on the display, also on printer, with possibility to save the results and generate reports from the results.



© 2017 "Milkotronic" Ltd All rights reserved



<u>Comparison between the direct microscope analysis (DMSCC) and</u> <u>analysis with LACTOSCAN SCC.</u>

A comparison test with measuring of 14 reference samples of raw cow milk treated with glycerin, thimerosal and dimethyl sulfoxide was conducted in the laboratory "Buluritest", Scientific Research Sector, University "Prof. d-r AssenZlatarov"- Burgas, Bulgaria. The total amount of somatic cells in each reference sample of cow milk has been analyzed by the direct microscopic analysis (DMSCC) and by LACTOSCAN SCC.

Direct microscopic analysis, as a control method, was conducted according to ISO 13366-1:2008 «Milk –Enumeration of somatic cells – Part 1: Microscopic method (Reference method)».

10 microliters of the analyzed milk samples were applied on a marked area of 1 cm^2 from the surface of the glass slide obtaining smear, which was dried. After drying, the smears were colored with a solution of methylene blue and propidium iodide fluorescent dye (PI). After that, the colored cells were counted with microscope. The amount of the counted cells in a specific surface (1 cm^2) was multiplied with the working coefficient so the amount of somatic cells in 1 ml to be determined.

The counting was made by counting 50 fields in one smear, while moving the lens of the microscope consequently vertical and right. Therefore from each milk testing sample were analyzed 10 smears, 50 fields were counted in each of them.

Determination of SCC by using **fluorescent-microscopic counter Lactoscan SCC** was conducted according to the company records.

The essence of the method is based on the fact that after demolition of the cytoplasmic membrane of the somatic cells in the milk, under the influence of the used by us dried lysogenic buffer reagent, the cores of the cells become accessible to the action of the used fluorescent dye SOFIA GREEN, which colors the nuclear DNA. The light signals

from the colored cells are detected by the detector of the so called CMOS-camera, i.e. on the elements of the CMOS-camera are formed images, where the emitted energy from every cell is displayed as an illuminated ball. By automatic moving of the mechanism on the axes X and Y, the device captures a maximum of 70 images. The images are analyzed with the help of the embedded software and in this way the amount of the somatic cells is determined. The whole process after placing a LACTOCHIP with a sample in the carriage is automatic.

Performing the analysis

Slowly stir the milk sample while avoiding foaming. The temperature of the sample may be in the range from 10 to 40 °C. With the automatic micropipette put 100 μ L of the stirred milk in an Eppendorf with dried lysis coloring buffer reactive. Slowly stir for a couple of seconds by a multiple pipetting with the automatic micropipette or by using Vortex stirrer. After a stay of about a minute, the content in the Eppendorf is stirred again and 12 μ L of the colored milk are being pipetted in the microfluidic camera of Lactochip x2.



the microfluidic camera can stop and after that put the loaded chip into the carriage of Lactoscan SCC. Now follows an automatic analysis of the sample up to visualizing of the results (expressed in thousands of cells in 1ml) on the display of the device. These results along with all captured images are automatically stored in the data base of the device.

Analysis of the results

The number of somatic cells in the analyzed milk sample is reported by the results, obtained on the display of the device and are expressed in thousands of cells/ml. Repeatability results of every milk testing sample are evaluated by the coefficient of variation (CV), calculated for each of the two methods SCC. The accuracy and repeatability of the system Lactoscan SCC and the average comparison and regression analysis of the SCC data between Lactoscan SCC and the direct microscope analysis (DMSCC) were also established on the base of the obtained results.

<u>Results</u>

1. Accuracy of the system Lactoscan SCC

The accuracy of our method for defining SCC with the help of Lactoscan SCC, based on the technical abilities of the device with its 2 versions for 4 and 9 microliters of measuring is:

Number of cells 1 ml	Lactosca mode microl	n SCC e 4 iters	Lactosca mode microl	an SCC e 9 iters	Minimum requirement of the standard ISO13366-1		
	Counted cells	CV%	Counted cells	CV%	Counted cells	CV%	
100000	400	5,0	900	3,3	100	10,0	
500000	2000	2,2	4500	1,5	400	5,0	
1000000	4000	1,6	9000	1,1	400	5,0	
1500000	6000	1,3	13500	0,9	400	5,0	

Image comparison of the direct microscopic somatic cell counts (DMSCC) and Lactoscan SCC

A) Stained somatic cell image with propidium iodide (DMSCC)

B) Stained somatic cell image with methylene-blue (DMSCC)



Somatic cell concentration: **700x10³/ml**

Magnification: 1:1000



Somatic cell concentration: **1200x10³/ml** Magnification: **1:1000**

C) Stained somatic cells with Sofia Green dye captured by Lactoscan SCC.



Somatic cell concentration: **1200x10³/ml** Magnification: **1:4**

For an analysis with version 4 microliters, the implied accuracy is significantly better than the minimum requirements of the standard. For a chip with a version 9 microliter of the tested milk sample, the result is the maximum possible accuracy for electronic devices.

The measurement results for 14 milk samples are summed up in the following table.

		Lactosca	n SCC		DMSCC						
Milk test samples	Average cell/ml	Standard deviation	Coefficient of variation %	Number of counted cells in 50 fields	Average cell/ml	Standard deviation	Coefficient of variation %	Number of counted cells in 50 field	Necessary visual fields		
1	48 200	2974	6.17	190	50 315	5701	11.33	7.5	667		
2	105 400	5015	4.76	411	108 009	10701	9.91	16.1	311		
3	148 000	6815	4.60	569	165 033	12329	7.47	24.6	407		
4	244 500	9180	3.75	952	258 954	12728	4.92	38.6	389		
5	338 300	12472	3.69	1320	348 178	20846	5.99	51.9	289		
6	453 400	15254	3.36	1768	461 554	30135	6.53	68.8	291		
7	559 000	18708	3.35	2189	554 805	28646	5.16	82.7	242		
8	677 100	20469	3.02	2604	621 891	24301	3.91	92.7	216		
9	776 070	22944	2.96	3019	777 531	37612	4.84	115.9	173		
10	978 900	26917	2.75	3871	989 524	47568	4.81	147.5	136		
11	1 091 000	26541	2.43	4275	1 104 242	43037	3.90	164.6	122		
12	1 140 200	19904	1.75	4502	1 179 379	40102	3.40	175.8	114		
13	1 485 000	24855	1.67	5844	1 463 154	61846	4.23	218.1	92		
14	1 668 900	30700	1.84	6598	1 661 730	53859	3.24	247.7	81		

The linear correlation between the values of SCC in the 14 milk test samples with a progressive increase of SCC from 50 000/ml to 1 650 000/ml (R 2, 0,995; 95% CI, 0,990 \sim 0,999; p <0.01), measured by DMSCC and Lactoscan SSC is presented in Fig. 2

In comparison with the data from SCC for the 14 milk test samples, measured by DMSCC, the system Lactoscan SCC shows an acceptable and similar repeatability and accuracy of the conventional device for SCC in the correlation analysis of Pearson.

At low values of SCC in the milk, it is necessary in the method of DMSCC to be watched much more than 50 visual microscope fields in order to fulfill the requirement of the standard. Like this, for example, for the tested milk with SCC 50 000/ml with the method of DMSCC were counted from 6 to 9 cells in 50 microscope fields, and with the use of Lactoscan SSC - from 171 to 212 cells from a single measurement.
Comparison results between Lactoscan SCC and DMSCC are presented in the graphs below:



Figure1



The common data for the repeatability of SCC in each of the 14 milk test samples with the use of DMSCC are summed up in the following table:

Concentration	1	2	3	4	5	6	7	8	9	10	Arithmeticaverage	Standard deviation	Coefficient of variation %
50 000	53 669	53 669	46 960	46 960	60 378	46 960	40 252	46 960	53 669	53 669	50 315	5701	11,33
100 000	107 338	100 630	100 630	120 756	93 921	120 756	114 047	107 338	120 756	93 921	108 009	10701	9,90
150 000	154 299	181 133	181 133	161 007	161 007	167 716	140 881	161 007	167 716	174 425	165 033	12330	7,47
250 000	261 637	268 346	248 220	261 637	254 928	268 346	268 346	261 637	268 346	228 094	258 954	12729	4,92
350 000	348 849	362 267	368 975	301 889	322 015	362 267	348 849	348 849	355 558	362 267	348 178	20846	5,99
450 000	483 022	462 896	509 857	436 062	456 188	409 227	442 770	442 770	489 731	483 022	461 554	30135	6,53
550 000	529 983	570 234	516 565	523 274	576 943	583 652	536 691	583 652	536 691	590 360	554 805	28646	5,16
650 000	603 778	597 069	597 069	630 612	657 447	630 612	623 904	630 612	590 360	657 447	621 891	24302	3,91
750 000	838 580	724 533	711 116	805 037	764 785	791 620	778 202	791 620	771 494	798 328	777 531	37612	4,84
950 000	1 059 965	979 461	925 792	992 879	919 084	999 587	979 461	1 046 548	1 033 131	959 336	989 524	47569	4,81
1 100 000	1 153 886	1 113 634	1 046 548	1 113 634	1 086 800	1 113 634	1 046 548	1 160 595	1 059 965	1 147 177	1 104 242	43038	3,90
1 150 000	1 113 634	1 160 595	1 153 886	1 160 595	1 133 760	1 227 681	1 214 264	1 194 138	1 220 972	1 214 264	1 179 379	40102	3,40
1 450 000	1 549 696	1 442 358	1 408 814	1 435 649	1 395 397	1 589 948	1 428 940	1 469 192	1 475 901	1 435 649	1 463 154	61847	4,23
1 650 000	1 650 325	1 744 246	1 596 656	1 697 286	1 737 538	1 603 365	1 697 286	1 616 782	1 636 908	1 636 908	1 661 730	53860	3,24

Milk samples	1	2	3	4	5	6	7	8	9	10	Arithmetic average	Standard deviation	Coefficient of variation %
1	43 000	46 000	48 000	46 000	50 000	54 000	50 000	47 000	49 000	49 000	48 100	2974	6,17
2	112 000	101 000	103 000	103 000	112 000	99 000	100 000	111 000	108 000	105 000	105 400	5015	4,75
3	142 000	161 000	147 000	157 000	142 000	145 000	151 000	150 000	140 000	145 000	148 000	6815	4,60
4	251 000	242 000	243 000	230 000	232 000	249 000	242 000	243 000	253 000	260 000	244 500	9180	3,75
5	343 000	335 000	330 000	366 000	322 000	339 000	335 000	328 000	350 000	335 000	338 300	12473	3,69
6	440 000	436 000	474 000	472 000	453 000	446 000	432 000	469 000	450 000	462 000	453 400	15255	3,36
7	571 000	563 000	548 000	530 000	543 000	584 000	558 000	584 000	539 000	570 000	559 000	18709	3,35
8	715 000	675 000	648 000	686 000	682 000	671 000	679 000	683 000	643 000	689 000	677 100	20469	3,02
9	800 700	800 000	802 000	774 000	773 000	763 000	733 000	750 000	789 000	776 000	776 070	22944	2,96
10	940 000	934 000	975 000	1 006 000	987 000	1 007 000	993 000	967 000	970 000	1 010 002	978 900	26918	2,75
11	1 102 000	1 071 000	1 064 000	1 077 000	1 083 000	1 080 000	1 081 000	1 083 000	1 115 000	1 154 000	1 091 000	26541	2,43
12	1 153 000	1 164 000	1 146 000	1 142 000	1 146 000	1 122 000	1 100 000	1 154 000	1 120 000	1 155 000	1 140 200	19904	1,75
13	1 517 000	1 494 000	1 466 000	1 505 000	1 497 000	1 501 000	1 466 000	1 484 000	1 489 000	1 431 000	1 485 000	24855	1.67
14	1 694 000	1 677 000	1 647 000	1 663 000	1 702 000	1 705 000	1 677 000	1 639 000	1 607 000	1 678 000	1 668 900	30701	1,84

The common data for the repeatability of SCC in each of the 14 milk test samples with the use of Lactoscan SCC in measuring of 4 μ L samples are summed in the following table:

As it is seen, the intervals are within those, recommended from the International milk association for somatic counters (to150 000 CV %< 10: from 150 000 to 250 000 CV %<7: from 250 000 to 400 000 CV%<6 and over 400 000 CV%<5). With the use of Lactoscan SCC, based on the microscope counting of somatic cells with a plastic chip for analyzing of the milk, the values of SCC are less variable than with DMSCC. Nevertheless, the obtained ratio of SCC measured by the system Lactoscan SCC, in the milk test samples to those, measured by DMSCC, vary from 0,897 ~ 1,089. Therefore, the values of SCC in the milk test samples, defined with Lactoscan SCC are acceptably repeatable.

Comparison of CV% for the 14 milk test samples for the results, obtained from the both methods of analysis – by the use of Lactoscan SCC and the direct microscope analysis (DMSCC) is presented in the graph below:



Figure 3

Comparison of the potential accuracy of the devices, based on Direct Fluorescent image low magnification microscopic recognition, starting from the declared data from the producers about the effective volume of the measured milk and the results from the Normal distribution of Poisson.

	Lactoscan SCC 4 microlite	C mode ers	Lactoscan SCC 9 microlite	c mode ers	C-Reader "A max 9 micro	DAM" liters	NucleoCounto 100 1 microli	er SCC ter	DELAVAL I 0.1 microl	DCC iters
Number of cells 1 ml										
	Total cells	CV%	Total cells	CV%	Total cells	CV%	Total cells	CV%	Total cells	CV%
50 000	200	7,07	450	4,71	450	4,71	50	14,14	5	44,72
100 000	400	5,00	900	3,33	900	3,33	100	10,00	10	31,62
300 000	1200	2,89	2700	1,92	2700	1,92	300	5,77	30	18,26
600 000	2400	2,04	5400	1,36	5400	1,36	600	4,08	60	12,91
1 000 000	4000	1,58	9000	1,05	9000	1,05	1000	3,16	100	10,00
1 500 000	6000	1,29	13500	0,86	13500	0,86	1500	2,58	150	8,16

Comparison on the base of some of the main features of the devices, based on Direct Fluorescent image low magnification microscopic recognition, starting from the declared data from the producers.

	Lactoscan SCC	DMSCC	<u>C-Reader "ADAM"</u>	NucleoCounter SCC <u>100</u>	DELAVAL DCC
Feature	Milkotronic Ltd.		DigitalBio Technology	Chemometec	DeLaval
Staining Method	Fluorescent dye	Fluorescent/Exclusion dye	Fluorescent dye	Fluorescent dye	Fluorescent dye
Reagent	Sofia Green	EtBr/MB	PI	PI	PI
Microscope	Fluorescent (4 x)	Light (1000 x)	Fluorescent (4 x)	Fluorescent (0,7x)	Fluorescent (4 x)
Sample Volume Needed	100uL	1 mL	100uL	100uL	60 uL
Dilution	no	no	1:2	1:2	no
Measured Volume	9 µL	0,149 μL	9 μL	1 µL	0,1 μL
Captured Images	max 70	~ 50	90	1	1
Machine Analysis Time	max 120 samples/hour	-	max 60 samples/hour	up 100 samples/hour	max 60 samples/hour
Cell Density Range	1 x 10E4 ~ <mark>20 x 10E6</mark>	1 x 10E4 ~ <mark>10 x 10E6</mark>	0 ~ 10x10E6	1 x 10E4 ~ <mark>2 x 10E6</mark>	1 x 10E4 ~ 2x 10E6
Image Processing	Fluorescent image recognition	Human	Fluorescent image recognition	Fluorescent image recognition	Fluorescent image recognition
Calibration	No	No	No	No	No
Applications	Count, Cell size	Count, Cell size, Morphology	Count, Cell size	Only Count	Only Count
Process Time	Fast	Very slow	Fast	Fast	Fast

113 Operational manual

© 2017 "Milkotronic" Ltd All rights reserved

Personal Error	N		Y	Ν	N	N	
Method	Similar Standard Method		Standard Method	Similar Standard Method	Similar Standard Method	Similar Standard Method	
Repeatability	4 µl 9 µl 100 000 CV(%)- cell/ml 5% 400 000 CV(%)- cell/ml 3% 1000 CV(%)- cell/ml 3% 1000 CV(%)- 000cell/ml 2% 1% 1%		<150 000 cell/ml CV(%) - 10% 150 000-250 000 cell/ml CV(%) - 7% 250 000-400 000 cell/ml CV(%) - 6% > 400 000 cell/ml CV(%) - 5%	100 000 cell/ml CV(%) - 3% 400 000 cell/ml CV(%) - 2% 600 000 cell/ml CV(%) - 1%	100 000 cell/ml CV(%) - 10% 400 000 cell/ml CV(%) - 5% 1 000 000cell/ml CV(%) - 3%	100 000 cell/ml CV(%) - 12% 400 000 cell/ml CV(%) - 8% 1 000 000cell/ml CV(%) - 7%	
Sample Prep.	Simple			Complex	Simple	Simple	Simple
Calibration		N		-	Ν	N	N
Cell Types	Sor	matic Cell	Most Cells		Somatic Cell	Somatic Cell	Somatic Cell
Throughput		High		Low	High	High	High
Device Price		Low		-	Low	Low	Low
Peripheral Device	Chip, dry Stain/ Lysis reagent			Microscope, Reagent	Chip, Stain Solution	Cassette, Lysis Buffer	Cassette, dry Stain/ Lysis reagent
Device Size	Small			-	Small	Small	Small
Weight	5 kg			-	10 kg	3 kg	4 kg
Reagent				Manual prep	Stain Solution	Coated in Cassette, Lysis Buffer	Coated in Cassette, Lysis Buffer

It should be noticed that Lactoscan SCC increases the images 4 times, and such magnification is the most suitable and this provides an opportunity for measuring the size of the particles and respectively for more correct differentiation, which has a determining influence for the accuracy of the measurement. For comparison NucleoCounter SCC 100 has only 0.7 magnifications which don't allow a visual differentiation of the cells.

It is necessary to add as well the large volume of measurement = up to 9 μ L, which provides an opportunity for measurements with values of CV, that are significantly better than the standard requirements. The guaranteed high accuracy of measurement with Lactoscan SCC in combination with the competitive price of the device and the lowest on the market price of consumables / the price of consumables was limiting till today the use of this kind of devices / makes it currently extremely attractive.

Comparison on some of the main features of Lactoscan SCC with the devices, based on fluorescent flow cytometry, starting from the declared data of the producers:

	L	actoscan SCC		Somacount 150	SomascopeTM	
		E				
Feature	1	Milkotronic Ltd.		Bentley Instruments	Delta Instruments	
Staining Method	F	luorescent dye		Fluorescent dye	Fluorescent dye	
Reagent		Sofia Green		EtBr/0,083 mg	DAPI	
Microscope	Fl	uorescent (4 x)				
Sample Volume Needed	100uL			3,5 ml	2 ml	
	No			1:2		
Measure Volume	max 9 µL			max 10µL	max 10 μL	
Captured Images	max 70			-		
Machine Analysis Time	max 120 samples/hour			max 150 samples/hour	max 400 samples/hour	
Cell Density Range	1 x 10E4 ~ 20 x 10E6			0 ~ 10x10E6	0 ~ 10x10E6	
Image Processing	Fluorescent image recognition			Fluorescent flow cytometry	Fluorescent flow cytometry	
Calibration		No		Y	Y	
Applications	C	ount, <mark>Cell size</mark>		Count	Count	
Process Time		Fast		Fast	Fast	
Personal Error		N		N		
Method	S	imilar Standard Method		-	<u> </u>	
Repeatability	4 μl 9 μl 100 000 cell/ml CV(%)-5% CV(%)-3% 400 000 cell/ml CV(%)-3% CV(%)-2% 1 000 000cell/ml CV(%)-2% CV(%)-1%			100 000 cell/ml CV(%) - 5% 300 000 cell/ml CV(%) - 3% 500 000 cell/ml CV(%) - 2%	100 000 cell/ml 400 000 cell/ml 1 000 000cell/ml	
Sample Prep.		Simple				
Calibration		Ν		Y	Y	
Cell Types	Somatic Cell			Somatic Cell	Somatic Cell	
Throughput	High			High	High	
Device Price	Low					
Peripheral Device	(Chip, dry Stain/ Lysis reagent				
Device Size		Small				
Weight		5 kg		33 kg	10 kg	



Figure 1

Comparison between Lactoscan SCC and FOSS device based on disc cytometry-fluoride-optoelectronic counter with rotating disc, based on the data declared by the producer.

	La	ctoscan SCC		Fossomatic 5000TM
Feature				
		Milkotronic		FOSS
Repeatability	100 000 cell/ml 400 000 cell/ml 1 000 000cell/ml	4 μl CV(%)-5% CV(%)-3% CV(%)-2%	9 μl CV(%)-3% CV(%)-2% CV(%)-1%	100 000 cell/ml CV(%) - 7% 300 000 cell/ml CV(%) - 5% 500 000 cell/ml CV(%) - 4%



Figure2

<u>Conclusion</u>

A method for defining the number of somatic cells in samples of cow milk with the help of Lactoscan SCC has been developed.

With the help of LACTOSCAN SCC can be AUTHENTICALLY defined the SCC in the raw or preserved with bronopol, boric acid, sodium azide, potassium dichromate and other milk preservatives.

The device can REPRODUCIBLY be used for defining the number of somatic cells both in milk testing laboratories and in dairy farms and milk processors.

Work with Lactoscan SCC is with HIGH REPRESENTATIVENESS (by a single measurement can be defined the SCC in 4 or 9 microliters of milk, while for DMSCC this quantity is 0.149 microliters).

The repeatability indexes for Lactoscan SCC doesn't inferior to other devices with fluorescent microscopy of the other globally presented companies (Foss Electric, Bentley Instruments, and Delta Instruments etc.)

With the device Lactoscan SCC it can be worked easily, simple, fast, accurately and reliably and its portability makes it extremely comfortable to use not only in laboratories but also in the farms with aims to improve the quality of raw milk.

Detection of Yeast mastitis in raw milk samples using LACTOSCAN SCC

The topic of the Yeast mastitis and in general the presence of yeast cells in fresh raw milk is very important for all who work in the dairy sector: dairy farmers, veterinarians, analytical laboratories, scientific researchers, manufacturers of equipment for somatic cell counting and others, which is though barely illuminated in researches and practice.

Its importance is determined by a number of reasons like:

1. Increased and uncontrolled use of antibiotics which stimulates the growth of yeasts in the udder;

2. The presence of yeast mastitis has increased over the last decade;

3. In some countries, more than 10% of the animals suffer from such kind of mastitis.

The main cause of yeast mastitis are: Candida albicans, Candida catenulata , Candida glabrata. Usually these yeasts are also found in milk from animals that have not developed the symptoms of yeast mastitis and at least in the standard ISO 13366-1 IDF 148-1 Milk — Enumeration of somatic cells — Part 1: Microscopic method (Reference method) should be pointed what to be done if such cells come under the microscope. This is important because they are with the size of 5-6 micrometers and should be listed by the size criteria which says that all leukocytes, epithelial cells and cell fragments of a size greater than 4 microns are counted.

"INTERNATIONAL STANDARD ISO 13366-1 IDF 148-1

Milk — Enumeration of somatic cells — Part 1: Microscopic method (Reference method)

8.2 Determination

8.2.1 Reading optimization

Using the microscope (5.3), count the cell nuclei in the obtained smear (8.1.1 or 8.1.2) of fields, entirely filled with milk smear only. Choose the best magnification (from 500x to 1 000x), in order to have an average maximum number of 20 cells in each field.

The cells possess a stained nucleus. The cells generally are 8 μm or larger. Do not count cells less than 4 μm (see Figure 1). Count fragments only if more than 50 % of

nuclear material is visible. Count cell clusters as one, unless the nuclear unit(s) is (are) clearly separated.

See also Figures 2 and 3



Macrophage 8-30 µm

Figure 1 — Examples of cells



PMN 10-14 µm

The relation between is big. Phagocytosis, antigen presentation, secretion chemoattractants





Lymphocyte 5-10 µm

The relation between cytoplasma/nucleus is small. Nucleous intensively stained T helper T suppressor B cell



Epithelial cell 10-14 µm

Nucleus round. Cytoplasma weakly stained

Figure 2 — Examples of cells from bulk cows' milk (1 000 x magnification)



Cell lengths: $L1 = 9,79 \ \mu m$ and $L2 = 2,77 \ \mu m$

Figure 3 — Examples from cells from bulk cows' milk



(500x magnification) (1 000x magnification) Cell length: L22 = 9,08 μ m; L23 = 8,27 μ m; L24 = 4,95 μ m; L25 = 7,39 μ m; L26 = 6,37 μ m and L27 = 3,58 μ m. In the example of a cluster, as shown in Figure 3, five cells have to be counted 1.27

In the example of a cluster, as shown in Figure 3, five cells have to be counted. L27 is omitted because its diameter is less than 4 μ m."

Nothing is said about the existence of yeast cells in milk, which looks similar to the leucocytes and attention has to be paid in order not to be added to the total number of somatic cells. And if in a referent laboratory the laboratory technician, who has never, been taught to recognize the yeast cells from leucocytes, can count and add their number to the total number of somatic cells in reference milk samples, which then to be used for calibration of electronic counters.

The images below show yeast cells dyed with the same reagents used in the ISO 13366-1 IDF 148-1:





Yeast cells dyed in methylene blue in a milk sample



Yeast cells dyed in propidium iodide in a milk sample

It is clear that they can be counted as somatic cells and can affect on accuracy in determining the number of somatic cells in milk from there and on determination of the mastitis' type: healthy animal, sub-clinical bacterial mastitis or clinical bacterial mastitis, from there to affect on the manner and the type of treatment /treatment with antibiotics assist the propagation of the yeast/, on milk price and others.

Until now no one indicated the number of yeast cells per ml milk that does not lead to yeast mastitis as until now IDF has not advised methods /reference DMSCC or devices/ for differential cell count in milk, by which you can determine which cells are with leukocyte origin, which are with epithelial, which with yeast origin and differential counting of leucocytes /determination of neutrophils and macrophages separately as a criterion for the development of mastitis/.

Only in 2016 talks started about such requirements for the reference method and the new reference methods, possessing such features.

In IDF documents:

1. SCSA S15 Improvement of the reference method for SCC:

"Thoughts on a new method

Properties

- Clear identification of cells
- Measurement of the same cells as today's reference method (no change in definition or measurement level)
- Fair repeatability
- Improved performance
- Easy to use for characterization of reference material
- Ideally applicability for differential SCC"

and

On the other hand, the electric methods for somatic cell counting are monopolized by 2-3 companies. Basically they are using flow-cytometry in its oldest version, compared with the medical flow-cytometry, without determination of cell size and their differentiation as it is in medical fluorescent flow-cytometers. In practice, they can measure the yeast cells, but manufacturers do not deepen in the problem. They set apparatus not "see" most of the yeast cells by setting threshold /level of the separation of the signals from the fluorescent cells is higher than the level which emit yeast cells/, as they are small than 4-6 μ m, and possess less DNA. Because these devices are like big black boxes, the operators have no visual control as monitoring under a microscope and the presence of yeast goes unnoticed.

In LACTOSCAN SCC, for the first time there is a possibility to automatically count the number of yeast cells in milk samples when determining the number of total somatic cells with reference accuracy.

Below there are images, downloaded from LACTOSCAN SCC devices, located in different parts of the world. When in the milk there are only leukocytes, induced by bacterial mastitis or only a normal level of leukocytes / macrophages, neutrophils and lymphocytes / the average cell size is in the range 9-14 μ m.



Image 1. Milk sample from a healthy animal



Graphic from Image 1. Minimum cell size set to 3 µm

When the cell size range is set from 3 to 20 μ m, the counted cells are 4, equivalent to 41×10^3 . The cells size is between 13 and 14 μ m. The average cell size is 13.75 μ m and there are no small cells with size under 9 μ m. The cell size peak is at 14 μ m.



Graphic from Image 1. Minimum cell size set to 8 µm

When the cell size range is set from 8 to 20 μ m, the counted cells are 4, equivalent to 41×10^3 . The cells size is between 13 and 14 μ m. The average cell size is 13.75 μ m and there are no small cells with size under 9 μ m. The cell size peak is at 14 μ m. According to the graphics, it can be concluded that there are only somatic cells in the milk sample. There is no change in the final result when the cell size range of the counted cells is set to from 3 to 20 μ m and when set from 8 to 20 μ m.





Graphic from Image 2. Minimum cell size set to 3 µm



When the cell size range is set from 3 to 20 μ m, the counted cells are 497, equivalent to 5021×10^3 . The cells size is between 4 and 15 μ m. The average cell size is 7.48 μ m. Most of the cells are with size under 9 μ m. The cell size peak is at 7 μ m.



Graphic from Image 2. Minimum cell size set to 8 µm

When the cell size range is set from 8 to 20 μ m, the counted cells are 186, equivalent to 1879×10^3 cell/mL. The size of the cells which are included in the final result is between 8 and 15 μ m. The average cell size is 9.40 μ m. The cell size peak is at 8 μ m.

If we compare the results from both graphics, we can conclude that on the image there are 67% yeast cells with size between $4 - 8 \mu m$ and 37% somatic cells. The higher number of yeast cells, 311 counted cells = 3142 x10³ cell/mL, is a confirmation for presence of yeast mastitis.

The number of somatic cells, 186 counted cells = 1879×10^3 cell/mL, is a confirmation of presence of clinical mastitis.

According to the results, we can conclude that the analyzed milk sample is from an animal with yeast mastitis and clinical mastitis at the same time.



Image 3. Milk sample from an animal with yeast mastitis





When the cell size range is set from 3 to 20 μ m, the counted cells are 238, equivalent to 2074x10³ cell/mL. The cells size is between 5 and 13 μ m. The average cell size is 6.88 μ m. Most of the cells are with size under 9 μ m. The cell size peak is at 7 μ m.



Graphic from Image 3. Minimum cell size set to 8 µm

When the cell size range is set from 8 to 20 μ m, the counted cells are 43, equivalent to 375x10³ cell/mL. The size of cells which are included in the final result is between 8 and 15 μ m. The average cell size is 8.56 μ m. Most of the cells are with size 8-9 μ m. The cell size peak is at 8 μ m.

If we compare the results from both graphics, we can conclude that on the image there are 82% yeast cells with size between 5 – 8 μm and 18% somatic cells. The higher number of yeast cells, 195 counted cells = 1699x10³ cell/mL , is a confirmation for presence of yeast mastitis.

The number of somatic cells, 43 counted cells = 375×10^3 cell/mL, is a signal for possible development of sub-clinical mastitis.

According to the results, we can conclude that the analyzed milk sample is from an animal with yeast mastitis and possible development of sub-clinical mastitis at the same time.



Image 4. Milk sample from an animal with yeast mastitis





When the cell size range is set from 3 to 20 μm , the counted cells are 2221, equivalent to $18679 x 10^3$ cell/mL. The cells size is between 5 and 14 μm . The average cell size is 7.60 μm . Most of the cells are with size 6-8 μm . The cell size peak is at 7 μm .



Graphic from Image 4. Minimum cell size set to 8 µm

When the cell size range is set from 8 to 20 μ m, the counted cells are 942, equivalent to 8137×10^3 cell/mL. The size of the cells which are included in the final result is between 8 and 14 μ m. The average cell size is 8.94 μ m. Most of the cells are with size 8-10 μ m. The cell size peak is at 8 μ m.

If we compare the results from both graphics, we can conclude that on the image there are 66% yeast cells with size between 5 – 8 μ m and 44% somatic cells. The higher number of yeast cells, 1279 counted cells = 10542 x10³ cell/mL, is a confirmation for presence of yeast mastitis.

The number of somatic cells, 942 counted cells = 8137×10^3 cell/mL, is a confirmation for presence of clinical mastitis.

According to the results, we can conclude that the analyzed milk sample is from an animal with yeast mastitis and clinical mastitis at the same time.



Image 5. Milk sample from an animal with yeast mastitis





When the cell size range is set from 3 to 20 μ m, the counted cells are 141, equivalent to 1213×10^3 cell/mL. The cells size is between 6 and 17 μ m. The average cell size is 8.89 μ m. Most of the cells are with size 6-11 μ m. The cell size peak is at 7 μ m.



Graphic from Image 5. Minimum cell size set to 8 µm

When the cell size range is set from 8 to 20 μ m, the counted cells are 77, equivalent to 662x10³ cell/mL. The size of cells which are included in the final result is between 8 and 17 μ m. The average cell size is 10.78 μ m. Most of the cells are with size 8-9 μ m. The cell size peak is at 9 μ m.

If we compare the results from both graphics, we can conclude that on the image there are 55% yeast cells with size between 6 – 8 μ m and 45% somatic cells. The higher number of yeast cells, 64 counted cells = 551 x10³ cell/mL, is a confirmation for presence of yeast mastitis.

The number of somatic cells, 77 counted cells = 662×10^3 cell/mL, is a confirmation for presence of clinical mastitis.

According to the results, we can conclude that the analyzed milk sample is from an animal with yeast mastitis and clinical mastitis at the same time.

In the results that we received from customers in Columbia, we found a lot of animals from one herd with average cell size 7-8 μm . We did researches of such animals and the researches emphatically have shown that the increased number of small cells /4-7 μm / always is due to the existence of yeasts in milk. The yeasts have such cell size.

In researches for India and Brazil it is said that up to 12% of the animals have yeast mastitis. Probably the percentage for the other part of the world is close to 12% in some areas.

The big problem is that if these cells are counted as it is required in the standard for somatic cell counting /all cells with size above 4 μ m are counted for somatic cells/, though the boarder of 500 000 cells indicating the existence of bacterial mastitis will surpass many times not due to somatic cells, but mainly due to the counted yeast cells, which proliferate very fat in milk and even the antibiotic treatment helps for this proliferation.

Very easily the vet can err for such milk samples and to accept that they have to treat mastitis for bacterial mastitis but it could have been yeast mastitis. For this reason, LACTOSCAN SCC is set the range for somatic cell counting from 8 to 20 μ m for cow milk. This way cells with cell size 4-6-7 μ m are divided from somatic cells. Of course, this way cells, mainly small lymphocytes with similar cell size, are not included in the final result but generally they are only 1% from the total number of lymphocytes and body cells entered in milk.

Below there are photos of Petri dishes with such milk sample (average cell size is 6-7-8 μ m), cultured in a special culture medium for yeast, showing the proliferation of yeast in them. The microbiological cultures and the results are clearly visible from the photos:



Photo 1. Milk from a healthy animal

Photo 2. Milk from an animal with yeast infection



Photo 3. Milk from animal with yeast infection



Yeast colony

Photo 4. Milk from an animal with yeast infection



Yeast colony

II. ULTRASONIC MILK ANALYZER

If you are in the SCC mode, by pressing the icon in the red rectangle with picture of the device will switch to the milk analyser mode:

E	Main	🖺 Graph	Options	Advanced					<u>8</u> V 08-	10 17:04
		Νο	Image				Start			
			Ū		-	ID Client Milk type Average cells per mL				
	ID	Chip	Date	Client	Milk	Dilution	n Cells p	er mL	Avg. cell size [µm]	
Empty										Details Summary
	actoscan.com									2016-08-10 17:04

Specification of LACTOSCAN COMBO's Ultrasonic milk analyzer

Switching Adapter

• Input:

100-240 V ~1.6 Amax. 50-60 Hz

Output: +12 V---3 A min.
Output power: 36-42 W

Measurement modes

- cow milk
- sheep milk
- UHT milk
- goat milk
- buffalo milk
- cream 25%
- cream 45%
- whey
- recovered milk
- other /pasteurized milk/

In	the	table	below	the	configuration	of	the	milk	analyzer	including	spares	and
aco	cesso	ories, v	vhich ar	e de	livered on cust	om	ers r	eques	t are liste	d:		

Nō	Description	Item Nº	pcs	
	a) included in the set: $igtimes$			/
	b) not included in the set (may be additionally bought):			
1.	1 sample measurement time	60 sec		
		30 sec		
2.	RS232 Interface Cable - Analyser-IBM PC	LSS006		
3.	pH measuring system	LSS009	1	
4.	pH probe with cable and holder	LSS010	1	
5.	Buffer solution Ph 60 ml (pH7.00±0.01/20°C)	LSS011	1	
6.	Buffer solution pH 60 ml (pH4.00±0.01/20°C)	LSS012	1	
7.	Milk conductivity measuring system	LSS013	1	

8.	Buffer solution conductivity 50 ml (5.02 (±5%) mS/cm (18±0.1°C)	LSS014	1	
9.	RS232 Interface Cable - Milk Analyser – Serial Printer/IBM PC	LSS018	1	
10.	High-fat measurement function	LSS020	1	\square
11.	Spare O-ring for the pH probe		1	

1. FUNCTION

The function of the milk analyzer is:

- to make quick analysis of milk on fat (FAT), non-fat solids (SNF), proteins, lactose and water content percentages, temperature (°C), freezing point, salts, total solids, as well as density of one and the same sample directly after milking, at collecting and during processing.
- Total integration of the milk sample's measured parameters with the processes for storing, processing and sending the results towards the corresponding receiver.

2. TECHNICAL PARAMETERS

2.1. Working modes characteristics:

The program of the milk analyzer has four working modes.

2.1.1. Measurement mode milk / dairy product – first type

2.1.2. Measurement mode milk / dairy product – second type

2.1.3. Measurement mode milk / dairy product – third type

These modes have been calibrated on customers' request for 3 milk types from the following: cow, sheep, UHT, buffalo, goat, camel milk, cream, ice cream mixtures, whey, recovered milk, etc. before leaving the production facilities and the text on the display will be for the corresponding types.

2.1.4. Cleaning

2.2. Measuring range:

Fat	from 0.01% to 25%
SNF	from 3% to 15%
Density **	from 1000 to 1160 kg/m ³
Proteins	from 2% to 7%
Lactose	from 0.01 % to 6 %
Water content	from 0 % to 70 %
Temperature of milk	from 1°C to 40°C
Freezing point***	from – 0,4 to – 0,7°C
Salts	from 0,4 to 1,5%
PH*	from 0 to 14
Conductivity *	from 3 to 14 [mS/cm]
Total Solids*	from 0 to 50 %

* Option, on customers' request

** Density data are shown in an abbreviated form. For example 27.3 have to be understood as 1027.3 kg/m³. To determine the milk density, write down the result from the display and add 1000.

Example: result 21,20; density = $1000 + 21,20 = 1021,2 \text{ kg/m}^3$

The abbreviated form of the density is used also when entering data for samples in working mode **Recalibrate**, for example:

If the measured sample density is 1034.5 kg/m3, then in the menu for entering the samples parameters used for calibration, across the parameter Den =, you have to enter 34.5.

**** Please, carefully read Appendix Freezing Point.

2.3. Accuracy:

Fat	± 0.10%
SNF	± 0.15%
Density	\pm 0.3 kg/m ³
Proteins	± 0.15%
Lactose	± 0.20%
Water content	± 3.0%
Temperature of milk	± 1°C
Freezing point	± 0.001°C
Salts	± 0.05%
РН	±0.05%
Conductivity	±0.05
Total solids	± 0.17%

The difference between two consequent measurements of one and the same milk could not exceed the maximum permissible absolute error.

2.4 Correct ambient conditions:

Maximum permissible absolute error is guaranteed in case of normal ambient conditions:

Air temperature	from 10°C to 40°C
Relative humidity	from 30% to 80%
Power supply	220V (110V)
Extent of contamination at normal environme	ntal conditions2



Maximum permissible absolute error values in point 2.3 are in dependence on the correctness of the corresponding chemical method, used for component content determination. In point 2.3. are used the following reference methods: Gerber - for fat, gravimetric – for SNF, Kjeldahl – for protein. The boundary for maximum variation of repeatability when the power supply voltage is from +10 to -15% from the nominal voltage values (220 V) have to be no more than 0.8 accuracy according point 2.3. The analyzer is used in conditions free of outer electrical and magnetic fields (except the magnetic field of the Earth) and vibrations.

2.5. Continuous working time:

.....non-stop

2.6 Milk sample volume per one measurement:

All rightsreserved

3. Qualification of raw milk, thermally treated milk, other dairy products and derivatives

3.1. Taking samples and preparation for analysis

In order to receive reliable results in qualification of milk, dairy products and derivatives are needed: precise samples taking; correct samples storing (in need to be preserved); correct preparation before making measurement. The rules and requirements for this are described in details in *Appendix Preparing Samples*.

3.2. Making the measurement.

3.2.1. Preparing the analyzer for working mode

3.2.1.1. Put the analyzer on the working place, providing good ventilation and not in the vicinity of heat providing devices or sources. The temperature in the premises has to be in the boundaries $10-30^{\circ}$ C.

3.2.1.2. Switching on the analyser is done with a switch at the rear panel of the analyser. The tablet, controlling the analyser is powered by its power switch. Wait for Windows to be loaded. The operator has an access to all the resources of the system, where can be set, if necessary - printer (for example standard A4 printer), WiFi or connection to other peripheral devices (for example mouse and keypad, flash memory etc.). There's an USB Hub, which ports are mounted on the back panel of the analyser. These are standard USB ports, connected to the tablet (controlling computer).

There are 2 ways of control of the tablet/software:

A. By direct touch on the buttons on the display (using the possibility of the touch screen of the tablet) and buttons on the front panel (if the version is with such buttons);

or

B. With the mouse and/or the keyboard, working like a standard computer.

Version A is suitable for daily work for measuring and cleaning the analyser, while version B is recommended for setting and servicing.

The software can be started directly after Windows has been loaded or by pressing an icon, in dependence of the custom's requirement. After the software is loaded, the main menu appears on the display:


Description of the main screen:

Position 1. System menu – the upper one on the display. It is described with the file Supplier Tools LS-Menu Help. Serves for direct control of the settings/test/service mode of the analyser.

Position 2. System Status – includes groups Supplier, Liters and information tab, imitating the display of the standard (analyser without a tablet) Lactoscan (upper right corner). Serves for entering the ID of the deliverer, liters/ kgs which were delivered and shows the current information, describing the exacution of the specific commands sent from the tablet towards the PCB of the analyser.

Position 3. Tab for measurement results – Fat=00.00, SNF=... Information tab and could not be changed by the operator. It serves for displaying the measurement results. The last measurement results stay active till the next measurement is started.

Position 4. Control buttons. The buttons at both ends are with fixed names – Print and Clean – for starting the corresponding commands. The rest 3 are with changeable names, because they show the name of the corresponding calibration (type of milk). By pressing them, a measurement on the chosen calibration is started. If the analyser is with outer keyboard, mounted in the box of the device (a number of buttons with the corresponding designations), by pressing the buttons left/right is moving the choice of the calibration. Using this keyboard, the chosen calibration is started by pressing the button Enter.

Position 5. Archive tab, type List Box at the lower part of the display. It shows the last 10 measurement results. The meaning of the results is described in the list over them.

Position 6.Information line. It is at the bottom of the display. The group Save File shows the name of the file, where the results will be saved. It is generated automatically, based on the current date. If needed, it ca be changed. The group Meas Count – shows the number of the measurements already done. At the right of this line is displayed the version of the PCB (Main PCB – on which the measurement system is based) and software in it.

3.2.2. Working with the analyser in measurement mode.

This is the main working mode. It is needed a milk to be poured in the glass, to place it in the recess of the analyser, to choose the needed calibration and start it. The sample is sucked, display is changed as follows:



The calibration used at the moment remains on the display. The previous results are hidden. First appears the temperature of the sample. A clock/timer appears in the middle of the display, showing in seconds the remaining time till the appearance of the measurement results. After measurement completion, the results appear in the tab for the measurement results as well as in the archive tab i.e. the initial start display with valid results appears. The results are archived to a text file and can be used for later result analysis. If there is a internet connection(e.g. WiFi, 3G modem, etc.) the results can be send to a cloud, where they can be stored or further analysed



If there's serial printer connected (loaded with paper, with closed cover), the results are printed out. The analyser is ready for the next command (it is in Idle mode) i.e. ready for the next measurement. Button Print is active – by pressing it a copy of the results can be printed out (it is advisable one printout to remain in the milk collecting center, while the second one to be given to the deliverer).

4. CLEANING THE ANALYZER

This procedure prevents gathering milk fat residues and milk stone on the sensor. The milk stone consists of milk solids, calcium, iron, sulphates, magnesium, etc. All these substances form layer on the pipe and sensor's walls, which leads to deviations in the measurement results and blocking up the piping.



The company-producer recommends usage of the chemicals, supplied with the analyser – alkaline and acidic (Lactodaily and Lactoweekly). You may order them separately or together with the analyser. Try to use only these chemicals for cleaning the analyser.

In case you missed to order these chemicals, the alternative is to use alkaline and acidic cleaning solutions for dairy equipment by one the companies, producing such chemicals, as for example:

http://www.diversey.com

http://www.ecolab.com

http://www.calvatis.com



Do not use chemicals not intended for usage in the milking systems or vessels in the dairy sector. Pay special attention to the concentration of the acidic chemical. **Increased concentration may damage the measuring sensor.**

4.1. Periodically cleaning (rinsing) the analyzer

It is done in the process of routine work of the analyzer. Its aim is to prevent drying

up and adhesion of different milk components in the milk analyzer's measuring

system.

4.1.1. Periodical cleaning frequency.

It is easy to understand what is the period on which the rinsing could be done as the analyzer reminds you when it is necessary. This is done by a sound signal in 1-second cycle after the set time intervals elapse:

• 55 min. after switching on the power supply of the analyser, but idle work;

• 15 min. after the last measurement of real milk sample.

*Idle Mode is that part of the standard working mode, when the analyser is not making measurements. There's embedded in the analyser system for measurement of the idle time. The idle time is measured starting from the last action of the operator. In dependence of it (what the operator last did), are taken decisions regarding the cleaning.

There are 2 options:

Option A: If the analyser:

- 1. Was only switched on but was not started in measurement mode,
- 2. Or the last action was cleaning,
- 3. Or the last action was measuring sample with very low Fat (similar to water)

Then the signal for cleaning is started after 55 min.

Option B: If the last thing done with the analyser was measurement of normal milk sample, the signal for cleaning is started after 15 min.

After cleaning completion, new measurement takes place in above described time intervals.

4.1.2. Making the rinsing

After above message is received put in the recess of the analyzer a sample holder with alkaline cleaning solution or water.

Press Clean to start the rinsing mode.

In this mode the analyzer makes 8 cycles and stops.

Already used solution is poured out of the analyser. Now the device is ready for the next measurement. In case of doubt that the analyzer is still not well cleaned, the procedure Cleaning may be executed repeatedly.

4.2. Complete cleaning

4.2.1. Complete cleaning frequency

This cleaning is done after finishing the work with the analyzer at the end of the working day or if it is obvious that the measuring system of the analyzer is contaminated in case of intensive work with it. It is done with alkaline cleaning solution.

Preparation of 3 % alkaline solution of Lactodaily for circulation cleaning in the milk analyzer:

- 1. Take the package 100 g concentrated chemical Lactodaily
- 2. In appropriate vessel (for example bucket) pour 1 l water.

3. Add the powder and then again water up to 3 l.



For a single cleaning cycle you need only 25 ml cleaning solution. We recommend you to prepare working solutions of cleaning chemicals, enough for normal work for 1 week, because, during their stay unused, the working solutions loose their strength and also is difficult to store them.

Then follow the instruction for milk analyzer cleaning.

4.2.2. Cleaning

4.2.2.1.Rinsing milk residues

Fill in the glass with water. Put it in the recess of the analyser and start command Cleaning from the main menu. After finishing it pour out the contaminated water.

4.2.2.2.Cleaning with alkaline cleaning solution

Fill in the glass with warm (50-60 C) alkaline cleaning solution. Put it in the recess of the analyser and start the command Clean. After finishing it, pour out the contaminated liquid.

4.2.2.3.Rinsing with water

Fill in the glass with water. Put it in the recess of the analyser and start command Clean. After finishing it pour out the contaminated water. Now the device is ready for work.

4.2.2.4. Cleaning with acidic solution

It is recommended to be done every day.

Preparation of 3 % acidic solution of Lactoweekly for circulation cleaning in the milk analyzer:

- 1. Take the package 100 g concentrated chemical Lactodaily
- 2. In appropriate vessel (for example bucket) pour 1 l water.
- 3. Add the chemical and then again water up to 3 l.

Labels for the cleaning chemicals



The following procedure is executed:

1. Rinsing the milk residues:

Fill in the glass with water. Put it in the recess of the analyser and start command Clean. After finishing it pour out the contaminated water.

2. Cleaning with acidic solution

Fill in the glass with warm (50-60 C) acidic cleaning solution. Put it in the recess of the analyser and start the command Clean. After finishing it, pour out the contaminated liquid.

3. Rinsing with water

Fill in the glass with water. Put it in the recess of the analyser and start command Clean. After finishing it pour out the contaminated water. Now the device is ready for work.



Please, pay attention that, when the analysers gives a signal for need of cleaning 15 min after the last measurement of real milk samples or 55 min. after being powered and not used, cleaning is made ONLY with alkaline solution in concentration 1-3%.

During the basic/final cleaning consequence is: alkaline solution – water – acidic solution - water



By pressing the button Clean, a cleaning procedure is started.



After the cleaning cycles are finished, the analyser returns "Lactoscan Cleaned" message and "Ready to start"



5. Menus' description.



These functions are NOT for everyday use. They are used in case some parameters of the analyser to be changed, after careful acquaintance with these commands and how to use them

Service mode menus are selected from the main menu system at the top of the screen. For easier work with the service menu we recommend using a mouse and a keyboard. Structure of the menu system:

5.1 Menu File

System Setup

Purpose: To select the communication ports. Used in production conditions, or if necessary in the service conditions

Reset System

Purpose: To reset the hardware of the device - the specialized board. It is equal to switching the power Off / On.

Exit

Purpose: Exit from the software, controlling the analyser and returning to windows desktop. It is used in cases when the tablet needs settings to be madefor example connecting to a WiFi network, connecting an external printer, etc.

ShutDown

Purpose: Exitfrom the software and shutting down the tablet. It is used to turn off the system. After pressing shutdown, turn off the pcb using the switch on the back panel.

5.2 Menu Supplier

Add Supplier

Purpose: To enter the details of a new supplier of milk. It is used in case of building or using a report system.

Edit Supplier

Purpose: To edit / change data for already introduced suppliers of milk.

5.3 Menu Tools

LSCal

Purpose: To start an external software tool, used for the calibration of the device. The way of work with the tool LSCal.exe is described separately in the service documentation. After completion of the calibration, to resume work with the program press the button "Finish".

Calibration

The program LSCal.exe is used to calibrate Lactoscan Milk Analyzers- WLS Version. Two samples are used – with high and low fat, the third measurement group is measurement with water. Measurements are taken in the following order: high fat sample, low fat sample, water.

Note:

Before starting the calibration procedure, the operator has to be acquainted with the analyzer's manual, mainly in its part preparation of samples, described in service manuals.

Program control:

After the program is started the display shows the following:

🕶 Untitled - LSCal	
<u>File Edit View H</u> elp	
Calibration sample - parameter values 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	COM Port © 1 © 2 © 3 © 4 Start Measure Calibration © 1 © 2 © 3 Download
Ready	Meas: 0 / 0 //

Radio button:

COM Port

- 1 COM 1 connection
- 2 COM 2 connection
- 3 COM 3 connection
- 4 COM 4 connection

Calibration

- 1 Calibration cal 1
- 2 Calibration cal 2
- 5 Calibration cal 3

Control buttons:

for work	Start	- Starts the communication, this means the device is ready
	Measure	- starts a measurement in Calibration mode.
	Download	- saving the new calibration in the device

Edit boxes:

Group Calibration sample – parameter values - set the sample parameters values for each parameter at a time.

Colomns N1, N2, Temp - the results of the measurements in calibration mode, used for making the calculations of the new calibration.

Workflow

This program receives the control from the main program by selecting the Tools menu and choose LSCal.



It opens the software tool for calibration. The operator chooses the Com Port, and presses the button Start. After that inputs the parameter values in the fields. With the button Measure starts a new measurement of the sample in calibration mode.

Details of the calibration process are described in the service documentation.

After finishing the calibration, close the LSCal program and press button Finish of the main program to return to normal use of the analyser.

5.4 Menu LS-MENU

Special Modes

Cycle

Purpose: Starts the unit in cyclic measurement of the ultrasound. The results appear in the results field, in place of Fat, SNF ... The results are remembered in a file. The mode is used for industrial purposes or for service work by qualified professionals. This type of measurement is not used in everyday laboratory work.

	Supplier Number:	Kilos:		Getting Ready Ready to start Mode Cycle
	51 co	ТЕМО	09.7	
	51 se		20.7	
		Stop		
N1 N2 Temp	Ser N Ampl	11		
Save File:	20160810n1.txt	Measure Count:	Reset	Version: 59 5/8/16

5.5 Menu Correction

Temperature Correction

Purpose: To correct the measured sample temperature. Used by trained service specialists. The measured temperature can be varied within + - 9.9 degrees.



Correction

All Corrections

Purpose: Provides information about all existing corrections of the system parameters. The window is informative, for editing adjustments proceed to the selection of calibration and selection of parameters to be adjusted (see below)

Sections		
Calibration 1	Calibration 2	Calibration 3
Fat:-0.95	Fat:0.29	Fat:0.13
SNF:-4.75	SNF:0.35	SNF:-0.75
Density:-4.75	Density:-0.45	Density:-1.15
Lactose:-0.95	Lactose:-0.36	Lactose:0.67
Salts:0.95	Salts:0.04	Salts:0.33
Protein:0.95	Protein:-0.06	Protein:-0.16
Water:-9	Water:3	Water:-3

Corrections

Calibration

Parameter:

Purpose: The menus are used for selecting which parameter of which calibration will be corrected.



5.6 Menu Settings

Set Calibration Names

Purpose: To edit the names of the calibrations. The operator can introduce a free string of 8 alphanumeric characters. After starting, the screen displays the following, the operator must follow the instructions of the program:

Calibration 1	Cow
Calibration 2	Buffalo
Calibration 3	Mix Milk
Cancel	Set

Using the keyboard, the operator inputs the names of calibrations. Memorize them by pressing the Set button or refuse corrections by pressing button Cancel.

Save Calibration

Purpose: To save the calibration to a file in the computer. Select the calibration to be saved



Restore Calibration

Purpose: Restore an already saved to the tablet calibration

*Note: Each calibration is specific to the analyser. You can save a calibration and restore it only to the same machine

Autoprint Enable/Disable

Purpose: To select whether the device will automatically print the results after finishing the measurement



Format Results:

Purpose: To select the format of the results – main or large. The device comes with factory settings for Large results. The machine should be set for operating with Large Results

High Fat Control

Purpose: To enable/ disable high fat measurement for a specific calibration for machines with High Fat Option. When selected from the menu, choose a calibration:

File	Tools	LS-Menu	Measurem	ient He	lp							
		S	upplier Nur	nber:			Kilos:		Ready Pleas for hig	y to star e select jh fat	t calibrat	ion 🖡
		Calibr	ation 1 Hig	h Fat:		Enable	d					
	Cancel		Enable					Disab	ole			
Time	Sup N	N Kg	Cal Tem	p Fat	SNF	Den	Pro	Lac A	d W 🛛 Fr F	Cond	pН	Salts
•						m						•
								Reset	Versio	on: 59 5	/8/16	

High Fat Speed For Calibration

Purpose: For machines **with** High Fat Option. Turn on High Fat speed for calibration, before attempting a calibration with high fat samples

		S	Supplier	Numb	er:			Kilos:		R	eady t	o start		*
			Turn	On					Turn	OFF				
Time	Sup N	Kg	Cal	Temp	Fat	SNF	Den	Pro	Lac	Ad W	Fr P	Cond	рН	Salts
•						T	1		Reset	Ve	ersion	: 59 5/8	8/16	Þ

5.7 Menu Tests

Test Pump

Purpose: To test of the pump unit. Used in service conditions by qualified specialists. To end the test press the button "Finish".

Ultrasound

Purpose: To set the measuring system of the device. Used in service conditions by qualified specialists. The procedure is similar to that used in standard devices without integrated tablet. Current messages are displayed in the

field	System	Status	(top	left	corner	of	the	display)
Test	US - Pha	se: 1		÷				
					-			
						N	lovt	
							ICAL	
						XXX		

To finish the procedure, press the button "Finish".

Set Amplitude

Purpose: To set the metering system of the device. Used in service conditions by qualified specialists. The procedure is similar to that used in standard devices without integrated tablet. Current messages are displayed in the field System Status (top left corner of the display). When the procedure is finished apparatus brings a message that the amplitude is found and asks the user to set this amplitude or cancel the procedure

Amplitude Found Set or Cancel Ready to Start	Supplier	Litres
	Cancel	Set

5.8 Menu Cleaning

Starts cleaning of the analyser. It is equal to pressing the button Clean from the main menu.

5.9 Menu pH & Co Meter

pH Meter En/Dis

Purpose: To turn on / off pH measurement in the process of measuring other parameters. It is available for devices with manufactured hardware option pH.



pH Test

Purpose: To set a pH measurement system of the device (if the device is manufactured with an optional pH). Used in service conditions by qualified specialists.

pH Measure

Purpose: For off-line measurement of pH. In this case the analyser works only as pH meter, without measuring other parameters of the sample.

pH Meter Calibration

Purpose: To calibrate the system for measuring the pH. After starting, the screen displays the following, the operator must follow the instructions of the program:



The operator inputs the values of the buffers used for calibration. Selection of buffers and principles of work with the pH meter are described in the document WLS-TechManual. By pressing the Set button to go to the actual calibration. Reports of this process is displayed in the System Status, in the upper right corner of the screen. By pressing the Cancel button procedure is terminated and transferred to the main screen.

Co Meter En/Dis

Purpose: To turn on / off the conductivity measurement in the process of measuring other parameters. It is available for devices with manufactured hardware option conductivity.



Co Meter Test

Purpose: To set up the measuring system to the conductivity of the device (if the device is manufactured with an option conductivity). Used in service conditions by qualified specialists.

Co Meter Calibration

Purpose: To calibrate the system for measuring conductivity. After starting, the screen displays the following, the operator must follow the instructions of the program:



The operator enters the values of the buffer that will be used for calibration. Selection of buffer and principles of working with the system for measuring conductivity are described in the document WLS-TechManual. Pressing the Enter key passes to the actual calibration of the measurement system. Reports of this process is displayed in the System Status, in the upper right corner of the screen. By pressing the Cancel button procedure is terminated and transferred to the main screen.

Sample calibration

Purpose: To make a precise calibration of the Conductivity Meter. For the calibration are used 5 conductivity buffers

The procedure is as follows:

Enter first sample value First Sample Value: Cancel Enter

-	Put a Sample and press start		_		
		Littes	-	Put new sample and	ш
				press Start	
	First Sample Value:	2			
		otart			

- Repeat the same steps for the rest of the buffers. For each sample five measurements are made.

5.9 Menu Pump Control

Purpose: The following menus are used for test of the pump unit. Used in service conditions by qualified specialists.

Pump In

Pump Out

Stop

5.11 Menu Measurement

Discard Last measurement

Purpose: Discard the last measurement from the daily report

5.12 Menu Help

LS Identity

Purpose: Displays information about the hardware of the device - serial number version of the program board.



About

Purpose: Displays information about the program in the tablet control hardware appliance.

Admin Mode

Purpose: Enter admin password for enabling admin menu. For receiving your admin password please contact your Milkotronic representative



5.13 Menu Admin Mode

Restore Factory Settings

Purpose: Return the machine in the state of which it was produced.

NOTE! This procedure will erase all user data such as calibration corrections, new calibrations and will restore the machine to the state in which it was received from Milkotronic.

6. Peristaltic pump service

Peristaltic pump



7. Possible malfunctions and error messages, troubleshooting

In the table below are described the possible malfunctions during the milk analyzer's exploitation and ways for their repair/remedy. If the problem persists after all recommended measures are taken, please, connect the nearest service center for help. Do not forget to tell the analyser's identity.



To receive the analyzer's identity, refer to point 3.2.1.3.

Error message	Possible problem /cause	Repair/remedy					
2 MA		Immediately switch off the analyzer.					
overheated	Overheated	Pay attention the analyzer to be situated away from direct sunlight or heating devices.					
a continuous sound signal	milk analyzer	Wait 5-10 minutes the device to cool down or to be normalized the ambient temperature and switch i on again.					
3 Empty Camera	Insufficient quantity of the milk sample sucked in the system or air in the sample	 The analyzer is ready to measure the next sample. In order to avoid the future appearance of the same error message, please, check the following: The sample is prepared according the instructions and there aren't air bubbles in it. There is a real suction of the sample after starting measurement, i.e. it is obvious that the level of the milk sample in the sample holder decreases. In other case – there is damage in the suction system. Avoid the end of the suction pipe to be above the surface of the liquid (not dipped enough). Avoid curdling of the milk sample. Clean immediately if there is a sample curdled in the system. In mode Measurement, after starting the measurement, remove the sample holder and see if there is no milk poured back in the sample holder. 					
4 Sample Overheat	Sucked overheated	The analyzer is ready to measure the next sample. In order to avoid the future appearance of the same					

sample	error message, please, check the following:
	-The sample is prepared according the instructions and its temperature does not exceed the maximum permissible sample's temperature.
	-Complete the procedure for checking the analyzer in case of error message Empty Camera .

8. Making corrections and recallibration of the device

In the process of work with the analyser there is a possibility the results to start differing between the data for some of the measuring parameters when measured with the milk analyzer and the corresponding reference method of analysis (Gerber for fat, Kjeldhal for proteins etc). In order to establish the possible discrepancy and to correct the readings of the milk analyser do the following:

8.1. Taking samples and preparation of samples for checking the accurracy of the milk analyser, making corrections and recalibration

This is a basic moment for the correct checking the accuracy of the analyser and for making correct and precise correction and calibration. It is accomplished according Appendix Sampling and preparation of samples for verification the accuracy of the milk analyzer, making corrections and recalibration.

8.2. Determination the type of the discrepancy:

8.2.1. Making measurements

Make measurements with different samples (not less than 3) with known values of a separate parameter (for example fat content), determined by the known reference methods of analysis (for example Gerber's method for determination of fat content). For more accuracy it is recommended among these samples to be also such with values, close to the lowest and highest bounds for the measured parameters.

Make 5-time measurement for each of the samples. Calculate the average value for each sample parameter, without taking into consideration the first measurement for each sample.

8.2.2. Analysing the measurement results

Make comparison between the values of the parameter from the reference sample and measured with the analyser. Make analysis of the difference received.

8.2.2.1. If the received differences are relatively constant value for samples with different content of the analysed parameter, it is necessary to make correction.

For example

M% of the reference samples:	2,20	3,00	3,80	4,60	5,20
M%average when measuring					
with the analyser:	<u>2,38</u>	3,174,01	L4,795,4	2	
Difference:	0,18	0,17	0,21	0,19	0,22

Conclusion: the difference is relatively constant value and correction is possible to be done with -0.2 % (see Corrections, p6.3.3)

8.2.2.2. If the differences are not a constant value it is necessary recalibration to be done.

For example.

M% of the reference **Conclusion**: It is obvious that the difference is variable value and recalibration have to be done.

8.3. Making corrections

8.3.1. Possible corrections, limits and changing steps

Every parameter from each calibration may be separately corrected. Below is the table with possible corrections, limits and changing steps:

Parameter	Increasing	Decreasing	Step
FAT	0.95%	0.95%	0.01%
SNF	4.75%	4.75%	0.05%
Density	4.75%	4.75%	0.05%
Lactose	0.95%	0.95%	0.01%
Salts	0.95%	0.95%	0.01%
Proteins	0.95%	0.95%	0.01%
Added water	9.00%	9.00%	1.00%
Sample's temperature	9.90°C	9.90°C	0.1°C

Corrections are described in a dedicated software tool.

8.3.2. Making verification

After the corrections are made put the milk analyser in working mode and make several times measurement of reference samples with known values of the corrected parameter. If the difference between the values of the parameter from the reference methods and milkanalyser are in the limits for the parameter it may be considered that the correction is successfully made.

If the discrepancy between the measurements from the milk analyser and classical methods is bigger than is necessary to make second correction according above described way.

If after the second correction the results are unsatisfactory we recommend making a calibration of the analyser.



When making corrections or calibrations be 100% sure in the accuracy of the reference methods result.

8.3.3 Recalibration of the device

*Note: Before starting calibration please read appendix 1 of this manual. To start the procedure go to Menu -> Special Modes -> Calibration The window opens:

File						
					Getting Re Ready to s	əady ^ start
		FAT	3.21	Salts:	0.63	
		SNF	7.62	Lactose:	4.19	
		Density	27.07	Protein:	2.79	
Cancel	1		Measure	Einisk		Delete
Fat SNF	Den Lac	Salts Pro	N1 N2	Ampl		Delete
	00470407					•
Save File:	20170425	.txt Mea	asure Count:	Reset	Version	: 60 25/1/17

Input the parameter values for the high sample. By pressing the "Measure" button start a measurement

File											
										Getting Ready to	Ready o start
				FA	T	3.21		S	Salts:	0.63	
				SN	IF	7.62		Lac	tose:	4.19	
52				Dens	sity	27.07		Pro	otein:	2.79	
55	sec										
											Delete
Fat	SNF	Den	Lac	Salts	Pro	N1	N2	Ampl			
•											,
Save File:		20170)425.txt	t	Mea	sure Co	unt:		Reset	Versio	n: 60 25/1/17

Finished measurements are added to the table

File											
										Getting F Ready to	Ready start
				FA	λT	3.21		S	alts:	0.63	
				SN	١F	7.62		Lact	tose:	4.19	
	0			Dens	sity	27.07		Pro	tein:	2.79	
	0 sec)								
Car	icel					Measur	e		Finish		Delete
Fat	SNF	Den	Lac	Salts	Pro	N1	N2	Ampl			
3.21	7.62	27.07	4.19	0.63	2.79	61246	61813	17102			
3.21	7.62	27.07	4.19	0.63	2.79	61244	61810	17180			
3.21	7.62	27.07	4.19	0.63	2.79	61243	61806	17395			
3.21	7.62	27.07	4.19	0.63	2.79	61241	61807	17321			
•						III					Þ
Save File: 20170425.txt Measure Count: Reset Version: 60 25/1/17											

Stir well the sample with high fat value. Temperature of the milk sample has to be between 18-22°C.

You have to measure high fat sample 5 times. Results are shown on the computer's display.

Delete the first two measurements – mark the value you want to delete and press Cut. Check if the results in the columns N1 and N2 are stable and in the limits 60000-65000. For stable results are taken those with difference 2x-3x one from another. For example:



If results for N1 and N2 are with bigger difference than 3x, check milk's acidity (it has to be normal one); stir the sample and measure once again 5 times; in case you receive again not stable results, address our service.

If the results are normal, write in fields **Calibration sample – parameter values** the new values of low fat sample.

Stir well the sample with low fat value. Temperature of the milk sample has to be between 18-22°C.

You have to measure low fat sample 5 times. Results are shown on the computer's display. After the fifth measurement of the low fat sample

Fil	е											
											Getting F Ready to	Ready 6 o start
					FA	T	1.21		S	Salts:	0.75	
					SN	IF	7.62		Lac	tose:	4.98	
					Dens	sity	32.6		Pro	otein:	3.32	
	() sec										
	Canc	el					Measur	е		Finish		Delete
	Fat	SNF	Den	Lac	Salts	Pro	N1	N2	Ampl			
1	1.21	7.62	32.6	4.98	0.75	3.32	61241	61807	17093			
	1.21	7.62	32.6	4.98	0.75	3.32	61240	61806	17256			
	1.21	7.62	32.6	4.98	0.75	3.32	61242	61807	17163			
	1.21	7.62	32.6	4.98	0.75	3.32	61240	61806	17285			
	1.21	7.62	32.6	4.98	0.75	3.32	61242	61807	17256			
•							III					۱. Electric de la construcción de la const

Delete first 2 measurements. Check the stability of the results in columns Nº1 and Nº2

In the field **Calibration sample – parameter values** enter 0, because the next sample will be water.

Pour water. Measure water 5 times without changing it. The water has to be opaque. This is normal for calibration with water.

Check stability of the results in columns N $^{0}1$ and N $^{0}2$

Delete the first 2 measurements with water.

After finishing water measurement check carefully the results and if they are stable, click **Finish**

Confirm by pressing **the calibration channel** and the program will download the new calibration coefficients in the analyser.

When the analyzer is connected the next time it will be ready for work with those milk types you've just calibrated it for.

Checking the calibration

Switch on the calibrated device.

Make sure it shows the same serial number as this already

calibrated. For checking use the third sample with medium

FAT content. Measure the milk 5 times in the mode you've calibrated it.

In case that the device is not connected towards printer write down the results. Ignore the first two results. The rest three could not differ one from another more than 0,05% FAT, 0,07% SNF, 0,7% Density.

9.Basic printer characteristics

9.1. Basic printer characteristics

printing method	Thermo
Type of paper	Termorulony sensitive to the heat of the
	party outside
Recommended types of paper	From 55 g / m2 to 70 g / m2
Width	57,5 mm ± 0,5 mm

Printer control panel



- 2. Paper roll compartment
- 3. Printout button
- 4. Cover opening button

Changing the roll

To change the paper rolls proceed as follows:


1. Press the button as it is shown on the figure.



- 2. Position the paper roll making sure it unrolls in the proper direction.
- 3. Tear off the paper and close the cover

Sequence of actions when printing

The embedded printer is automatically switched on with switching on the device. After finishing the measuring procedure the results are automatically printed out. If there is a need of repeated printing of the same results the operator has to press the button Print on Main Display (each pressing the button Print - the results will be printed out). If you do not need to print out the results after the measurement is completed, just open the lid of the printer.

10. Additional possibilities of the analyser

10.1. Connecting to **12** V DC power supply.

If there is a need the analyser to work on place without electrical supply available, then it could be powered by car battery or other 12 V DC external power supply. Use the 12 V power supply cable (art. number 30030).

10.2. Measuring high fat samples (option).

The standard device measures samples up to 25% fat.

On customer's request, the device could be produced with possibility to measure samples up to 45% fat. The customer can choose which calibration to have this possibility and which not, as well as during the process of exploitation to change the measuring mode i.e. to pass from measuring normal fat percentage towards high and vice versa.

What the operator sees during these passes is the difference in the speed of sucking the sample. For that purpose, the high-fat sample has to be preliminary heated up to 30C + - 3C.

If the device is with High Fat option, the operator can select which calibration to run at a high fat speed. By going to Menu->Settings->High Fat control the operator can select which calibration to run at high fat speed.



Select the calibration you need so you can edit the status

Supplier Number:		KG:	high fat Please select c high fat	alibration for	• •
Cancel	UH		Sheep		
Time Sup N Kg Cal Tem	p Fat SNF	Den Pro Lac	Ad W Fr P	Cond pH S	alts
•	m				F F

By clicking on the calibration the details windows appears and the operator can enable or disable high fat measurement for this calibration

Supplier Number:	KG	: high f Pleas high f	at e select calibrat at	ion for
Calibration 3 High F	at: Enable	d		
Cancel]	Disable		
Time Sup N Kg Cal Temp	Fat SNF Den	Pro Lac Ad	W Fr P Cond	pH Salts
<[III.			•

By pressing the appropriate button, the operator can select the type of measurement and exit the menu. If you change the type of measurement of the calibration to do calibration for the corresponding speed. When calibrating for measuring samples with high fat content, before proceeding with the calibration procedure, the operator must select the menu ->Settings ->High Fat speed for calibration



And the operator should press on "Turn On"

Supplier Number:	KG	:	high fat Please select calibration for high fat	
Turn On		Т	urn OFF	

With this the device is ready for calibration with high fat samples. The user should start calibration procedure by going to Menu->Special Modes->Calibration

10.3 Saving and Archive of measurement data

After each measurement is finished, the results are stored in file with nameyyyyMMdd.txt, in C:/LSdata/import, where yyyy is the year, MM is the month and dd is the day, in which the measurements were made. This means there is such file for each day in which measurements were made. The structure of the file is as following:

Date 'tab' time 'tab' supplier 'tab' litres 'tab' serialNumber 'tab' calibration number 'tab' temperature of the sample 'tab' fat 'tab' snf 'tab' density 'tab' protein 'tab' lactose 'tab' Added water 'tab' salts 'tab' conductivity 'tab' ph 'tab' freezing Point,where 'tab' is separator between values

At the same time the results can be sent to a cloud server via Internet. So if the client provides protocol for communicating with their cloud server, the program will automatically send the results, when the measurement is made. At the moment the machine sends the results to a test cloud server using HTTP REST API, so if given a proper REST API address, the results can be easily redirected to this address. After finishing the measurements for a shift or group of measurements of operators choosing, the inbuilt data collection program can be started, which comes installed on every tablet – LSAN-DB (username:1, password: 1), which manual LSANDB_UserGuide.pdf can be found in c:/LSdata

10.4 Possibility of more than 3 calibrations

Lactoscan MCCW allows the user to have more than the 3 calibrations that are inbuilt in the device.

At any moment the analyser has 3 active calibration, which names are written on the middle 3 buttons of the Control Buttons (Operation manual section 3.2.1)



Each Calibration made in this device can be saved to the tablet, which is operating the analyser(operation manual, section 5 – Settings – Save Calibration)

File	Menu	Tools	Meas	uremen	t Hel	р								
	Speci	al Mod	es 🔸	ber:			KG:							*
	Corre	ections	•											
	Settir	ngs	•	Set C	alibrat	tion N	ame							-
	Tests		۲	Save	Calibr	ation						I AC 0	0.00	
	Clear	ning		Resto	ore Cal	ibrati	on							\prec
	Final	Clean		Auto	print E	nable	/Disabl	e		0 (_	ADW 0	0.00	
	pH C	lo Mete	er ⊦	Form	nat Res	ults						_		
	FrP 00	.00		High	High Fat Control						C	OND (00.00	
				High	High Fat Speed For Calibration				า					
	Print		С	ow		U	HTass			Sheep		C	Cleanir	ig i
Time	Sup N	Kg	Cal	Temp	Fat	SNF	Den	Pro	Lac	Ad W	Fr P	Cond	pН	Salts
•							11							

At any moment the user can Restore a calibration on one of the 3 calibration channels (operation manual, section 5 – Settings – Restore Calibration)



This way the user can have any number of calibrations saved to the tablet and when needed Restore(replace) one calibration with another

Appendices

Appendix 1: Preparation of samples for milkanalysers' calibration

For calibration are needed samples of cow milk with the following parameters:

		Low Fat	High Fat	Middle	
1	Cow	2,2%	5,2%	3,6%	

For the calibration are needed:

- 1. Distilled water
- 2. Min. 3 milk samples with known values for fat, SNF, protein, density, lactose, salts.

Calibration samples have to be with low, middle and high values of the analyzed components. Samples have to be representative for given milk type. Volume of the sample has to be enough for making min 5 measurements for each sample – not less than 1,00 l. Changes in the analyzed parameters in the samples have, if possible, to cover the whole measuring range – i.e. used samples to be with low, middle and high content of the analyzed components.

Methods of milk samples preparation for calibration.

For milk sample with middle value of the analysed components we recommend to use milk taken from not less than 10 animals from most common in the region breed.

Sample with low and high value are prepared on the following way:

- 1. Pour the fresh milk with FAT at about 3.7% in a separating funnel.
- 2. Leave the funnel with the milk in refrigerator for 12 hours at temperature +5-+8 ° C.
- 3. Draw the substratum of the separated milk in a vessel, mix it well, pour it and heat it in water-bath up to 20°C.
- 4. Pour the upper layer in another vessel.
- 5. Determine the concentration of the measured components (FAT, protein, SNF, density, lactose, solids) by using certified methods.



The analyser's accuracy depends only on the correctness of the chemical analysis of the components in the samples and the normal acidity during calibration! It is recommended the first cow milk sample with low fat content to be with the following parameters:

2-2,3% FAT; 8.7-9% SNF; 3,3-3,5 % Protein; 4,8-4,9% Lactose; 0,75 Salts; 1030-1033 kg/m3 Density.

The second cow milk sample with high fat content to be with the following parameters:

5-5,3% FAT; 8.4-8,79% SNF; 3,1-3,2% Protein; 4,6-4,7% Lactose; 0,7 Salts; 1028-1029 kg/m3 Density.

If, after milk's separation you do not obtain samples in the requested range, then, by adding milk with high fat value into the low fat milk sample you can obtain necessary value-2,3%

Analogous to this, by adding low fat milk sample into a milk sample with high fat value you may receive 5,3%

Samples with medium values are received by mixing low fat and high fat samples in necessary proportion.

If there is a need of longer sample storing they have to be preserved; the most commonly used preservative is potassium dichromate (K2Cr2O7) - 1 g for 1 000 ml.

When using samples, stored shortly, preliminary pour the sample from one vessel to another in order to distribute the milk components evenly paying attention not to form foam in the sample.

When the samples are stored for a longer period it is recommended to warm it up to 35-45 °C, and the vessel to be shaken carefully. In case that there is a cream stuck on the vessel's surfaces – remove it. The sample is poured from vessel to vessel several times and is cooled down (advisable to 20 °C /.



If there is separated liquefied fat or white particles with irregular form on the vessel's walls reliable results could not be received.

Because it is very difficult both lactose and salts to be measured but they are substantial and influence in great extend when determine added water. That's why it is better both lactose and salts to be calculated by using SNF results.

The milk must be for sure without added water.

If you are unable to make the analysis of milk in certified methods in a pinch you can use the following formulas:



DETERMINATION OF THE BASIC PARAMETERS IN THE MILK SAMPLE BY USING FORMULAS IS NOT AS PRECISE AS USING THE ARBITRARY METHODS, BUT IS SUITABLE FOR USAGE IN FIELD WORK.

1. Determination some of the parameters by formulas

There is dependence between the different parameters in milk and its density, which may be expressed with mathematical equation. On this base different formula, tested and confirmed by the classical laboratory methods for analysis, are developed. We recommend the following:

2. SNF determination.

For determination of SNF the correlation dependence exists between the milk's density, fat and SNF in the milk. When the density and the fat are known, the SNF can be calculated.

There are several formulas with different applicability.

A/ When the salts and fat are known

SNF is calculated by subtracting the fat percentage from the salts.

SNF = Salts - F(%)

Where

Salts – salts in (%),

F – fat content in (%),

This formula is used for determination of SNF in whey, buttermilk, and cream.

B/ Known quantity of fat and density (most commonly used method when maximum accuracy is needed).

We recommend the following formula:

$$SNF = \frac{0,075 * F\% + 100 - 100 / density}{0,378}$$

This is a universal formula and actual for milk of almost all kind of cows and sheep all over the world.

3. Determination of lactose content

We recommend the following formulas:

A/ for cow milk

Lact. = SNF * 0,55 (%)

Where

SNF – content of SNF in percentages (%),

0,55 – constant coefficient.

B/ for sheep milk

Lact. = SNF* 0,45 (%)

Where

SNF –solids-non-fat content in percentages (%),

0,45 – constant coefficient.



This is an actual coefficient for sheep breeds on the territory of the Balkan Peninsula.

4. Determination of salts content

We recommend using the following formulas:

A/ for cow milk

Salts = SNF* 0,083 (%)

Where

SNF – solids-non-fat content in percentages (%),

0,083 – constant coefficient.

B/ for sheep milk

Salts = SNF * 0,075 (%)

Where

SNF – solids-non-fat content n percentages (%),

0,075 – constant coefficient.



This is an actual coefficient for sheep breeds on the territory of the Balkan Peninsula.

5. Determination of total proteins content

We recommend using the following formulas:

A/ for cow milk

Protein = SNF * 0,367 (%)

Where

SNF - solids-non-fat content in percentages (%),

0,367 – constant coefficient.

B/ for sheep milk

Protein = SNF * 0,475 (%)

Where

SNF – solids-non-fat content in percentages (%),

0,475 - constant coefficient.



This is an actual coefficient for sheep breeds on the territory of the Balkan Peninsula.

Appendix 2. Freezing point determination

1. Methods for determination.

The milk analyzer determines the freezing point of each sample and the quantity of added water. The milk analyser does not measure the freezing point, but calculates it from the components it depends on. The basic components in the milk are water, solids, lactose, FAT, proteins, minerals (salts) and acids. The freezing point depends only on the diluted in the milk components and quantity of the solvent (in the milk it is water). The ultrasonic technology allows direct measurement of FAT, proteins, lactose + salts (the soluble components, only influencing the freezing point), and the quantity of the solvent in % is determined by 100 % – total solids %, total solids = lactose % + FAT % + proteins % + salts % + acids %.

Without understanding the meaning of the freezing point – determined or shown from the milk analyzer added water result easily may lead to a mistake for the value of this parameter.

2. The basic freezing point.

Milk freezes at lower temperature than water. The average freezing point of the raw milk in the most regions is at about -0,540°C. The average reading for your region is called "basic" freezing point.

The freezing point of milk is a "physiological constant". This does not mean that it will not vary. In fact feed, breed, season, time of lactation, climate, whether the sample is taken at the beginning, middle or end of lactation – all these factors will have an effect on the freezing point of the individual sample. This means that there is an average value of all these numbers. The more samples used in obtaining this average, the more reliable it is as a base. Or the basic freezing point is an average of freezing points of milk, taken from many cows. When a laboratory checks a producer, it is only comparing the average of the producer's cows against a larger area average.

The Health authorities establish the basic freezing point or agriculture departments in some regions, sometimes by universities, separate dairy producers, or their associations. Frequently, tolerances have been established on top of a basic freezing point to allow some variations in the milk as well as device or operator variations.

Withoutmentioningthebasicfreezingpoint, the Association of Official Analytical Chemists now recommends an upper limit freezing point at $-0,525^{\circ}C$ (2,326 standard deviations above the most recently determined North American average of $-0,5404^{\circ}C$), below which there will be at 95% confidence that will show 99% of all freezing point determinations on unwatered milk:

"if the freezing point is -0,525 °C or below, milk may be presumed to be free of water or may be confirmed as water free by tests, specified below. If the freezing point is above -0,525 °C, milk will be designated as "presumptive added water" and will be confirmed as added water or added water free by tests specified below. Evaluate extreme daily fluctuations in the freezing point of herd, pooled herd, or processed milk for presence of added water".

"Presumed added water", as described above, must be "confirmed" by means of tests on authentic milk samples obtained as specified in the AOAC METHODS.

After determination the freezing point of your sample via the milk analyzer, the added water is calculated using the following formula:

$$AddedWater = \frac{FrPo \operatorname{int}_{Base} - FrPo \operatorname{int}_{Calc}}{FrPo \operatorname{int}_{Base}} * 100[\%]$$

Where:

FrPointBase is the basic freezing point

FrPointCalc is measured freezing point

Note:

If the freezing point is not correctly determined, the result for the added water is not valid. In this case results for FrPoint and AddWater are not shown on the display and on the printout from the printer. If the density of the measured sample is 0, the result for AddWater is not valid and is also not shown on the display and the printouts.

Sample:

First variant

If you've entered for milk analyzer basic freezing point -0.520°C (according article 5.9 of the EU Milk Hygiene Directive 92/46/EEC), measured freezing point -0.540°C, using the above pointed formula you'll receive -3,8%. Because it is not possible the added water to be negative value, the milk analyzer indicates 0% added water. The reason for this is the tolerance in the basic freezing point, reasons for which are described below.

If in the same milk we add 3,8% water, and the basic freezing point is the same, the milk analyzer will measure freezing point -0.520°C, and will indicate again 0% added water.

Second variant

If you've entered for the device basic freezing point -0.540°C, measured freezing point -0.540°C, the milk analyzer will indicate 0%. When you add 3,8% water, the device will indicate 3,8%-added water.

From the above mentioned follows that it is very important to enter correct basic freezing point in the device.

The device's results for added water may give information about doubt of added water in the milk and the exact value of this added water may be determined after a

"cowshed sample" is taken and the result for the freezing point, measured by the milk analyzer of the "cowshed sample" is entered as basic freezing point in the formula for calculation of added water.

Appendix 3. pH measuring

1. General information

PH probe is a unit, measuring the solution acidity or alkalinity degree. It is measured on scale of 0 to 14. The term pH is derived from "p", the mathematical symbol for the negative logarithm, and "H", the chemical symbol of Hydrogen. The formal definition of pH is the negative logarithm of the Hydrogen ion activity.

2. pH Electrode

For pH measurement the milk analyzer needs a combination electrode, compatible with most pH electrodes that have BNC connectors and zero potential (the pH where the mill volt output of the electrode equals 0) near 7 pH.

2.1. Electrode part

The electrode is the most important part of the pH measurement. The electrode glass membrane is fragile and must be handled with care. To protect the glass membrane and to maintain activation, a protective rubber cap containing a suitable storage solution covers the glass membrane.

2.2. Electrode care & Electrode maintenance

pH Electrodes are susceptible to dirt and contamination and need to be clean regularly depending on the extent and condition of use. At no time should one touch or rub the glass bulb as this causes the build-up of electrostatic charge.

2.3. Storage

For best results, always keep the pH bulb wet. An optimal storage solution for combination electrode is pH 4 buffer with 225 grams of KCl per liter. Table salt, NaCl, can be used if KCl is not really available. Other pH buffers or tap water are also acceptable storage media, but avoid storage in de-ionized water. The protective rubber cap filled with the buffer solution provides ideal storage for long periods.

2.4. After Use

After measurement is completed, follow the sequence below for storage.

- Wash the electrode and reference junction in de-ionized water.

- Close the refilling hole by returning its rubber sleeve or stopper cap. (Necessary for only refillable electrode).

- Store the electrode as mentioned above (see section Storage).

2.5. Electrolyte Replacement (for refillable electrode only).

The reference electrolyte needs to be refilled when the electrode has been used for a long period, or when the internal electrolyte has dried up. To accomplish this, follow the procedure described below.

- Remove the protective rubber cap or sleeve;
- Remove the protective rubber sleeve to expose the filling port of the electrode;
- Remove the old reference electrolyte with a syringe;
- Fill the new reference electrolyte.

2.6. New electrolyte preparation:

- Open the KCl container;
- Add in de-ionized water until it reaches the level of 20 ml;
- Close the container and shake it to dissolve the KCl;
- Add in fresh electrolyte until it reaches the level of the refilling port. The reference electrolyte used should be 3M(Mol) KCI;
- Replace the rubber sleeve.

2.7. Re-use the electrode.

- Rinse the liquid junction with de-ionized water.



If these steps fail to restore normal electrode response, you may attempt to rejuvenate it (See: Electrode Rejuvenation).

2.8. Electrode cleaning

Electrodes which are mechanically intact can often be restored to normal performance by one or combination of the following procedures.

- Salt deposits:

Dissolve the deposit by immersing the electrode in tap water for ten to fifteen minutes. Then thoroughly rinse with de-ionized water. Wash the electrode pH bulb in a little detergent and water. Rinse electrode tip in with de-ionized water.

- Oil/Grease films:

Wash electrode pH bulb in a little detergent and water. Rinse electrode tip with deionized water.

- Clogged Reference Junction:

pH electrodes have junction, which allows the internal fill solution of the measuring electrode to leak out into the solution being measured. The junction can become clogged by contamination in the solution. If a clogged junction is suspected it is best to clear the junction.

Heat up the diluted KCl solution to 60-80°C. Place the sensing part of the pH electrode into the heated KCl solution for approximately 10 minutes. Allow the electrode to cool while immersed in some unheated KCl solution.

- Protein Deposits

Prepare 1% pepsin solution in 0.1 M HCI. Allow the electrode to stand in this solution for five to ten minutes. Rinse the electrode with de-ionized water.

2.9. Electrode activation

Generally, if the procedure of storage and maintenance had been closely followed, the electrode can

be used immediately. However, should the electrode response become sluggish, it may be possible that the bulb has dehydrated.

The bulb can be dehydrated by immersing the electrode in an ideal storage solution (e.g. buffer pH 4 solution) for 1-2 hours. If this fails, the electrode may require reactivation. If the above procedure does not reactivate the electrode to acceptable status, try rejuvenation the electrode by following the procedure outlined below.

2.10. Rejuvenation Procedure

Dip and stir the electrode in freon or alcohol for 5 minutes.

Leave the electrode in tap water for 15 minutes.

Dip and stir the electrode in concentrated acid (HCI, H_2S_4) for 5 minutes.

Leave the electrode in tap water for 15 minutes.

Dip and stir in strong base (NaOH) for 5 minutes.

Leave the electrode in tap water for 15 minutes.

Test with standard calibration solution.

Finally, test with standard calibration buffer solution to see if the electrode yields acceptable results. You may repeat again for better response (maximum 3 times). If the response does not improve, then the electrode has completed its useful life. Replace with a new electrode.

2.11. Electrode Lifespan

pH electrodes have a finite lifespan due to their inherent properties. How long a pH electrode will last will depend on how it is cared and the solution it is used to

measure. Even if an electrode is not used it still ages. Electrode demise can usually be characterized by a sluggish response, erratic readings or a reading, which will not change. When this occurs an electrode can no longer be calibrated. pH electrodes are fragile and have a limited lifespan. How long an electrode will last is determined by how well is maintained and the pH application. The harsher the system, the shorter the lifespan. For this reason it is always a good idea to have a back-up electrode on hand to avoid any system down time.

3. Buffer Solutions

Buffers are solutions that have constant pH values and the ability to resist changes in that pH level. They are used to calibrate pH measurement system.

PH buffer solution description (Pharmacopoeia standard)

Description	pH 7.00±0,01/20°C	pH 4.00±0,01/20°C
Composition	Potassium dihydrogen phosphate, Di-sodium hydrogen phosphate	Borax, Sodium hydroxide solution
Temperature parameters	10ºC - 7.06	10°C - 4.00
	25°C - 6.99	25°C - 4.00
	20ºC - 7.00	20°C - 4.00
	30°C - 6.98	30°C - 4.00
	40°C - 6.95	40°C - 4.00
	50°C - 6.91	50°C - 4.05

Use only this types standard buffers for calibration!

4. pH Electrode Calibration

pH Electrodes are like batteries; they run down with time and use. As an electrode ages, its glass changes resistance. For this reason, electrodes need to be calibrated on a regular basis. Calibration in pH buffer solution corrects for this change.

Calibration is an important part of electrode maintenance. This assures not only that the electrode is behaving properly but that the system is operating correctly.

Usually pH meters require calibration at 3 specific pH values. One calibration is usually performed at pH 7, second and third are typically performed at pH 4 and pH 10.



It is best to select a buffer as close as possible to the actual pH value of the sample to be measured. Use standard calibration buffers that the temperature and the sample solution are the same.

Use the operation manual for the corresponding pH meter.

To start the calibration procedure select Menu -> pH & Co Meter ->pH Meter Calibration



Input the values of the buffers, that are going to be used and press "Set"



Follow the instructions, received in the field

File Menu Tools Measuremen	t Help			
Supplier Number:		KG:		Ready to start Put High buffer and press Start
Buffer I	ligh:	7.4		1
Buffer	Low:	4.68		
Cancel		Start		

Wait for couple of measurements and press "Stop". The buffer should be returned in the glass.

After that water should be placed in the machine

File Menu Tools Measuren	nent Help			
Supplier Number:		KG:	2.010 Put Water and Press Start Ready to start	* *
Buffe	ər High: 7.4			
Buff	er Low: 4.68			
Cancel	Sta	irt		

Next step is the second buffer:

File Menu Tools Measurement Help)	
Supplier Number:	KG: Put second buffer and press Start Ready to start	•
Buffer High:	7.4	
Buffer Low:	4.68	
Cancel	Start	

Wait for couple of measurements and press "Stop". The buffer should be returned in the glass.

Again place water in the machine

After the rinsing the machine returns "Procedure completed message" and the operator should press "Finish" button to complete the procedure

File Menu Too	ls Measurement Help		
Supp	lier Number:	KG:	Ready to start Procedure Completed Ready to start
	Buffer High:	7.4	
	Buffer Low:	4.68	
	Cancel	Finish	

For Sensorex pH electrodes, originally supplied with the milk analyser read the following information:

Temperature compensations

The output of pH electrodes varies with temperature in manner, predicted by theory. When needed, Sensorex can supply electrode holders with build-in automatic temperature compensators. The need of automatic compensation depends on the temperature variation, the pH value being measured. At pH of about 7 there is no error due to temperature and, of course, at a constant temperature there is no error. As shown in the following table, the pH error due to temperature is a function of both the temperature and the pH value being measured. At a pH of about 7 there is no error due to temperature and the pH value being measured. At a pH of about 7 there is no error. The more the temperature and, of course, at a constant temperature there is no error. The more the temperature changes from the ambient calibration temperature and the more the pH departs from 7 the greater is the pH error.

°C	рН										
	2	3	4	5	6	7	8	9	10	11	12
5	.30	.24	.18	.12	.06	0	.06	.12	.18	.24	.30
15	.15	.12	.09	.06	.03	0	.03	.06	.09	.12	.15
25	0	0	0	0	0	0	0	0	0		0
35	.15	.12	.09	.06	.03	0	.03	.06	.09	.12	.15
45	.30	.24	.18	.12	.06	0	.06	.12	.18	.24	.30
55	.45	.36	.27	.18	.09	0	.09	.18	.27	.36	.45
65	.60	.48	.36	.24	.12	0	.12	.24	.36	.48	.60
75	.75	.60	.45	.30	.15	0	.15	.30	.45	.60	.75
85	.90	.72	.54	.36	.18	0	.18	.36	.54	.72	.90

pH temperature error table

0 pH Error Range

Less than .1 pH Error Range

5. PH helpful hints

For greatest accuracy in pH measurement, follow these guidelines:

Use the same technique to measure samples, which was used for calibration.

Be consistent with stirring rates, times and conditions.

Calibrate with buffers, which are close in temperature to that of the sample.

Calibrate the pH electrode regularly, e.g. once an hour for accuracy to within 0.01 pH, or once a day for accuracy to within 0.1 pH.

Use fresh buffers for calibrations. Avoid contamination of the stock buffer solution and do not use it beyond the expiry date.

Keep all connections dry.

Immerse the electrode far enough into the solution to insure the reference junction is below the surface.

Allow adequate time for the electrode to stabilize in standards and samples before taking a reading.

Clean the electrode periodically. Allow more time for aged electrodes.

Do not use the pH electrode in solutions of fluoride ion at low pH. This will etch the glass membrane.

Sulphide vapors can permeate the electrode wick and contaminate the reference element. Minimize contact in such environments and change the reference electrolyte frequently.

Lactic acidity and pH This table shows the relationship between the values of pH

and°⊤

۰T	Deviations	рΗ	°T	Deviations	Нα
•		average	•		average value
		value			5
Raw milk			Pasteurized milk		
16	6,74-6,70	6,72	16	6,68-6,64	6,66
17	6,69-6,65	6,68	17	6,63-6,58	6,61
18	6,64-6,58	6,62	18	6,57-6,52	6,55
19	6,57-6,52	6,55	19	6,51-6,46	6,49
20	6,51-6,46	6,49	20	6,45-6,40	6,43
21	6,45-6,40	6,43	21	6,39-6,35	6,37
22	6,39-6,35	6,37	22	6,34-6,30	6,32
23	6,34-6,30	6,32	23	6,29-6,24	6,26
24	6,29-6,24	6,25	24	6,23-6,19	6,21

6. Preparation for pH measurement

When the analyzer is with pH measuring option, it is received from the customer with pH probe packed separately and there's a stopper on its place. If you need to measure pH follow the procedure below:

- 1. Loosen the nut anti-clockwise (1).
- 2. Pull up the stopper (2)
- 3. Carefully place the pH probe paying attention not to remove the sealing O-ring (3)

4. Place the probe with the nut in the hole (C & D) and tighten it.

Preparation for work with the pH probe





It is very important to close the nut tightly, paying attention not to allow air to enter the system.

Placement of the pH probe





If you are working with the analyser regularly (each day) do not remove the probe after work.



If you 'll not use the analyser more than 2 days, you must take out the probe and to place the stopper back.

The pH probe must be stored separately as per the instructions of point 2.3 Appendix 3.





Designed for calibrating pH meter. For this purpose, use 2 reference buffer output to the screen as a Low buffer (eg 5.00 pH) and High buffer (eg 7.00 pH).

Use this procedure only if you have a sufficient number of buffers for calibration, since they can not be reused.

If you do not have enough buffers, place the probe in the containers near the unit as shown in the Fig.

7. PH measuring.

Measuring pH is an additional feature of the analyser and is optional.

Remove the protective rubber cap of the pH electrode. Take care to handle it appropriate in order not to be damaged. Use de-ionized or distilled water to rinse the electrode before usage. Fill in the sample holder with milk, put it in the recess of the analyser and dip the pH electrode into the milk sample, ensuring complete dip of the electrode in the sample. Stir gently for homogenization of the sample.

Measuring can be done in two modes:

Off line by starting the menu **pH & Co Meter | Measuring**, when the analyser works only as a pH meter.

On line automatic pH measuring, when measuring the rest of the sample's parameters.



When starting work with pH meter first connect the probe/sensor, and then the power supply of the device.

Having in mind the characteristics of the process of pH measuring it is necessary to dip the pH probe in the sample and then start measurement.

Appendix 4. Conductivity measuring

1. Method of determination.

Conductivity (or Electrolytic Conductivity) is defined as the ability of a substance to conduct electrical current. It is the reciprocal of the resistance.

In a healthy animal*, the mean value of electric conductivity is:

Milk type	Conductivity values
Cow milk	between 4 to 6 mS/cm (18°C);

Sheep milk	between 3 to 5 mS/cm (18°C);
Buffalo	between 2,5 to 5 mS/cm (18°C);

*These values depend on the geographical region, the breed and on other factors.

Milk conductivity changes on the concentration of ions in the milk:

Added water, sugar, proteins, insoluble solids	Decrease the ion's concentration. Milk conductivity decreases.				
Added salts	Increase the ion's concentration. Milk conductivity increases. Increase the ion's concentration. Milk conductivity increases. Often the milk is falsified by adding salt: towards milk with good characteristics: fat 4%, SNF 8,8, conductivity 4,5 are added salt and water. Then the results are changed to 3,2 and 8,8, conductivity 10. In other words adding water regulates the increased value of SNF and density till normal (within the boundaries/parameters) and even the fat is normal. By the values of these parameters may be determined if the sample is falsified, but the only characteristic, proving this is conductivity, which is out of boundaries nevertheless added water. But be careful, as the falsification is not the only possible reason for conductivity increasing. The other possibility is mastitis that's why we recommend using another (chemical) method for checking it.				
Significantly extreme value (6,5 - 13,00 mS/cm (18°C)	Should indicate the development of mastitis. Infections damage the tissue of the udder. This allows sodium and chlorine ions from the blood to be released into the milk. The concentration of ions in the milk is thereby raised, and it can more easily conduct an electrical current - the conductivity of the milk increases.				

Milk conductivity can be used as tests for degree of water evaporation in condense milk production.

Milk conductivity change notifies of powder (dry) milk solution rate.

2. Conductivity measurement

Conductivity measurement is additional possibility of the analyser and is delivered on customers request/

3. Corrections in conductivity measurement.

Please, follow the instructions on the display of the tablet.

4. Conductivity calibration buffer preparation

In order a standard buffer for conductivity measuring to be prepared follow the instruction below:

- Take the packet with the powder buffer.
- Carefully shake the packet in order to gather the powder at the bottom.
- Cut one end of the packet.
- Empty its content in a measuring mug with 1 l volume, paying attention all its content to be emptied.

For standard buffer: 5,02 ms – 3,056 r

- Add 600-700 ml distilled water, which was preliminarily deaerated in vacuum dryer or boiled and then cooled down to 20 °C.
- Shake the mug till the powder is fully dissolved.
- Add distilled water to the mark.

Basic concepts

User - any individual who is a user of the device and / or the software

Disclaim

Important legal information

The information contained in this user manual is subject to a disclaimer and is copyrighted. It is subject to a personal data protection system.

Discharge clause

The information contained in this manual is for informational purposes only and Milkotronic Ltd. is not responsible for them.

The information contained in the instruction does not constitute professional or legal opinion.

However, Milkotronic Ltd. does its best to distribute accurate and up-todate information, Milkotronic Ltd. can guarantee this. If we are notified of errors, we will do what is necessary to correct them.

Milkotronic Ltd. is not responsible for the use of data and information contained in this document. Users are required to take all necessary precautions before using this information, which they do at their own risk.

Milkotronic Ltd. is not responsible for the subjective perceptions and interpretation of the accuracy, completeness and usefulness of the information resources.

Milkotronic Ltd. shall be entitled to compensation for all damages, costs and claims of third parties as a result of violation of these General Terms and / or unauthorized use of the Services.

General Terms of Use

The devices purchased come with Remote Support included, in the form of a Remote desktop application. Remote maintenance is provided only with the prior consent of the User at a convenient time for the User and the representative of Milkotronic Ltd.

If the user becomes unwilling to use remote support, the user can uninstall / delete the Remote desktop application. As a result of this action, the Customer can not have any claims against Milkotronic Ltd. for maintenance.

By buying the device you agree to the terms of use and the software provided.

Copyright

© Milkotronic Ltd., 2019

Reproduction is authorized provided the source is acknowledged, unless otherwise specified.

If prior permission is required for the reproduction or use of text and multimedia data (pictures, etc.), that permission will void the general permission referred to above and clearly indicate any restrictions on use.

Warranty card

LACTOSCAN COMBO

The warranty is valid for a period of 1 (one) year. The incorrect working, transporting and storage make the warranty invalid.

Serial Nº

Date of purchase:

Password:

Distributor:

Signature:

Stamp:

",Milkotronic" Ltd preserves the right to choose the methods for checking the device in order to establish the validity of the warranty. Devices with expired warranty are not subject to a free warranty service.

Service card

Customer:

Service report:

Service entry date	Damage	Date of receipt	Signature

Contacts

For more information or technical support, visit our website <u>www.lactoscan.com</u>or contact us:

Headquarters and Service:

4, NarodniBuditeli Street 8900 Nova Zagora BULGARIA tel/fax: + 359 457 670 82 office@lactoscan.com