

**ULTRASONIC MILK ANALYSER**  
**DIRECT MEASUREMENT OF COLD MILK SAMPLES,**  
starting at 5°C  
No need of periodical calibration

# Lactoscan Standard

**60 sec, 30 sec measurement**

*Inbuilt printer (option)*

*Inbuilt keypad (option)*

*In-flow pH meter (option)*

*Multilingual display*

**Results**  
F=03.02 S=07.93  
D=28.55 P=02.90  
L=04.35 W=04.23

**结果**  
脂03.32 固07.69  
密27.38 蛋02.81  
糖04.21 水07.30

**Жир=03.26 1**  
СОМ=07.68  
ПЛТ=27.38  
Бел=02.81



**USB flash drive (option)- unlimited storage**  
of measurement results in Excel file,  
ready to be processed in different reports

## SPECIFICATIONS Lactoscan S:

| Parameter             | Measuring range                     | Accuracy                |
|-----------------------|-------------------------------------|-------------------------|
| Fat                   | from 0.01% to 45%                   | ± 0.06%                 |
| SNF                   | from 3% to 40%                      | ± 0.15 %                |
| Density               | from 1000 to 1160 kg/m <sup>3</sup> | ± 0.3 kg/m <sup>3</sup> |
| Protein               | from 2% to 15%                      | ± 0.15%                 |
| Lactose               | from 0.01 % to 20%                  | ± 0.20%                 |
| Added Water content   | from 0 % to 70 %                    | ± 3.0%                  |
| Temperature of milk   | from 5°C to 40°C                    | ± 1°C                   |
| Freezing point        | from - 0,400 to - 0,700°C           | ± 0.005 °C              |
| Salts                 | from 0,4 to 4%                      | ± 0.05%                 |
| pH (option)           | from 0 to 14                        | ± 0.05%                 |
| conductivity (option) | from 2 to 14 mS/cm                  | ± 0.05 (mS/cm)          |

Milk analyzer Lactoscan makes quick analysis of milk and liquid dairy products:

Cow milk (25%)  
Sheep milk  
Buffalo milk

Whey  
Cream (up to 45%)  
Skimmed milk (0,01% FAT)

Concentrated milk (up to 1160kg/m<sup>3</sup>)  
Recovered milk  
Etc.

And can be calibrated by the customer with specific samples of:  
Yoghurt, Flavoured milk, Ice-cream mixtures etc.

The measurement accuracy  
not dependant on milk's acidity  
High-end ultrasonic technology  
for analyzing any kind of milk

- Easily readable results
- Possibility of connecting pH probe
- Integrated pH and conductivity measuring
- Self calibration without computer

## KEY FEATURES:

- User-friendly: simple in operation, maintenance, calibration and installation
- Portable and compact design
- Very small quantity of milk required
- Low power consumption
- No use of hazardous chemicals
- One year full warranty
- Adjustments of the measurement accuracy can be done by the RS 232 interface
- ESC POS Printer Support
- Two samples self-calibration



## MEASURING PARAMETERS:

Fat  
Solids-non-fat (SNF)  
Density  
Protein  
Lactose  
Milk sample temperature  
Added water  
Salts  
Freezing point  
Total Solids (option)  
pH (option)  
Conductivity (option)  
Ion meter (option)



**ENVIRONMENTAL  
CONDITIONS:**

Ambient air temperature  
10°C -40°C (option 43 °C)

Milk temperature 1°C -40°C

Relative humidity 30% - 80%

**ELECTRICAL  
PARAMETERS:**

Switching Adapter

Input:  
100 - 240 V ~1.6 A max.  
50-60 Hz

Output: +12 V 4.17A min.

Output power: 50 - 65 W

**MECHANICAL  
PARAMETERS:**

Dimensions  
(W x L x H)  
130x230x220 mm

Weight < 3 kg

Stainless steel  
cover box



**Lactoscan Standard**

60 sec, 30 sec measurement




pH measuring system (degree of acidity) (option)  
In-flow pH measuring system (degree of acidity)  
All parameters are measured by sucking one 25 ml milk sample.

Integrated pH meter



1. Power switch
2. 12 V output
3. pH probe input
4. 12 V input
5. Serial interface (RS232/ printer)

**MILK DATA Program:**  
milk collection accounting and management software.



Internet Remote Firmware  
Update (Download center)



Embedded function real date  
and time - registers the exact  
time for delivery of milk in  
the milk collecting center

Conductivity - built-in  
conductivity sensor gives high  
performance freezing point  
analysis - avoid poor quality  
milk due to added water;  
detects mastitis and  
falsifications with added salt

Visualisation  
out of range measurement

Milk conductivity measuring system (option)

Milk conductivity changes depending  
on concentration of ions in the milk.

Milk conductivity can be used for:

- test of udder health (detection of subclinical mastitis);
- control of water evaporation grade in condense milk production;
- determination of dry milk solubility rate.



ECS POS Serial Printer (option)  
Prints the result of the current  
measurement.

buy online on  
www.lactoscan.com



**MILKOTRONIC LTD**

Bulgaria, 6000 Stara Zagora  
167, Tsar Simeon Veliki, A, 25  
Tel: + 359 457 670 82,  
Fax: + 359 457 670 83

e-mail: salesdpt@lactoscan.com  
Eng, Spa, Port: +359 886 841 389  
e-mail: sales2@lactoscan.com  
Eng, Rus: +359 886 841 395  
e-mail: mainservice@lactoscan.com  
Eng, Rus: +359 888 811 017

e-mail: Prokopiev@lactoscan.com  
Eng: + 359 885 597 232  
e-mail: milkanalyzers@lactoscan.com  
Eng, Rus: +359 896 753 803

distributors: