


## DESCRIPTION

Independent, battery powered thermostatic device with a set for determination of antibiotics and inhibitors in the collected milk.

Prevents the risk a dairy processor to receive milk, in which milk acidic bacteria to develop (which is especially dangerous in case of yoghurt production).

Such device may be installed in the milk collection center or on the milk transportation tank vehicle.



## INTENDED USE

Test for the detection of inhibitors and antibiotics in milk. This is a rapid and simple method to check whether milk contains antibiotics in a concentration exceeding the Maximum Residue Limits (MLR).

For more information about the Eclipse Test for Detection of inhibitory substances in milk see:

[Eclipse50](#)

[Eclipse 100](#)