






## LACTOSCAN THERMOSTATIC DEVICES

Thermostatic device type	Picture	Dimensions	Standard set	Type of test / tube
LACTOSCAN TBS		W x L x H: 90x157x56 mm  Weight: 0,570 kg	Switching adapter: Input: 100-240V ~1.6 A max. 50-60Hz  User manual Red and Black cable for connecting to car/truck lighter.	BetaStar®
LACTOSCAN TCH		W x L x H: 90x157x59 mm  Weight: 0,720 kg	Switching adapter: Input: 100-240V ~1.6 A max. 50-60Hz  User manual Red and Black cable for connecting to car/truck lighter.	Charm® CowSide II
LACTOSCAN TDC		W x L x H: 90x157x65 mm  Weight: 0,910 kg	Switching adapter: Input: 100-240V ~1.6 A max. 50-60Hz  User manual Red and Black cable for connecting to car/truck lighter.	Delvotest®,
LACTOSCAN TET		W x L x H: 90x157x47 mm  Weight: 0,500 kg	Switching adapter: Input: 100-240V ~1.6 A max. 50-60Hz  User manual Red and Black cable for connecting to car/truck lighter.	Twinsenser, Eclipse 3G
LACTOSCAN TET combo		W x L x H: 90x157x52 mm  Weight: 0,650 kg	Switching adapter: Input: 100-240V ~1.6 A max. 50-60Hz  User manual Red and Black cable for connecting to car/truck lighter.	Twinsenser, Eclipse 3G, Eclipse Farm 3G
LACTOSCAN TSO (SPECIAL ORDER)**	On customers' request we can produce for you THERMOSTATIC DEVICE for ANY TYPE of TEST or COMBINED FOR SEVERAL, for ANY TYPE OF TUBE/VIAL or STRIP, with YOUR LOGO on the front panel.			